



Hill-Rom P8010 Surgical Stretcher - Fully Refurbished

Hill-Rom's flexible surgical stretcher (P8010) is ideal for any surgical setting, including outpatient, OR or freestanding surgery centers. The articulating head section and wrist rests are designed to suit a variety of surgeon preferences, regardless of the setting or type. A surgical rail allows easy attachment of gas delivery and drape support devices, PACU extenders and arm boards, and surgical hand tables.

Features:

- 700 lbs weight capacity
- 26" width
- Low height for safer stretcher ingress and egress
- Articulating head section for maximum flexibility
- Surgical rail for easy attachment of gas, extenders, and hand tables
- Dual-locking, four-wheel brakes for a strong, stable platform
- Brake/Steer pedals – head and foot
- Next generation Steering Plus™ system
- Siderail system (tuck-away)
- Power Washable
- New Mattress

Options/Accessories:

- Enhanced comfort mattresses
- IV poles
- Wrist Rest
- Armboards
- Restraint straps



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.