



Mindray M5 Ultrasound - Fully Refurbished

The M5 Ultrasound delivers a full range of functionality in a portable, laptop sized package. Weighing about 13 pounds, the Mindray M5 Ultrasound system is capable of producing brilliant 1D color images with enhanced contrast and spatial resolution through tissue harmonic imaging. Thanks to its easy format, the Mindray M5 can be easily transported and can provide a wide range of application services to meet more medical occasions including in the operating room, out in the field, in ambulances, etc.

Features:

- High-resolution 15" non-glare screen
- i-Touch - one button optimization
- Trapezoid
- Dual and Quad screen imaging
- Duplex / Triplex
- MBP- Increases resolution and frame rate by 4 times
- FTO- Increases soft tissue differentiation
- AVI- Distinguish tissue movement from blood flow
- i-Station- robust patient management system
- 80 GB hard drive (100k still images)
- DVD-R/RW
- 2 USB ports
- S-video output
- AC adapter & battery (4500mAh)
- Traveling case
- B-Mode/M-mode/PW/HPRF/Color Doppler/Directional Power Doppler
- Applications: OB/GYN, Urology, Vascular, Abdominal, Small parts, Basic cardiac



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.