



auxoMEDICAL

NUVO VisitOR Pro LED Light - New

The VistOR PRO Procedural Light features a large 16-inch, highly-polished faceted reflector which keeps shadowing to a minimum and maximizes the light output. This pre-focused light also combines a great depth of field, smooth, drift-free positioning and an easy-to-change bulb that doesn't require the use of tools.

Features:

- 64,500 Lux light intensity from a single quartz halogen bulb
- Color correct, extremely cool light output
- Large pre-focused pattern
- 36-inch depth of field
- Large faceted reflector provides shadow free illumination
- Color temperature: 4,000 K
- Depth of Field: 36 inches
- Bulb Life: 2,000 hours
- Bulb Voltage/Wattage: 12VAC/90W
- Available as a Ceiling Mount, Wall Mount, Freestanding Mobile Unit or Track Mounted on the Chuttle Track.

Available:

- Single Light Head
- Single Light Head + Monitor
- Single Light Head + Camera
- Dual Light Head
- Dual Light Head + Monitor
- Dual Light Head + Camera
- Triple Light Head
- Triple Light Head + Monitor
- Triple Light Head + Camera



Why purchase equipment from Auxo Medical?

At Auxo Medical, we don't take any decision lightly. Whether we focus on quality of new or refurbished medical equipment, repair savings, or improved operational cost-efficiency—our eye is always on the smallest details that can affect the level of care provided and your budget. Our consultants pride themselves in helping you plan and develop the most appropriate facility for your changing needs.