



Philips iE33 Ultrasound - Fully Refurbished

Keeping pace with the clinical needs of managing patients with cardiac disease, such as heart failure, valvular and coronary artery disease, and congenital heart disease, calls for exceptional levels of image quality, quantification and clinical performance. The iE33 intelligent echo system is addressing this growing requirement with complementary 2D and Live 3D imaging, quantification tools, and connected work ow solutions that help you answer the major questions related to cardiac disease management: structure, efficiency, size, and function.

Features:

- xMATRIX cardiac transducer with 3000 elements
- Automated Stress Echo
- Premium Image Quality
- Full Doppler functions and Live 3D TEE
- 20-inch high resolution articulating flat panel display
- Fully articulating control panel
- Live Xplane multiplane imaging
- Tissue Harmonics and Tissue Harmonic SonoCT
- XRES, iSCAN one-touch optimization, iFOCUS one-button focal zone optimization, iOPTIMIZE one-button optimization multifunction auto image optimization
- Integrated CD/DVD burner and export to AVI and JPG
- DICOM 3.0 and USB port for import/export of images
- On-Board workstation data management
- 8-speaker high-fidelity stereo audio
- Panoramic Imaging
- Channels scalable up to 57,000 channels
- Multiple triplex mode operation
- High Resolution interactive color touch panel with ambiance adjustments
- Pull-out alphanumeric keyboard, 2D/M-Mode, Color Power Angio, CW/PW/Color Doppler
- Duplex CW Doppler, Cine Loop, User-Defined presets and calculations



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.