



auxoMEDICAL

Powerheart AED G3 Plus - Fully Refurbished

Designed for professional and medically trained responders, the G3 Pro features the industry's first color ECG display, manual override, and 3-lead ECG monitoring capabilities. With its RHYTHMx analysis algorithm, the G3 Pro lets first responders proactively monitor for life-threatening arrhythmias. This is a major benefit. The earlier a responder recognizes cardiac arrest, the shorter the time to defibrillation, and the better the chance for survival.

Features:

- Color ECG display, manual override, and 3-lead ECG monitoring capabilities.
- Our patented Rescue Ready technology to self-test all critical components (battery, pads, and electronic systems) daily.
- Non-committed shock, pacemaker detection, and synchronized shock features.
- A built-in metronome (which can be switched off) to set the pace for CPR compressions.
- Escalating STAR biphasic technology for customized patient therapy.
- Allowable Energy Range: Escalating Variable Energy (VE) 95J to 351J
- RescueCoach
- Metronome
- Daily, weekly, and monthly Rescue Ready self-testing
- Fully or semi-automatic mode
- Synchronized shock
- Pacemaker pulse detection
- Non-committed shock
- Download patient event data
- Pads tested for presence and functionality



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.