



Samsung Accuvix A30 Ultrasound - Fully Refurbished

Samsung's Accuvix A30 Ultrasound System establishes a benchmark for imaging. Offering features such as the proprietary EZ Exam and ElastoScan, this system facilitates clinical decision-making and boosts diagnostic confidence. Its 21.5-inch large screen LED monitor combines with enriched 3D performance to deliver uncompromising image quality for increased detection rates. Its customizable interface and progressive ergonomic design provide flexibility and ease of use. This is complemented by advanced automation and imaging tools which transform separate exams into a single streamlined procedure for faster, more relevant results.

Features:

- 2D, M-Mode, Color M-Mode
- Color/PW/Power/Spectral/Power Doppler
- 3D/4D imaging
- Directional Power Doppler Imaging
- 21.5" wide-screen High Definition LED flat-panel monitor 1,920 x 1080 pixels
- 9-inch LED touchscreen control
- Cardiac Measurements
- ECG
- SonoView image storage and retrieval
- EZ Scan
- DICOM
- Panoramic Imaging
- Elastoscans elastography
- Auto IMT automatic Intima-Media Thickness measurements
- HDVI high definition volume imaging
- Speckle Reduction Imaging
- Spatial filter omnidirectional imaging
- Advanced 3D facial detection and smart filter for volume imaging



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.