

Skytron Aurora II LED Lights - Fully Refurbished

Skytron's Aurora II is a LED surgical lighting system equipped with high-performance hybrid LED optics that provides surgeons with direct control of intensity and focus. Aurora II's features include manual focus of high intensity LEDs and a centrally positioned ring of fixed-focus LEDs, which can enhance intensity while minimizing eye fatigue. Color temperature is selectable for every procedure, including bright (4,500K) or soft white (4,000K). The system provides bright, cool and shadow-free illumination without heat or glare, even during extended surgical cases, allowing the entire surgical team to remain comfortable while under bright illumination.

Features:

- Direct Surgeon Control of Intensity & Focus
- Brighter Intensity, Deep View Illumination
- Cool Beam Temperature, Less Heat
- Shadow Free Optics
- Ultra Energy Efficient, Ultra Low maintenance
- Eco-friendly LEDs (last over a decade)
- Precise, 360° Positioning at all points of articulation, without drift
- Future Ready Architecture
- HD Camera Ready



Options:

- Single, dual, and triple configurations
- Extended arms for hybrid, imaging, etc that need an arm reach up to 10'
- Ergon Skybox Center and Center Tandem Mounts
- Additional display arms
- Mobile stand light model

Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.