



Skytron ST9 Exam Lights - Fully Refurbished

ST9 Exam Light provides high quality, high performance illumination within an extremely durable light weight polymer light head, creating the perfect balance between performance, maneuverability, and durability. The ST9 features advanced high intensity, shadow free optics without glare. Applications include LDR, LDRP, OB/GYN, Exam, Emergency Room & Special Procedures.

Features:

- High Intensity Illumination with up to 45,200 lux (4,200 fc) at a working distance of 42"
- Advanced VSRD (Vertically Segmented Reflector Design) Optics produce a myriad of rich and vibrant lighting segments blending into a high quality, soft white column of light
- Radically Durable, Featherlight Polymer Light head Construction provides for optimum maneuverability and durability
- Cool, Color Corrected Light provides natural Tissue rendering with a Color Rendering Index of 93
- Non-Sterile Positioning Handle is provided on the light head for positioning control
- Sealed Light head provides Smooth Surfaces for easy wipe-down between cases
- Fully Adjustable Friction Points assure smooth Movements and drift free performance
- Light head Mounted Intensity Control provides Infinite range of intensity adjustments
- Total system output is to be 45,200 lux at a color temperature of 4,000° K
- Configuration options include wall mounted radial arm, ceiling mounted radial arm, or portable stand



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.