



Steris AMSCO Examiner 10 Lights - Fully Refurbished

The AMSCO Examiner 10 Examination Lighting System illuminates with brilliant intensity for accurate diagnosis and treatment in physician offices, clinics, emergency rooms, nursing homes, and veterinarian offices.

Features:

- Virtually eliminates drift that can interrupt procedures
- 34,500 lux
- Produces a large light pattern
- Superior color rendition for accurate diagnosis and treatment
- Reduces glare and provides a cool beam of light
- Easy to position
- Variety of mounting options
- Intensity: 34,500 lux
- Pattern size: 10" (254 mm)
- Color temperature: 4400K
- Mobile base size, height, weight (all include 20' power cord)
 - 28 x 28" (711 x 711 mm), 55" (1397 mm), 41 lb.(18.6 kg)
 - 34 x 34" (864 x 864 mm), 55" (1397 mm), 51 lb.(23.2 kg)
 - 34 x 34" (864 x 864 mm), 67" (1702 mm), 55 lb.(25 kg)



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.