

Steris AMSCO Gemini Light - Fully Refurbished

Steris AMSCO Gemini Surgical Light allows full range of mechanical and optical flexibility required for surgery. Ideal for minor surgical procedures, obstetrics, and specialized examinations. The 22" (559 mm) diameter lighthead provides cool, shadow-reduced, color-corrected light at the operative or examination site. The touch switch control panel on the lighthead allows selection of small or large beam pattern; and low, medium, or high beam intensity. Ceiling mount, track mount, and mobile floor stand models are available.

Features:

- Allows full range of mechanical and optical flexibility required for surgery
- Easy, touch switch operation
- Well-defined lighting with no hot spots and excellent shadow control
- Provides cool, shadow-reduced, color-corrected light
- Low, medium, or high beam intensity
- Available in ceiling mount, track mount, or mobile floor stand
- Light Head: 22" (559 cm)
- Floor Stand Mount: 175 lbs (79 kg)
- Single Ceiling Mount: 165 lbs (75 kg)
- Dual Ceiling Mount: 245 lbs (111 kg)
- Drift free stability
- Excellent optical performance
- Easy touch switch operation
- Color temperature: 4200 K
- Voltage requirements: 100/120/220/240 volts 50/60hz
- Lamp voltage: 200w



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.