

## Steris AMSCO QDJ05 Warming Cabinet - Refurbished

---

The Steris AMSCO Single Compartment Warming Cabinet is the perfect size when space is limited and convenient access is required. This model provides 3.1 cu. ft. (18" depth) or 4.2 cu. ft. (24" depth) capacity and offers a compact size to fit any desired space. Designed to be more than a blanket warmer, the Single Compartment Warming Cabinet warms sterile intravenous fluids (IV fluids), surgical irrigation fluids, or linens and blankets.

### Features:

- LED readout for instant temperature confirmation
- Superior durable materials to withstand years of use
- Adjustable shelving for flexible organization and storage
- A design to meet AORN Recommended Practices for Safe Environment of Care, including technology to warm IV fluids
- Reliable heating with uniform and efficient heat transfer
- Visual access to inventory with dual-pane tempered glass or traditional solid steel doors
- Temperature lock-out for consistency and control
- Celsius and Fahrenheit selection for convenience and accuracy
- Single heating chamber with up to 4.2 cu. ft. total capacity
- Optional keypad door lock and electronic data reporting with USB port
- Free-standing or recessed installation
- Steel or dual-pane tempered glass doors
- 18" (457mm) or 24" (610 mm) depth



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.