



Steris V116 Steam Sterilizer - Fully Refurbished

Steris V116 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V116 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles.

Features:

- Florescent display for easy reading
- Ink on paper impact printer
- 16" x 16" x 26" chamber
- Rack and 2 shelves
- Freestanding unit
- Vertical Door
- Pre-vacuum Steam Sterilizer
- Fast & efficient
- Emergency manual exhaust valve
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Communication interface compatibility w/ most PC devices
- Electronic water conservation valve

Specifications:

- Chamber Size: 16" x 16" x 26"
- Overall Exterior Dimensions (including door lock-wheel): 74.5"H x 30"W x 35.75"D (with door closed)
- Operating Temperature: 250° to 270°
- One (1) Rack, Two (2) Shelves
- Touch Screen Display

Additional Options: Cabinet Package, Loading Cart, Integral Steam Generator, Water Filtration Options



Why an Auxo Medical refurbished sterilizer?

Auxo sterilizers are refurbished to the highest standards in the industry. We completely disassemble each sterilizer to its core and completely rework or replace all parts and return the autoclave to OEM standards functionally and cosmetically. Each autoclave is reassembled, tested and is warranted against defects.