



Crest Ultrasonic Digital Bench-top Cleaner- Fully Refurbished

Crest Ultrasonics' patented ceramically enhanced transducers provide greater cleaning power and reliability, when compared to conventional wafer-type transducers — and they have transferred that high performance technology to our Powersonic™ line of benchtop ultrasonic cleaners. The tabletop cleaners ensure uniform cleaning throughout the tank by sweeping the ultrasonic frequency 3 kHz, which creates overlapping ultrasonic waves. This approach eliminates inconsistent cleaning due to hot spots in the cleaning solution.

Features:

- Tank and wrap of Stainless Steel
- Internal tank dimensions of 19.5 x 11.75 x 8" deep
- Overall Dimensions - 21 x 12.75 x 14.5"
- 7 gallon tank capacity
- 38 lb. ship weight
- 300 average Wattage
- Frequency 42 kHz- 45 kHz
- Lifetime Heater Warranty
- 2 year Warranty on parts and Labor
- Main on/off switch
- Cover included with all units
- All models meet FCC requirements
- All models are ETL/CSA approved
- All models are available in 117 volts and 230 volts
- Warranty: 2 years
- Basket Sold Separately
- Deep-drawn stainless steel cleaning tank with rounded corners for easy cleaning
- Industrial bonded transducer assembly
- Intrinsically safe heater with dry-run protection
- Digital power control of the ultrasonic generator provides consistent cleaning results regardless of tank temperature



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.