



## DuraSonic Ultrasonic Cleaner Digital 5.25G- Fully Refurbished

The new benchmark in benchtop cleaning! DuraSonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price. The DuraSonic ensures uniform cleaning throughout the tank by sweeping the ultrasonic frequency 3 kHz, creating overlapping ultrasonic waves. This eliminates inconsistent cleaning due to hot spots in the cleaning solution, caused by the unvarying or "standing" waves produced by conventional ultrasonic power supplies.

### Features:

- Main on/off switch
- Three configurations available: Analog model with timer, Analog model with timer and adjustable thermostatically controlled heater, or Digital model with timer, heat, power control and degas
- All timers have continuous-operation capabilities
- Cover included with all units
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency: 40 kHz center, 42 - 45 kHz sweep
- Adjustable, thermostatically controlled heaters (ambient to 80°C)
- Precise control of time, heat, power and degas on digital units
- Drains with flow control valve on all 3/4 gallon and larger units
- Basket is sold separately
- Tray Size : 19.5"x11.75"x 6"
- Outside Dimensions (LxWxH) : 21"x12.75"x14.5"
- Sonic Power Watt Avg. : 240
- Meet FCC requirements and are ETL/CSA approved
- Available in 117 volts and 230 volts
- Warranty: 1 year parts and labor, lifetime heater warranty
- Timer: 0 - 99 Minutes



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.