



## Ferno 35A Ambulatory Stretcher- Fully Refurbished

The 35-A Mobile Transporter X-Frame Ambulance Cot provides exceptional value for x-frame users. The 35-A weighs only 74 lbs, and has a record of dependability with over 20 years of use by EMS services all over the world. A full range of options and accessories enhance patient care and help reduce situations that may lead to patient or operator injury.

### Features:

- X-Frame loading style.
- Fold down drop frame.
- Foot end lead handle.
- Easy-pull release handle.
- Side release handle.
- Lubrication free.
- Easy roll wheels.
- 8 bed positions.
- 2 loading heights.
- Backrest adjusts 0-75°.
- Gas assist backrest.
- One hand swing down sidearms.
- Safety U bar and hook.
- Tubular aluminum construction.
- Tested to industry standards.
- Built-in shock frame.
- 175 fastener compatible.
- Length Max: 79 in (200 cm).
- Length Min: 63 in (160 cm).
- Height Max: 40 in (101 cm).
- Height Min: 9 in (23 cm).
- Loading Height: 28/32 in (71/80 cm).
- Weight: 74 lbs (34 kg).
- Load Limit: 500 lbs (227 kg).
- Width: 24 in (61 cm).
- Warranty: 6 month



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.