



auxoMEDICAL

Gendron MRI 1000 Gurney- Fully Refurbished

The Gendron 1000MR Bariatric MRI Transport Stretcher is designed specifically for use in the MRI suite, making it ideal for hospitals and medical centers everywhere. The MRI Transport Stretcher is constructed of non-magnetic materials making it safe for use in the suite, however it should not be utilized directly inside of the MR system. Proudly made in the USA, the 1000MR MRI Stretcher features side rails that swing down for easy access, a telescopic IV holder and a manual fowler backrest and mattress. Ideal for transferring patients from one department to another, the 1000MR has a weight capacity of 350 lbs.

Features:

- Side rails that swing down for easy access
- Telescopic IV holder and
- Manual fowler backrest and mattress
- Ideal for transferring patients from one department to another
- Weight capacity of 350 lbs.
- Warranty: 1 year parts only
- Safe for use in and around the MR suite
- Constructed of non-magnetic materials
- High capacity for the needs of the obese patient
- Height: 31 Inch
- Top Dimensions: 28 W X 76 L Inch
- Four 8 Inch Casters, 2-Locking
- 2" hinged litter pad
- Six position fowler backrest



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.