



Hausted 5E8 power Eye stretcher- Fully Refurbished

The 5E8 is a multi-position stretcher intended for treatment, transport, or recovery of a patient. With the uniquely designed articulating head section, this unit can be used for specialty procedures which require precision head positioning, such as ophthalmology, plastic surgery, or ENT procedures.

Features:

- Articulating head piece with dual operating knobs
- Power hand held controls:
 - Height Adjustment
 - Trendelenburg/Reverse Trendelenburg - Backrest Adjustment
 - Knee Flex Adjustment
 - Auto Contouring Adjustment
- Nurse Control Panel at foot-end
- Emergency Backrest Release lever to override electric/battery
- Battery back-up allows power functions when not plugged in



Specifications:

- Trendelenburgs: 0-16°.
- Height Range:
 - High: 02" (81 cm).
 - Low: 22" (61 cm).
- Overall Width: 29.8" (75.8 cm).
- Patient Surface Width: 24.5" (62.2 cm).
- Overall Length: 88.5" (225 cm).
- Patient Surface Length: 81.5" (207 cm).
- Backrest Degree of Movement: 0-80°
- Knee Flex Degree of Movement: 0-60°
- Retracto Rail Height: 14" (35.5 cm).

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.