



## Hill-rom GPS 880 Series Stretcher- Fully Refurbished

The Hill-Rom GPS General Procedure Stretcher is ideal for transport, recovery and minor procedures. It Features hydraulic adjustable height, with controls located on both sides. The Side rails collapse and can be tucked away when not in use, allowing zero transfer gap. All four wheels brake from one pedal located at all four corners of the stretcher.

### Features:

- Overall length: 84"
- 8" casters
- Weight capacity: 500 lbs.
- Trendelenburg/Reverse Trendelenburg
- Dual-end brake and steer
- Sleep surface: 78" (H) × 26" (w)
- Height range: 23.5" to 28.5"
- Four IV rod locations
- One permanent IV rod
- Integrated storage bin
- Full perimeter bumper
- Stow-away side rails
- Radiolucent surface
- Wide variety of room accessories
- Oxygen bottle holder
- Pneumatic head section
- Pneumatic thigh gatch
- 5th Wheel steering
- Service contracts for preventative maintenance



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.