



## Lifepak 20E Defibrillator- Fully Refurbished

The LIFEPAK 20e defibrillator/monitor is intuitive, yet includes sophisticated features that make a real difference. For difficult-to-defibrillate patients, the 20e uses Physio-Control ADAPTIV™ biphasic technology for the maximum range of energy settings—up to 360 joules. For patients who need additional shocks, this has been shown to be a better strategy for terminating shockable heart rhythms in patients who need more than one shock.

### Features:

- Hospital Grade Power Cord
- 3-Lead ECG Cable
- Quik Combo Cable
- Spo2 Adult Reusable Finger Sensor
- Test Plug
- New Battery
- Set of Adult Pads
- Biocertified
- Warranty:1 year
- In-hospital first responders equipped with a LIFEPAK 20e defibrillator/monitor can make the lifesaving difference for victims of cardiac arrest.
- Highly intuitive to use and adapts to various patient environments.
- Easy for infrequent AED-trained responders to quickly understand and use.
- More advanced monitoring parameters such as ECG, external pacing and pulse oximetry are displayed clearly through color-matched waveforms and values.
- Capable of escalating energy up to 360J for the broadest therapeutic dosage for patients who fail to respond to previous shocks.
- Automated daily self-tests, and a battery gauge provide assurance your device is ready
- Two display options
- Seven main operating modes
- AC line operated device with an internal battery as backup.
- Can easily print a CODE SUMMARY™ report.
- Capable of transferring data records by IrDA version 1.0



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.