



Lifepak 500- Fully Refurbished

Owning a LIFEPAK 500 AED puts you in great company with the other tens of thousands of owners of this device. . Most EMS agencies in the United States rely on Physio-Control for pre-hospital defibrillation equipment. The features and configurability of the LIFEPAK 500 make it user-friendly. Its technology makes it a proven life-saver. Like all AEDs, this one talks you through the rescue. But that's just the start. It also has an LCD display so you can read the instructional prompts if background noise prevents you from hearing them.

Features:

- The LIFEPAK 500 AED is available with one of two defibrillation waveforms: monophasic or biphasic.
- The LIFEPAK 500 AED uses disposable QUIK-COMBOTM pacing/defibrillation/ECG electrodes
- Infant/Child Reduced Energy Defibrillation Electrodes
- Voice prompts
- Tones
- Flashing LEDs
- Screen messages
- Recertified by trained biomed technicians
- LCD screen that delivers text prompts along with voice prompts
- Escalating Bi-phasic energy
- Programmable CPR options
- Adult pads included
- You can change the timing of CPR interval
- You can program an instruction to perform CPR while the AED is charging
- Medical professionals can activate an asystole detector, too



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.