



MRI Fowler Back Stretcher- Fully Refurbished

Fowler positioning comfortably accommodates patients for prolonged periods. Ideal for patient transport to and within the MRI suite, the MRI Fowler Back Stretcher is non-magnetic, all aluminum, and designed for heavy-duty use. The Fowler back is infinitely adjustable and the stretcher comes complete with mattress, restraining straps and retractable guard rails. Five-inch ball bearing swivel casters provide easy, controlled mobility and lock firmly to secure the stretcher in position as needed.

Features:

- 500-lb patient capacity
- Non-magnetic, all-aluminum, heavy-duty use design
- Includes mattress and retractable guard rails
- Completely non-ferrous
- Wheels: 5" locking casters
- Finish: Powder coat
- Shipping Weight: 170 lb (77.1 kg)
- Warranty: 1 Year
- Fowler back positioning

Specs:

- Overall Dimensions: 74.5" L x 27.75" W x 34" H
- Tabletop Dimensions: 72" L x 24" W
- Fowler Back: 0° to 80° infinitely adjustable, manual
- Side Rail Height from Stretcher Surface: 7.9", from Floor: 41.7"
- Compliance: MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System
- Weight: 110 lb.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.