



Midmark 500 Baby Crib- Fully Refurbished

Midmark 500 pediatric stretcher/crib, with its unique oval design, rolls easily on 5 oversized locking casters. It features manual Fowler, optional storage tray and pad. One handed operation lowers and raises side rails with side Rail Intermediate Stop Position. The 500 Pediatric Crib Stretcher offers full-function capability and it's unique 2-post support system design allows 97% access to the patient.

Features:

- Steel caster provides increased maneuverability.
- 5" casters with central lock and/or steering with controls from 2 corners.
- Trendelenburg/Reverse Trendelenburg.
- 2 IV pole receptacles.
- White powder coat finish.
- Unique rounded corners.
- Rails with 2" bar spacing.
- Full function fowler.
- Litter Top: Solid pan.
- Latex-free.

Specs:

- Patient Surface: 26" x 56.25" (66 cm x 143 cm).
- Overall Length: 60" (152 cm).
- Overall Width: 30.25" (77 cm).
- Height:
 - Fixed Litter Top: 33.25" (84 cm).
 - Top of Crib: 59.5" (151 cm).
 - With Security Top: 76.75" (195 cm).
- Litter Positioning:
 - Fowler Angles: -5°, 0°, 20°, 30°.
 - Trendelenburg: 0°, 5°, 10°.
- Side Rails Intermediate Stop: 13.75" (35 cm).
- Casters: 5" (13 cm).



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.