



Midmark 540 Stretcher- Fully Refurbished

The full-function 540 Stretcher is both versatile and economical. It features a patient weight capacity of 350 pounds, full perimeter bumper, hydraulic height adjustment, 8-inch casters, and four corner IV pole positions. The 540 is also latex-free.

Features:

- 500-pound patient weight capacity
- Fowler Position: 0-90 degrees
- Pneumatic assist positioning
- Precision casters lock from any corner
- Trendelenburg: 12 degrees positive
- Reverse Trendelenburg -7 degrees
- Hands-free hydraulic height adjustment: 24" - 33"
- Patient Capacity: 350lbs
- Tuck away side rails
- Dimensions: 29.5in W x 83in L
- Wheel diameter: 8"
- Conductive wheels
- Central brakes
- Transfer gap: 0
- Patient Surface: 77" x 27"
- 4 Wheel Locking
- Steering
- O2 Tank Holder



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.