



Sharpertek 11 Liter Ultrasonic Cleaner- Fully Refurbished

Sharpertek has developed, engineered and priced their ultrasonic cleaners based on the requirements and demand of the consumer. This American based company manufactures industrial sized ultrasonic cleaners which will function in a large range of industries. There is no limit to what these products can clean

Features:

- Free Parts' Basket
- Industrial grade model
- Side Handles for easy handling.
- Heating element for best cleaning results
- Heating element can stay on all day even if ultrasonic power is off.
- Stamina XP™ Series
- Digital Ultrasonic Timer
- Industrial Grade Transducer for best reliability
- Ultrasonic Frequency: 120KHz
- Adjustable feet provide stabilization and prevent rocking
- Quickly degases in 5 minutes
- Tank Material: Stainless Steel

Specs:

- Tank Capacity: 11 Liters (approx. 3 Gallons)
- Tank Inside Dimensions: 12" x 9" x 6" (L x W x H)
- Unit Dimensions: 13" x 11" x 16.5" (L x W x H)
- Power Supply AC 100 ~ 120V, 50 / 60Hz
- Power 850 W
- Weight: 35 Lb
- Warranty: 2 year



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.