



Sonic Systems, Inc. 400 Series Dual Tank- Fully Refurbished

In the hospital and healthcare environment, surgical instruments and devices require high level cleaning prior to sterilization and re-use. The use of ultrasonic cleaners for this critical step has been the accepted practice for generations. The 400 series allows you to add the convenience of just cleaning or rinsing and drying in a single console. SSI engineers and manufactures the 400 series in-house to ensure a quality product and customer satisfaction.

Features:

- All stainless steel cabinet
- Model 422/HL
- Tank dimensions of 15" x 24" x 12" d
- Safety interlocked front doors
- Tray Elevator
- Automatic fill and drain
- Piezoelectric ULTRABRAZE[®] transducers vacuum brazed to the all 14 gauge 316L radiating tank bottom.
- 1000 Watts ultrasonic
- Remote generator housed within the console
- Frequency of 40 kHz \pm 2 kHz random sweep
- Low level safety control
- Wash tank heater
- Automatic detergent injector
- Spray rinse
- Forced hot air
- 5.2kW
- 235 cfm
- PLC Controlled
- Rear Control Panel with Operator Interface Terminal (OIT)
- Selectable Liquid Level
- Aqualift[®] Power Operated Lid



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.