



auxoMEDICAL

Stryker Atlas 660- Fully Refurbished

The Stryker Atlas 660 Stretcher is the gurney solution for all departments. The 660 pound weight capacity with an extra-wide 31" patient top meets the needs of all patients. The Stryker Big Wheel technology is proven to offer easier mobility and control. It requires 50% less maximum startup force and 60% less turning effort. The Big Wheel cuts down on the possibility of back injuries to the operators. The 4" thick mattress provides greater support and comfort for patients. This state-of-the-art Atlas 660 Gurney is sold refurbished and includes pads.

Features

- 660-pound weight capacity
- Big Wheel
- Four Corner IV Receptacles
- Glideaway Siderails
- Rubber Bumpers
- Pop-up Push Handles
- Three-sided foot control hydraulics
- Hands-free Trendelenburg
- Integrated patient transfer boards
- 8" Casters with wheel covers
- Unilower Pedal



Specifications

- Overall Length: 84"
- Overall Width with Siderails Up: 34.5"
- Overall Width with Siderails Down: 31"
- Weight Capacity: 660 lbs
- Height Range: 22" - 33"
- Patient Surface: 31" x 76"
- Transfer System: 8" x 42"
- Knee gatch: 0° - 30°
- Side rails: 14 in (36 cm) x 54 in (137 cm)
- Trend/Reverse Trend: $\pm 15^\circ$

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.