



Stryker Trio Surgical Stretcher- Fully Refurbished

The Trio is a fusion of surgical stretcher, mobile operating table, and recovery platform in one unit. With the Stryker Trio the patient remains on the same surface for pre-op, surgery, and post-op. By doing so this greatly reduces patient handling and improves patient and staff safety. All while increasing case efficiency turnover time.

Features:

- 3" (8cm) Ultra comfort mattress with sealed seams
- Quadra Channel™ stabilizing lift system
- Siderails enable fluoro clearance
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activator
- Ultraview litter top for unobstructed fluoroscopic access
- Three-sided foot controls for lift/lower functions
- Retractable fifth-wheel steering
- Low profile base with thermoformed hood
- Full-length surgery accessory rail



Specs:

- Height (min): 26 (66 cm)
- Height (max): 40" (102 cm)
- Overall length: 80" (102 cm)
- Overall Width (siderails up): 33" (84 cm), (siderails down): 30" (76 cm)
- Weight Capacity: 500lbs (228kg)
- Backrest: 1° - 80°
- Foot Section: 0° - 50°
- Trend/Reverse Trend: +20° / -15° Trend
- Lateral Tilt: ±15°
- Patient Surface: 27" × 78" (69 × 178 cm)
- Siderails: 12" × 49" (30 × 124 cm)
- 6" (15cm) Omni Surface™ casters

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.