



Zoll M Series 3 Lead ECG, SpO2, NIBP- Fully Refurbished

ZOLL M Series® defibrillators offer a family of products to meet your clinical needs from BLS to ACLS. Products range from AED and advisory defibrillators to manual units and Critical CareTransport (CCT). Thanks to the consistency built into each and every unit, you train once and use everywhere. Standardization has never been so easy.

Features:

- The smallest, lightest, brightest defibrillator unit in its class
- Battery Management built-in AC power and battery charger
- The ZOLL Uniform Operating System
- ZOLL's proprietary Rectilinear Biphasic low-energy waveform.
- Non-invasive Temporary Pacing offers superior capture at lower mean thresholds
- Built-In Recorder
- High resolution display
- Includes soft pads, EKG cable
- Optional Hand Paddles.
- 3 lead, SpO2, NIBP, Biphasic
- Available using synchronized cardioversion, external pacing, ECG with advanced monitor



Specs

- Size: 6.8 in (17.3 cm) high x 10.3 in (26.2 cm) wide x 8.2 in (20.8 cm) deep.
- Weight: 11.5 lb. (5.23 kg) with Multi-Function Cable and battery; 13.5 lb. (6.14 kg) with paddles.
- Design Standards: Meets or exceeds UL 2601, AAMI DF-39, AAMI DF-2, and IEC 601-2-4.
- Waveform: Damped sinusoid or ZOLL Rectilinear Biphasic.
- Charge Time: Less than 7 seconds with a new fully charged battery (first 15 charges to max energy.) Depleted batteries will result in a longer defibrillator charge time.

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.