



auxoMEDICAL

Tuttnauer® 3870EA Autoclave - Fully Refurbished

With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer autoclaves fill, sterilize, exhaust and dry at the touch of a button.

Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- Hepa filter 0.2 µm air filter (replaceable).
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs
- Program and customize cycle parameters to suit your needs and maximize flexibility
- Integrated printer for cycle documentation (optional)
- Interdependent temperature and pressure monitoring
- Cycle information recovery in the case of power failure or cycle interruption
- Insufficient water detector with automatic shut-off
- Supplied with a 3 shelf stand (30"Wx36"Dx34"H)



Specifications:

- Chamber Dimensions: 15" x 30"
- Volume: 22gal/85L
- Overall Dimensions: 34.5" x 26" x 20.6
- Standard Cassettes: 15 Full
- Tray Dimensions: 26.6" x 11.3" x 1" and 26.6" x 13.8" x 1"
- No. of Trays: 2
- Standard Unwrapped Cycle Time: Cold 32 Min, Hot 23 Min
- Voltage/Frequency: 230V, 50/60Hz
- Power/Current: 3000W/13.6A
- Shipping Weight: 213lbs

Why an Auxo Medical refurbished autoclave?

Auxo autoclaves are refurbished to the highest standards in the industry. We completely disassemble each autoclave to its core and completely rework or replace all parts and return the autoclave to OEM standards functionally and cosmetically. Each autoclave is reassembled, tested and is warranted against defects.