



Hill-rom P8000 Procedure Stretcher - Fully Refurbished

The Hill-Rom P8000 Stretcher is an ideal choice for the ER, OR, PACU or Surgery Center. Hill-Rom gurneys have superior control around corners and can self-adjust over uneven floor surfaces. These engineering features prevent vibrations and discomfort for the patient. The reinforced frame is made of auto-grade steel for durability and quality. The refurbished Hill-Rom stretcher has a 700 lbs weight capacity, making it ideal for providing top care to the ever-changing needs of the community. The four-wheel brake system is easy to operate, and the ergonomic push handles provide an easy-to-use gurney for the hospital staff. The Hill-Rom P8000 Stretcher is sold refurbished and includes mattress.

Features:

- 700 lb weight capacity
- 26" and 30" width
- BackSaver Fowler (includes manual knee gatch)
- Brake/Steer pedals – head and foot
- Next generation Steering Plus™ system
- Integrated O2 tank holder with utility tray
- Siderail system (tuck-away)
- New Mattress

Options/Accessories:

- Manual knee catch
- Active hand brake
- Enhanced comfort mattresses
- IV poles
- Restraint straps



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.