



## Hill-Rom P8000 Stretcher - Fully Refurbished

With its standardized platform, Hill-Rom's P8000 stretcher is ideal for use in virtually any transport situation. Whether you are in triage, radiology, or attending to critical patients this standardized platform provides the functionality needed for trauma situations with the mobility required for transporting. A full-length x-ray cassette access provides ultimate imaging flexibility from head to foot and captures quality images right on this stretcher, without moving your patient.

### Features:

- Maximum weight limit: 500 lbs. (317.5 kg)
- Low position -- floor to top of surface: 24.5" (61.6 cm)
- High position -- floor to top of surface: 37.8" (95.9cm)
- Overall length: 83" (210.8 cm)
- Overall width (siderails up): 36" (91.4 cm)
- Overall width (siderails stored): 30.3" (76.8 cm)
- Siderail length: 47" (119.4 cm)
- Siderail height above mattress/ surface: 14" (36.8 cm)
- Mattress size: 26" (66 cm) x 75" (190.5 cm)
- Maximum head elevation: 90°
- Maximum Trend/Reverse Trend: 18°
- Floor to base clearance: 3.5" (8.9 cm)
- Caster size: 8" (20.3 cm)



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.