



auxoMEDICAL

**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

-

2020



Autoclaves

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.



ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: support@auxomedical.com

Table of Contents

Autoclaves..... 1-8



auxoMEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

New & Refurbished
Medical Equipment

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

Amsco Eagle 10 Plus

Amsco Eagle Ten Autoclave provides clinics with Safe, fast and dependable sterilizer. The Amsco Eagle Ten Autoclave is a large capacity, counter-top sterilizer perfect for any size practice. Amsco has been a world leader in sterilization equipment, and the Eagle Ten offers many sophisticated features. The spacious 10" x 16" chamber can hold two utensil trays per load for wrapped, unwrapped and liquid goods. The outside dimensions are compact, saving counter space in the work area. Simply select wrapped, unwrapped, packs or liquid sterilization with the touch of one button. Other features include optional programmable controls, self-diagnostic error codes and safe operation



Features:

- Safe, fast and dependable sterilization.
- Large capacity, compact design
- Microprocessor controls
- Self-diagnostic error codes
- Safe operation
- Fast sterilization cycles
- Pre-programmed or programmable cycle selection.

Condition	Refurbished
-----------	-------------

Cox 6000 Rapid Dry Heat

The COX Rapid Heat Transfer Sterilizer is a fast, simple solution to your sterilization needs. Because it utilizes heated air that moves at a rate of 2,500 feet/minute, it is able to sterilize objects in as little as six (6) minutes -- faster and easier than any other sterilizer available.



No longer will you and your staff be bothered with all the problems inherent with other forms of sterilizations. Since the dry heat system won't corrode or rust your instruments, you won't have to replace them so often -- thus saving you money.

Sterilization Cycles:

- Cycle 1: 6 minutes for unwrapped instruments
- Cycle 2: 8 minutes for hand pieces and medical drills
- Cycle 3: 12 minutes for wrapped instruments

Included with Sterilizer:

- One Drawer
- One Removable Instrument Basket
- Basket Removal Tool
- Cooling Rack
- Instrument Rack (for wrapped instruments)

Condition	Refurbished
-----------	-------------

Dri-Clave 3000 Dry Heat Sterilizer 3 Tray

The FDA approved technique of dry heat sterilization works without pressure, steam, or chemicals. Dry heat sterilizes in the shortest possible time and you never need to dry the instruments. Ideal for all medical and dental instruments, it sterilizes instruments, needles, syringes, glassware, etc. at 320F (160C). Sharp instruments remain sharp; no steam or moisture to corrode the finest edge. Both operate on 115V AC and are UL Listed and FDA Registered.



Features

- The best method for sharp, micro, delicate instruments: Dry heat sterilizer kills bacteria without moisture and will not rust or dull sharp or serrated instruments.
- Maintenance-Free: Unlike wet heat methods, dry heat sterilizers require no routine cleaning.
- Safe and simple to operate: Dry heat sterilization works without pressure, steam or chemicals, which increases safety and makes operating simple.
- Convenient: Unit has automatic timers with heat "hold" capability for continuous use in a busy practice. Utility trays can be sterilized with lids on and stored without contamination. Plus, use it to conveniently dry wet packs.
- High-Quality: Made of stainless steel - it will not only never rust - it cleans up in a jiffy - and will look great for years.
- Covered by a full one year warranty.
- Economical: The EC-3000 Dri-Clave not only saves money with a low purchase price but continue to save money with low operating costs - no distilled water to buy!

Condition	New
-----------	-----

Dri-Clave Dry Heat Sterilizer - 2 Tray

The FDA approved technique of dry heat sterilization works without pressure, steam, or chemicals. Dry heat sterilizes in the shortest possible time and you never need to dry the instruments. Ideal for all medical and dental instruments, it sterilizes instruments, needles, syringes, glassware, etc. at 320 °F (160 °C). Sharp instruments remain sharp; no steam or moisture to corrode the finest edge. The operate on 115V AC and are UL Listed and FDA Registered. The best method for sharp, micro, delicate instruments: Dry heat sterilizer kills bacteria without moisture and will not rust or dull sharp or serrated instruments.



- Maintenance-Free: Unlike wet heat methods, dry heat sterilizers require no routine cleaning.
- Safe and simple to operate: Dry heat sterilization works without pressure, steam or chemicals, which increases safety and makes operating simple.
- Convenient: Unit has automatic timers with heat "hold" capability for continuous use in a busy practice.
 - Utility trays can be sterilized with lids on and stored without contamination
 - Plus, use it to conveniently dry wet packs.
- High-Quality: Made of stainless steel - it will not only never rust - it cleans up in a jiffy - and will look great for years.
- Economical: Our EC-2000 Dri-Clave not only save money with a low purchase price but continue to save money with low operating costs - no distilled water to buy!

Condition	Refurbished
-----------	-------------

Midmark Ritter M11

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.



Features include:

- Exclusive automatic-opening door for ease of use and a quick drying time
- Large 11" by 18" chamber for increased capacity
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration
- Self Program Controls enable you to create different cycle parameters to meet specific instrument processing needs
- Pre-Programmed Controls for convenience
- Green LCD Display for easy view of indicators

Condition	Refurbished
-----------	-------------

Midmark Ritter M9

The Ritter M9 is an autoclave that couldn't be easier and is time tested for reliability. You load it set it and forget it. Select your cycle and the time is set for you. Press Start and the Midmark M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam to give you fast drying time for your instruments. Its real convenience!



Features:

- Fills with correct amount of water when you hit the start button.
- Door opens automatically to dissipate steam and give you the fastest drying time.
- Built-in safety features such as overheat protection, a pressure relief valve and flashing "Error" light.
- Self Program Controls enable you to create different cycle parameters to meet specific instrument processing needs
- Green LCD Display for easy view of indicators

Condition	Refurbished
-----------	-------------

New Style Midmark M9/M11 Printer and cable

New Style Midmark M9/M11 Printer and cable Cable and 1 roll of recorder paper included



Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Pelton & Crane Magnaclave

The Pelton & Crane Magnaclave is a fully self-contained steam sterilizer with a large chamber, making it ideal for busy group practice, emergency clinics, surgery centers and small hospitals. The chamber capacity allows it to hold up to 14 Hu-Friedy cassettes or equivalent. The unit is sold with on basket tray and one medium shallow tray. It is equipped with a safety valve which limits the maximum chamber pressure to 38 PSI. This unit is sold refurbished and includes a stand. Electrical Requirements: Dedicated 30 amp circuit/220-240 VAC or 208 VAC.



- Large chamber capacity (16" dia x 30" deep)
- Total Dimensions:28"W, 34.5"D, 23" Tall (53" Tall w/ Stand)
- Safety lock system protects personnel from accidentally opening the door under pressure.
- Completes a flash cycle in 18 minutes from a warm start, 35 minutes from a cold start. Drying phase not included
- Reservoir drain is conveniently located on the front for easy cleaning.
- Simple and easy door locking mechanism.
- Water reservoir has large capacity to adequately cool condensed steam, reducing the need for refilling.
- Water level indicator
- Preheat feature allows faster sterilization time.
- Optional Stand available

Condition	Refurbished
-----------	-------------

Pelton & Crane Validator 10

The Pelton & Crane Validator 10 Autoclave is a economical, light-capacity sterilizing unit designed to accommodate the needs of small medical, dental, and other practices requiring reliable and safe access to sterilization equipment. Incorporated are a variety of safety and efficiency features, such as a door lock, which prevents the door from being opened while the chamber is pressurized. It also features a set of microprocessor-controlled features to ensure effective and safe operation; a vent valve will automatically open if chamber pressure exceeds specification with alarm and secondary release. Similar functionality exists for heat and electrical levels, as well as a sensor to prevent chamber overfilling. Unique "Pressure Locking" Door Design. Provides maximum safety.



- Automatic shut off and automatic venting system.
- Temperature and pressure gauges. Allow you to always be in control of the status and operation of the sterilizer.
- Safety Valve - Prevents excessive pressure build-up in chamber.
- Common parts can be easily replaced.
- You may insert the instruments wrapped, or unwrapped

Condition	Refurbished
-----------	-------------

Pelton and Crane OCR

The Pelton & Crane OCR Autoclave, Tabletop Sterilizer is known by some as "old faithful". These units were manufactured for years by Pelton and Crane and known for their solid reliability. When they are properly maintained, they will last forever.



The Sterilizer has a 9? x 18? chamber. This unit is ideally suited for surgery and veterinary practices and other specialty practices.

The Pelton & Crane OCR Autoclave is ideally suited for oral surgery practice groups and other specialty practices. It features a stainless steel door that is secured at four points and incorporates a door interlock switch to ensure that the door is properly locked before the cycle begins.

This is a work horse of a sterilizer. It will work forever for you. Pelton & Crane OCR Sterilizer with 9? x 18? chamber. Replacement parts are readily available from both the manufacturer and other sources. It is ideally suited for oral surgery practices groups and other specialty practices. This sterilizer comes with auto shut off, and auto venting mechanism.

Though this steam pressure autoclave is typically used in oral surgery practice groups, it can also be found in a range of other settings, including physician's offices, clinics, hospitals, nursing homes, and laboratories.

Reaches 250 °F in approximately twelve minutes from a cold start or approximately 5 minutes from a warm start. Maybe used as a steam or dry-heat sterilizer. Optional dry-heat door is available. Equipped with ASM E-type safety valve to limit maximum chamber pressure to 35 PSI. Reservoir drain conveniently located at front.

Condition	Refurbished
-----------	-------------

Ritter M11 Autoclave, M11-001

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.



Features include:

- Exclusive automatic-opening door for ease of use and a quick drying time
- Large 11" by 18" chamber for increased capacity
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration
- Pre-Programmed Controls for convenience
- Red LCD Display for easy view of indicators

Condition	Refurbished
-----------	-------------

Ritter Midmark M9, M11, Old Style, Connection Cable Included

Printer, Midmark M9, M11, Old Style, Connection Cable Included



Condition	Refurbished
-----------	-------------

Ritter Midmark M9, M9-001

The M9 autoclave couldn't be easier and is time tested for reliability. You simply load it - set it - and forget it. Select your cycle and the time is set for you. Press "Start" and the M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you fast drying time for your instruments.



Refurbished M9 Autoclave Features:

- Fills with correct amount of water when you hit the start button.
- Door opens automatically to dissipate steam and give you the fastest drying time.
- Built-in safety features such as overheat protection, a pressure relief valve and flashing "Error" light.
- Pre-Programmed Controls for convenience

Condition	Refurbished
-----------	-------------

Scican Statim 5000 G4 Cassette

The fully automatic, Statim 5000 G4 utilizes four, single touch sterilization programs to optimize instrument cleaning, from start to finish.

The Statim 5000 G4 has an additional sterilization cycle and a roomier cassette in order to process larger instruments. It's reliable and compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures. In just 17.5 minutes the Statim 5000 G4 can complete a wrapped cycle and only requires 9 minutes to complete an unwrapped sterilization cycle. The Statim 5000 G4 series has been modernized with a new color touch screen and intuitive controls.



Built-in data logger. No more messy thermal paper, clumsy printers or possible failures with external data loggers. The Statim 5000 G4 series now saves all the data from each cycle. The data can be copied to the included USB drive, viewed on screen, or stored in the Cloud via email.

- A large 3.5" high-resolution touchscreen offers a vivid display of messages and current cycle information all with extraordinary clarity.
- Ultra-fast 6 minute unwrapped cycle
- 14 minute wrapped cycle
- Exceptional instrument turnover
- 10 minute Dri-Tecc Drying Technology
- Aseptic transfer of instrument directly to the point-of-use
- Reduced instrument investment
- Prolongs instrument life

Condition	Refurbished
-----------	-------------

Statim 2000 Cassette Sterilizer

The Statim 2000 is a reliable cassette sterilizer that's compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures - it can complete an unwrapped sterilization cycle in just six minutes and a wrapped cycle in 14 minutes.*



The Statim 2000 is fully automatic. Simply select one of the three sterilization programs and with a single touch of the keypad, begin the cycle. An advanced microprocessor controls all functions and maintains optimal sterilization conditions, from start to finish.

Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Statim 2000 G4 Cassette Sterilizer

The Statim 2000 G4 is the same renowned and trusted autoclave it has been for over 20 years but now boasts a new contemporary look and connectivity that is the first of its kind. The G4 technology will change the way you interact by providing a direct channel of communication through the Internet to you, or anyone you desire. Still powered by SciCan's signature steam technology to provide sterilization and dryness at speeds faster than conventional chambered autoclaves, the Statim has been drastically upgraded with a level of interactivity never seen before.



Condition	Refurbished
-----------	-------------

Statim 5000 Sterilizer

The Statim 5000 sterilizer has an additional sterilization cycle and a roomier cassette designed to process bulkier instruments and greater loads. It's a reliable cassette sterilizer that's compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures - it can complete an unwrapped sterilization cycle in just nine minutes and a wrapped cycle in just 17.5 minutes.*



The Statim 5000 sterilizer is fully automatic. Simply select one of four sterilization programs and with a single touch of the keypad, begin the cycle. An advanced microprocessor controls all functions and maintains optimal sterilization conditions, at the start to the finish.

SciCan's Statim cassette sterilizers have a biological effectiveness that has been proven in tests for hollow and solid instruments at many internationally recognized institutes in Canada, the United States and Europe.

Condition	Refurbished
-----------	-------------

Tuttnauer 1730 Autoclave

The Valueklave 1730 will satisfy all of your basic sterilization needs at an economical price. The compact size (7" chamber) makes it a perfect primary unit for small offices or a second autoclave for larger practices. The fast cycle time and heated drying delivers a quick turnaround for wrapped or unwrapped instruments. The simple controls and easy closing door make this unit a breeze to operate. High speed and reliability are built into every Valueklave 1730.



FEATURES:

- Fast Cycle Time. 11 minutes from a hot start* (Standard Unwrapped Cycle Time).
- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications

Condition	New or Refurbished
-----------	--------------------

Tuttnauer 2340EA EZ9 Autoclave

This table-top autoclave is designed for the sterilization of wrapped and unwrapped instruments, glassware and other related items found in dental, medical and veterinary clinics, tattoo parlors, first aid rooms, hospitals, laboratories and other related industries.



The enhanced closed door active drying system maintains the sterile field and ensures complete drying of packs and pouches. The efficient air pump forces hot air circulation after the sterilization cycle for faster, more efficient closed door drying. Features a bacteriological Hepa filter 0.2 µm air filter (replaceable), Four pre-loaded cycles for unwrapped instruments, wrapped/packs and porous loads, glassware or slow exhaust for liquids or additional drying cycle.

All program parameters can be customized and stored. This unit has a contemporary design with a bright blue easy to read digital control panel. An automatic shut off at the end of both the sterilization and dry cycles. There are audible and visual cycle interruption alerts, and a double safety locking device that prevents the door from opening while the chamber is pressurized.

- Features/Description:
- Automatic shutoff at the end of sterilization and dry cycles.
- Automatic Dry Cycle at end of Sterilizer Cycle
- Front panel controls.
- Electronic Digital Control
- Safety device to prevent premature opening while unit is pressurized.
- Front side drain valve for ease of releasing reservoir water.
- Dual safety thermostat to protect against overheating.

Condition	New or Refurbished
-----------	--------------------

Tuttnauer 2340M Manual Autoclave

The Tuttnauer manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.



Features:

- Automatic shutoff at the end of sterilization and dry cycles.
- Front panel controls.
- Safety device to prevent premature opening while unit is pressurized.
- Front side drain valve for ease of releasing reservoir water.
- Dual safety thermostat to protect against overheating.

Condition	Refurbished
-----------	-------------

Tuttnauer 2540M

In the digital age, most autoclave manufacturers have discontinued production of manual autoclaves. Thank goodness Tuttnauer did not follow this trend. There remains a significant number of clinics, veterinary practices and medical offices who are not at all interested in changing their autoclave. The simple design of the Tuttnauer 2540M makes it an ideal unit for many satisfied users.



The Tuttnauer manual autoclave comes with an automatic shut off at the end of the cycle and a dual safety thermostat that protects against overheating. This affordable and practical autoclave is a great investment for anyone looking for a workhorse sterilizer to run day after day with very little maintenance or down time. As the Tuttnauer manual autoclave is one of our best sellers, we always strive to keep one in stock offering same day shipping to anyone who places an order.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- Long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications

Condition	New or Refurbished
-----------	--------------------

Tuttnauer 2540MK

The Tuttnauer M Series is also available in a quick cycle model, the Kwiklave (MK). This sterilizer offers all the proven, dependable features of our manual autoclave, with the added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turn-around or a high volume office. Speed and reliability are accomplished without sacrificing load size. 1 year parts and labor warranty.



Model 2540MK requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5kVA) may be required.

Features

- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.

Condition	New or Refurbished
-----------	--------------------

Tuttnauer 3870E Autoclave

With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer 3870E autoclaves fill, sterilize, exhaust and dry at the touch of a button.



Features of Tuttnauer 3870E

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- Hepa filter 0.2 µm air filter (replaceable).
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs
- Program and customize cycle parameters to suit your needs and maximize flexibility
- Integrated printer for cycle documentation (optional)
- Interdependent temperature and pressure monitoring
- Cycle information recovery in the case of power failure or cycle interruption
- Insufficient water detector with automatic shut-off
- Printer Included

Condition	New
-----------	-----

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Tuttnauer 3870E

The Tuttnauer Model 3870E autoclave is designed for sterilizing laboratory and medical equipment, media and even liquids. The Tuttnauer autoclave provides fast, safe and dependable sterilization without any need for special plumbing or venting. Tuttnauer autoclaves are manufactured with 316L stainless steel sterilization chambers which are highly corrosion-resistant. These autoclaves feature a fast exhaust option which reduces cycle time for solid goods and a slow exhaust option which releases chamber pressure slowly and prevents boiling over of liquids. Cold air pockets are eliminated by continuous pressure purging which ensures uniform sterilization.

The autoclave includes a 3.25L water reservoir which allows complete control the purity of the sterilizing steam and allows repeated use without the need to refill very often. The integral drying cycle prevents corrosion of expensive instruments and also dries absorbent materials. The autoclave also includes a low water level alarm which prevents damage to the heating elements. Tuttnauer autoclaves feature a double-locking door that is insulated for high heat and protects lab personnel from accidental burns. The door cannot be opened when chamber is pressurized by means of a dead-bolt mechanism for added safety.

- Automatic, Digital Control; Chamber Diameter x Depth: 15 x 27-Inch; Capacity: 84-Liters; Includes support stand and 2 Stainless Steel Tray Shelves (14 x 26.1-Inch)
- Includes: 60-Minute Timer, 3.25-Liter Reservoir & Low Water-Level Alarm. Temperature Range: 105 to 137°C (221 to 279°F)
- Digital display of temperature & pressure. Heating time to 273 °F (134 °C): 30-Minutes
- Heat-insulated, double-locking door protects against accidental burns. Door cannot be opened while chamber is under pressure.
- Automatic low water level shut-off prevents damage to heating elements when there is insufficient water for the sterilization cycle. UL Listed, CSA Certified & ASME Stamped.

Condition	Refurbished
-----------	-------------



Tuttnauer 3870EA Autoclave

The Tuttnauer 3870EA autoclaves have a closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches. A high efficiency air pump circulates hot air after the sterilization cycle for fast drying. These steam sterilizers use a replaceable HEPA 0.2m air filter.

Features of Tuttnauer 3870EA

- 15 x 30" Chamber.
- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- Hepa filter 0.2 µm air filter (replaceable).
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs
- Program and customize cycle parameters to suit your needs and maximize flexibility
- Integrated printer for cycle documentation (optional)
- Interdependent temperature and pressure monitoring
- Cycle information recovery in the case of power failure or cycle interruption
- Insufficient water detector with automatic shut-off

Condition	New or Refurbished
-----------	--------------------



Tuttnauer 3870M

This is Tuttnauer's largest manual table-top autoclave. Believe it or not there is still a demand for the 3870m analog unit among many healthcare facilities and veterinary practices. This beefy, no frills autoclave covers all the basics needed for proper sterilization. It can always be counted on even under harsh conditions. The large 15" x 30" chamber accommodates longer instruments or a large number of smaller instruments, making it a perfect match for spay and nueter clinics, equine practices, busy animal hospitals and multi-practice medical offices.



This freestanding autoclave is easy to set up requiring a dedicated 230 Volt, 13.6 Amp power outlet. No special plumbing is needed as the autoclave reservoir is filled with distilled water that can be purchased at your local supermarket. It is only necessary to drain the autoclave when doing the routine chamberbrite cleaning procedure. The waste water is caught in a bucket or the sink.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating

Condition	New or Refurbished
-----------	--------------------

Tuttnauer EZ10

The Tuttnauer EZ10 Fully Automatic Medical Autoclave Sterilizer features a simple one-touch design and pre-programmed cycles for: unwrapped instruments, wrapped and packs, glassware and an additional drying cycle. The Tuttnauer EZ10 Fully Automatic Autoclaves fill, sterilize, exhaust and dry at the touch of a button. Shut off at the end of both the sterilization and drying cycles is automatic. In addition, all program parameters can be customized and stored. The Tuttnauer EZ10 also features a replaceable HEPA filter air filter, and a high efficiency air pump that circulates hot air after the sterilization cycle for fast drying.



The Tuttnauer EZ10 is the easiest and most hassle free autoclave choice from QuickMedical. With the capability of sterilizing 3 full trays and 3 half cassettes at once, this is an excellent choice for the small practice desiring a one touch design in their sterilization routine. The choice of an integrated printer allows for print-out reports of cycle information.

The Tuttnauer EZ10 continuously monitors temperature and pressure for sterilization and has a double safety locking device preventing the door from opening while the chamber is pressurized. Depending on the instrument to be autoclaved the normal temperature ranges are 250°F (121°C) or 273°F (134°C). See More Info tab for detailed sterilization temperature information.

- Closed door active drying system will maintain sterility and ensure efficient drying.
- Includes a replaceable 0.2µm Hepa air filter.
- Front panel programmable controls.
- Safety door automatically seals when pressurized and cannot open until depressurized.
- 6 Gallons (23 Liters)
- Holds Up to 4 trays
- Cycle time: Cold start: 30 minutes, Warm start: 14 minutes

Condition	New or Refurbished
-----------	--------------------

Tuttnauer EZ9 Plus

The newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing.



Innovative EZView Multi-Color Display offers real time cycle information that is easy to read, even from a distance. The EZView displays cycle parameters, temperature, and pressure in the chamber. During a cycle each stage is displayed with an active color wheel and text.

Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high efficiency air pump circulates filtered air through the chamber for fast drying. 0.2µm HEPA filter (replaceable).

- Comprehensive 2 year parts & labor warranty and 10 year chamber warranty.
- 9" diameter X 19.8" deep chamber and flexible loading options. The exceptionally deep 19.8" chamber works in a variety of loading options to suit any office:
 - Use the 3 full sized trays for packs/pouches
 - Use the tray rack (included) to accommodate 2 full plus 2 half sized cassettes...no additional cassette rack to buy!
 - Use the optional Tuttnauer pouch rack (model AR910) to encourage proper loading
- Front or top water fill options: accessible EZFill front or top water fill option to enable easier filling of the reservoir.
- USB Port enables users to easily transfer cycle data via a flash drive to a PC.
- Network Port allows for connection to a local network to provide remote monitoring of the system when used with optional RPCR Software. Network port is located in the rear of the unit
- Dynamic Air Removal Technology to effectively remove air from the chamber.
- Designed to meet the most current sterilization standards ANSI/AAMI ST55.
- Optional printer documents a detailed history of each cycle.

Condition	New
-----------	-----

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers
Washer Disinfectors
Surgical & Exam Lighting
Surgical & Exam Tables Stretchers
& Patient Transportation
Vital Sign & Patient Monitors
Blanket/Fluid Warmers
Scrub Sinks & Accessories
Defibrillators
Electrosurgical Units
Surgical Endoscopy
Imaging
Tourniquet Systems
Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product	Page
Amsco Eagle 10 Plus.....	1
Cox 6000 Rapid Dry Heat.....	1
Dri-Clave 3000 Dry Heat Sterilizer 3 Tray.....	1
Dri-Clave Dry Heat Sterilizer - 2 Tray.....	2
Midmark Ritter M11.....	2
Midmark Ritter M9.....	2
New Style Midmark M9/M11 Printer and cable.....	2
Pelton & Crane Magnaclave.....	3
Pelton & Crane Validator 10.....	3
Pelton and Crane OCR.....	3
Ritter M11 Autoclave, M11-001.....	3
Ritter Midmark M9, M11, Old Style, Connection Cable Included.....	4
Ritter Midmark M9, M9-001.....	4
Scican Statim 5000 G4 Cassette.....	4
Statim 2000 Cassette Sterilizer.....	4
Statim 2000 G4 Cassette Sterilizer.....	5
Statim 5000 Sterilizer.....	5
Tuttnauer 1730 Autoclave.....	5
Tuttnauer 2340EA EZ9 Autoclave.....	5
Tuttnauer 2340M Manual Autoclave.....	6
Tuttnauer 2540M.....	6
Tuttnauer 2540MK.....	6
Tuttnauer 3870E.....	7
Tuttnauer 3870E Autoclave.....	6
Tuttnauer 3870EA Autoclave.....	7
Tuttnauer 3870M.....	7
Tuttnauer EZ10.....	8
Tuttnauer EZ9 Plus.....	8

Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER

Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's
Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions

Surgery Center Services of America
Chesapeake Healthcare Planning
Sodexo USA



STERILIZATION/BIO MEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's
Daughters
Virginia Eye Institute

OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTROSTEAM
SUSSMAN

Washer Disinfectant

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPERTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBROPTICS
WELCH ALLYN

EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424
www.auxomedical.com
sales@auxomedical.com

 **auxo**MEDICAL