



auxoMEDICAL

**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

-

2020



ECG

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.



ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: support@auxomedical.com

Table of Contents

ECG..... 1-3



auxoMEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

New & Refurbished
Medical Equipment

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

Atria 6100

The Atria 6100 is a practical, portable ECG / EKG machine for hospitals and cardiology clinics. The Atria 6100 portable ECG helps streamline workflow with a friendly user interface and clear, readable color display for easy reading of ECG waveform data.



- Large color monitor to view 3-, 6-, and 12-lead ECG printouts in real-time
- Portable 6100 ECG Stores up to 300 electrocardiograms
- Manual or automatic 3-, 4-, 6-, and 12-lead ECG printouts
- Standard with wired Ethernet, wireless 802.11, and Bluetooth; complies with FAA requirements
- Bi-directional communication using HeartCentrix connectivity software
- AccuPrint feature restricts ECG printing if lead wires are not properly attached

Condition	Refurbished
-----------	-------------

Burdick 8300

The Burdick 8300 is a practical, no-nonsense ECG device. The Burdick 8300 ECG / EKG has a friendly user interface and clear, readable 5.7-inch monochrome display to make reading data easy, with protocol choices designed to streamline your work flow.



- ECG Device stores 50 ECGs
- Manual or automatic 3-, 4-, 6-, and 12-lead ECG printouts
- Standard with wired Ethernet and compatible with wireless 802.11.
- Bi-directional communication with our HeartCentrix connectivity software
- AccuPrint feature restricts ECG printing if lead wires are not properly attached

Condition	Refurbished
-----------	-------------

Edan SE-1010, PC-based

PC-based ECG SE-1010 provides you an affordable and portable 12-channel ECG solution with powerful data analysis functions. Optional advanced analysis functions secure clinical accuracy and maximum decision-making efficiency.




Features:

- Simultaneous 12 Leads acquisition handheld design
- Real-time ECG waveforms display on PC screen
- Auto/Manual measurement and interpretation
- WORD, PDF and JPG format reports transferred via E-mail
- Optional advanced analysis functions: Timed vectorcardiogram, VCG(3D), SAEKG, QT interval dispersion, HRT, HFECG, Frequency Spectrum cardiogram
- 30s-900s HRV (heart rate variability) sample and analysis
- Stress test and exercise blood pressure monitor (option)

Condition	New
-----------	-----

Edan SE-12 Express ECG

The twelve-channel ECG SE-12 Express is a special system of smart ECG series, which support both resting ECG and Exercise ECG. Alphanumeric keyboard and one touch operation design makes it the perfect choice for you.



- Features:
- Large high resolution color touch screen
 - Alphanumeric keyboard
 - One-touch operation
 - Enhanced weak signal detection
 - Comprehensive filters & anti-noise technology
 - Real-time ECG waveforms freezing and review
 - MAdvanced SEMIP algorithm for automatic measurements and interpretation
 - Extended printing
 - Large internal memory which can be enlarged by USB flash disk
 - Built-in data management with reports review

Condition	New
-----------	-----

Edan SE-1200 and SE-1200 Express ECG

Excellent cost performance with built-in high-resolution thermal printer. You can record 12 ECG traces on A4 size and US Letter-size paper.

SE-1200 Features:• Alphanumeric keyboard and one touch operation• Sleep mode to save energy and extend LCD life• Automatic measurement and interpretation tested with authoritative CSE database• Real-time waveform freezing• 120 seconds ECG waveform reviewable• Two-step exercise test with periodic recording• Built-in rechargeable Li-ion battery• 100 ECGs in internal memory• Optional: PC-based data management software (Smart ECG Viewer)SE-1200 Express Features:• 8.4" color touch-screen• Alphanumeric keyboard and one touch operation• Sleep mode to save energy and extend LCD life• Automatic measurement and interpretation tested with authoritative CSE database• Real-time waveform freezing• 120 seconds ECG waveform reviewable• Two-step exercise test with periodic recording• Built-in rechargeable Li-ion battery• 200 ECGs in internal memory• Optional: PC-based data management software (Smart ECG Viewer)



Condition	New
-----------	-----

Edan SE-601 Series ECG

Six-channel ECG SE-601 series : SE-601A, SE-601B, SE-601C, cost-effective solutions including hospital wards, routine check, home visit and emergency cases.

Features:• Compact and lightweight design• One-touch operation• Alphanumeric keyboard• 12 leads simultaneous acquisition• Enhanced weak signal detection• Comprehensive filters & anti-noise technology Models: SE-601A• 3.5 inch STN LCD, 192*64 dots
• 50 ECGs internal storage• SE-601B• 5.7 inch STN LCD, 320*240 dots
• 100 ECGs internal storage• SE-601C• 5.7 inch TFT color LCD, 640*480 dots
• 100 ECGs internal storage



Condition	New
-----------	-----

GE MAC 1200

The MAC* 1200 digital, 12-lead electrocardiograph system offers comprehensive ECG solutions with practical features that help meet the needs of hospitals, clinics, office based practices and clinical trials. Its advanced algorithm capabilities, seamless connectivity to the CS and MUSE* Cardiology Information Systems and easy-to-use features provide the quality and convenience users require.



Condition	Refurbished
-----------	-------------

GE Mac 5000 EKG

MAC® 5000, GE's premier ECG system, delivers advanced disease management capabilities through its vast collection of proven algorithms. The MAC 5000 system offers the sophistication required for advanced ECG applications, while its ease of use extends this level of performance to the broadest range of users possible. And it's part of the complete GE suite of networked, non-invasive testing solutions designed to maximize patient throughput and department productivity. Add the stress option, and you combine the industry's leading high performance resting ECG system with comprehensive, equally advanced exercise-testing capabilities in a single, compact system.



- Mac 5000
- Cam 14 Module
- Module Cable
- Battery
- Set of Leads
- Interpretive Electrocardiograph with 12-lead analysis
- Waveform Display
- Full-sized Paper Format
- Alphanumeric Keypad

Condition	Refurbished
-----------	-------------

GE MAC 5500

Built on GE innovation in ECG acquisition and analysis, the MAC 5500 HD Resting ECG System is a premier ECG system, delivering advanced disease management capabilities through one of the industry-leading collections of algorithms and advanced networking.



The MAC 5500 HD represents a new definition of ECG speed and accuracy, and offers:

- Better triaging tools for Acute Coronary Syndrome (ACS) detection
- Configurable Critical Values to help speed time to treatment
- More accurate pacemaker detection for advanced interpretation of paced rhythms
- Enhanced barcode workflow, with orders and ADT download capability, work together to help improve efficiency across the enterprise
- Connected ECG Workflow to help save time and avoid costly misdiagnoses

The MAC 5500 HD also meets the need for productivity across the enterprise, by empowering a wide range of users to precisely capture, analyze, and communicate ECG data. All to ensure the right data goes with the right patient to help reduce demographic and billing errors, and speed time to treatment.

Sophisticated tools. Easy to use. The MAC 5500 HD enables clinicians to focus on what's most important: their patients.

Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Mortara ELI 250 ECG Machine



Compact and lightweight, the ELI 250c Electrocardiograph provides comprehensive functionality in an easy-to-use, portable device. Multiple connectivity options ensure patient data transfers easily from the ECG to the location of your choice. Experience the benefits of bidirectional communication via internal modem, LAN or wireless LAN with Mortara's E-Scribe™ and Athena products, as well as third-party EMR, HIS, and PACS systems via XML, PDF and DICOM®. Best 10 technology helps you capture the best 10 seconds of ECG data from system memory, reducing the need for repeat ECGs. Wireless data acquisition helps increase workspace flexibility by gathering 12-lead data without the need for a tethered cable.

Features:

- High-res color display offers 12-lead ECG preview and post-acquisition review
- Choose between WAM™ wireless acquisition module or the AM12™ acquisition module
- VERITAS™ algorithm uses gender-specific, adult/ped criteria for interpretation
- Buttons for ECG acquisition, rhythm printing & ECG transmission/order retrieval
- Chest and limb lead reversal indicators alert clinician to possible errors
- Communicates with Mortara ECG Safe™ cloud-based service
- Supports pacemaker spike detection with 40,000 sps/channel data capture
- 1k sps/channel capture for signal assessment during analysis, display & printing

Condition	Refurbished
-----------	-------------

Welch Allyn CP150



The Welch Allyn CP 150™ Resting Electrocardiograph (ECG) with optional spirometry combines the best traditional box ECG features like full-size printouts and clinical decision with a simple touchscreen interface, fast one-button operation, and flexible connectivity options to help you keep your staff productive and your tests accurate. And with the optional spirometry module, the CP 150 even lets you perform two reimbursable tests with just one device. The end result is a device that will help your facility save time, improve workflow, and efficiently manage patient information.

- 7" Color Touch Screen
- Advanced filters for optimal ECG trace quality
- Three (3) lead ECG Preview
- Battery operation for ease of transport (up to 8 hours) - Lithium Battery
- Store up to 100 results on device or easily transfer to USB memory stick

Condition	Refurbished
-----------	-------------

Welch Allyn Schiller AT-2 Plus



SCHILLER's AT-2 light offers the economics of a low cost unit with many of the high end features needed to care for your patients. The AT-2 light is a 12-channel ECG system offering battery operation of up to 6 hours, full size (8 1/2 x 11) reports to fit patient records, LCD screen and keyboard for easy entry of patient demographics and interpretation.

Switch between EKG and spirometry with the push of a button with this advanced yet economical combination unit. All the features of the AT-2 plus above, including preview screen and full-page reports, along with Schiller's advanced spirometry technology. Export data to PC with RS-232 interface.

- Alphanumeric keyboard
- 6/12 channel representation of all 12 simultaneously acquired leads
- Pulmonary Function Testing for FVC, SVC, MV and MVV
- Large, high resolution back-lit monitor for easy preview of ECG quality
- Automatic and manual mode
- Full size EKG printouts

Condition	Refurbished
-----------	-------------



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers
Washer Disinfectors
Surgical & Exam Lighting
Surgical & Exam Tables Stretchers
& Patient Transportation
Vital Sign & Patient Monitors
Blanket/Fluid Warmers
Scrub Sinks & Accessories
Defibrillators
Electrosurgical Units
Surgical Endoscopy
Imaging
Tourniquet Systems
Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product	Page
Atria 6100.....	1
Burdick 8300.....	1
Edan SE-1010, PC-based.....	1
Edan SE-12 Express ECG.....	1
Edan SE-1200 and SE-1200 Express ECG.....	2
Edan SE-601 Series ECG.....	2
GE MAC 1200.....	2
GE MAC 5500.....	2
GE Mac 5000 EKG.....	2
Mortara ELI 250 ECG Machine.....	3
Welch Allyn CP150.....	3
Welch Allyn Schiller AT-2 Plus.....	3

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER

Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's
Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions

Surgery Center Services of America
Chesapeake Healthcare Planning
Sodexo USA



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's
Daughters
Virginia Eye Institute

OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTROSTEAM
SUSSMAN

Washer Disinfectant

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPERTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBREOPTICS
WELCH ALLYN

EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424
www.auxomedical.com
sales@auxomedical.com

 **auxo**MEDICAL