



auxoMEDICAL

**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

—
2020



Hospital Beds

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.

ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.

Table of Contents

Hospital Beds.....	1-6
--------------------	-----



auxo MEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

**New & Refurbished
Medical Equipment**

Alterra 1232 Hi-Low Bed

Medline Alterra 1232 Hi-Low Full Electric Hospital Bed provides safety and offers a multitude of comfortable positioning options. The maximum height of 26 inches (66cm) makes it easy to care for residents and the low height of 7.25 inches (18cm) makes it ideal for Alzheimer's patients and any other resident at risk of falling. The bed moves from its lowest to highest height without any horizontal movement, eliminating wall damage. Battery back-up (included) keeps bed functioning during power outages or evacuations. Features:

Built for ease of use and safety for the caregiver?The maximum height of 26" (66 cm) makes it easy to care for residents; the low height of 7.25" (18 cm) makes it ideal for Alzheimer's patients and any other resident at risk of falling?The bed moves from its lowest to highest height within 47 seconds and without any horizontal movement, eliminating wall damage?Battery backup keeps the bed functioning in the event of a power outage (battery back up is not available on the Alterra 1200 series)?Smooth performance with ultra-quiet, dependable DC motors - 2 for height, 1 for the knee hinge and 1 for the reclining head deck?Combination grid/rib deck for enhanced comfort and strength?Bed deck dimensions: 80" (2 m) long, 34.5" (88 cm); Overall length: 87.5" (2.2 m)? 500 lb weight capacity? Weight: 302 lbs.



Condition	Refurbished
-----------	-------------

Hill-Rom Advance Series

The Hill-Rom® Advance® Series hospital bed is perfect for hospitals & surgery centers. This bed offers a variety of features that are perfect for your patient care needs.

Specifications:

- Three motors
- High-low travel time: 21 seconds
- Sleep surface: 36" (w) x 82" (l)
- Overall dimensions: 36" (w) x 91" (l)
- Bed frame clearance with 5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central pedals for braking and steering
- Automatic night light
- Instant CPR release



Condition	Refurbished
-----------	-------------

Hill-Rom Advanta 2 Bed

The Advanta™ 2 Med-Surg Bed offers improved ease of use through innovation and proven performance. Designed to help caregivers accomplish their goals, the Advanta 2 Med-Surg Bed provides patients egress/ingress, less patient positioning, and more convenience through storage and other features.

- Smart Bed technology enabled for connectivity with NaviCare® Nurse Call, other nurse call systems and EMRs.
- SafeView® alerts system
- Bed exit monitoring system with Hill-Rom's patented three-level detection algorithms
- Nightlight
- Audible "brake not set" alarm
- SlideGuard® bed articulation technology with Auto Contour™ feature
- Electric vascular foot position
- Care Grip patient egress assist handle
- Integrated scale
- Point-of-Care® nurse controls
- One-button Dining Chair® feature
- IntelliDrive® powered transport or central 5th steering wheel
- Automatic battery backup
- Patient in-rail storage compartment
- Digital head of bed angle display
- 30° head of bed alarm
- Sidecom® nurse call and entertainment controls
- One-hand bed extender
- Integrated Line Manager
- Positioning controls accessible to the patient



Condition	Refurbished
-----------	-------------

Hill-Rom Advanta Series

The Hill-Rom® Advanta® Series hospital bed is perfect for hospitals & surgery centers. This bed offers a variety of features that are perfect for your patient care needs.

Specifications:

- Controls at the bedside
- Dual locking brakes in all four corners
- Point-of-Care controls with large buttons
- Ergonomic grip handles
- Robust bed exit alert system
- Nurse call and wireless connections
- Convenient patient controls
- Visual safety reminders
- One-button low chair position



Condition	Refurbished
-----------	-------------

Hospital Beds

Hill-ROM Affinity 3 Birthing Bed

In labor and delivery, time and safety are top of mind and every second counts. We have designed every detail of the Affinity 3 birthing bed to help make the birthing process safe for the patient and caregiver. Additional back saving features like the no lift Stow and Go™ foot section and ergonomic push handles make the Affinity 3 a smart choice for caregivers.

Hill-ROM Affinity 3 Birthing Bed- Battery back-up- Automatic night light- Automatic seat tilt- Tente locking casters- Fluid basin- OneStep activation-Instant CPR release- Instant labor grips- Siderail release- Trend-like positioning- Foot supports- Foot section removal- Central brake & steer

Specifications: Total Length: 90" (229 cm). Length: From Roller Bumpers to Break in Seat Section: 61.88" From Break in Seat Section to End of Bed: 28.13" (Maximum Width: Siderails Stored: 36" Siderails Up: 39" Headboard Height (Maximum): 52" Siderail Height (Maximum): 14.75" Minimum Under-Bed Clearance: 5" Wheel Base: 50" x 29" Casters: 6" Total Weight (Maximum): 480 lbs Head Section Inclination (Maximum): 65°. Seat Section Inclination (Maximum): 15°. Trend-Like Position (Maximum): 8°. Bed Height Range: 18 - 34" (46 - 86 cm). Bed Lift Capacity (Max Safe Working Load): 500 lbs Foot Section Lift Capacity (Maximum): 400 lbs Head Section Lift Capacity (Maximum): 200 lbs



Hill-Rom Affinity 4 Birthing Bed

The Hill-Rom Affinity IV (Model P3700) is intended for use as a birthing bed within the acute care Labor and Delivery Unit. Hill-Rom has continued to make refinements to this popular model birthing bed. The Affinity IV introduces new design upgrades and optional features not available on the Affinity III model such as the Stow-and-Go™ Retractable foot section - convenient, stow away foot section stores under bed, eliminating lifting and saving caregiver valuable space in birthing suite; and EasyGlide™ Calf supports - wider support with strong release lever that can be actuated from any angle with one hand. The Affinity IV is a versatile, easy-to-use birthing bed that has multiple articulation features, ergonomic hand grips, motorized bed positioning, manually adjustable foot stirrups and calf supports, standard light-weight lift-off foot section for caregiver access, fluid basin, and for patient comfort a two-piece head and foot section mattress or optional adjustable air mattress system.



Standard Features:

- 2-Section Foam Head and Seat, and Foam Foot section mattress set with mattress covers.
- Point-of-Care™ caregiver and patient side rail controls (Bed hi/lo, backrest up/down and foot up/down)
- Dual ergonomic comfort foam gripped labor delivery handles
- Dual comfort foam gripped Trendelenburg handles
- Gas-assist Trendelenburg
- Bed positioning control lock-out switch
- Headboard with integrated push/steering handles
- Four-wheel braking system, and 1-wheel steering caster with new design dual centralized brake/steer pedal controls
- CPR release: emergency lowering
- Automatic 15-degree pelvic tilt
- Rechargeable battery backup with LED indicators
- OneStep™ side rail release (tuck-away side rails)
- Zero gap transfer
- Fluid basin
- Adjustable foot supports
- Two IV pole sockets (headend)
- 6" (15 cm) Tente™ casters
- Night light
- Hospital grade power plug
- 500 lb (227 kg) safe working weight capacity; Foot section 400 lbs

Condition	Refurbished
-----------	-------------

Hill-rom Affinity 4 Birthing Bed

The labor and delivery experience can be unpredictable for everyone involved. To help manage this, Hill-Rom offers a solution that can adapt and deliver throughout the birthing process.

The Affinity® 4 Birthing Bed is designed to support mom, baby, and you throughout the labor and delivery experience.

No-Lift Stow-And-Go Foot SectionAutomatic Seat Tilt, Gas-Assist Trendelenburg PositioningFluid Basin, Easy-Glide Calf SupportsOneStep Siderail ReleaseInstant CPR ReleaseCentral Brake & SteerLong-Life DC MotorMattress Included With Purchase



Condition	Refurbished
-----------	-------------

Hill-Rom Affinity II & III

The Affinity® Three Birthing Bed builds upon Hill-Rom's long-standing success of its premier Affinity® Birthing Bed line and improves hospital flexibility, efficiency, positioning, and deliver with a single product. And, its wide range of innovative features make it state-of-the-art in birthing beds. The Affinity® Three Birthing Bed is an ideal solution for improved comfort and satisfaction.

Today, many women have their first experience at a hospital through the obstetrics department. With more options for labor than ever before, you need to create a positive experience. The Affinity® Three Birthing Bed is made to accommodate virtually any choices your mothers prefer ... and it's built for superior comfort no matter what decisions they make. At the time to break down the bed for use, the intuitive design of the Affinity® Three Birthing Bed makes for a fast and easy transition. Its simplified OneStep™ features allow a single caregiver to perform tasks that previously required multiple steps. When it becomes necessary to transport during childbirth, the Affinity® Three Birthing Bed is one of the industry's most mobile beds. The easy-to-maneuver bed is designed to be the lightest birthing bed available. So, while the mother is in transit, caregivers can offer continued, confident care depending on the bed to improve mobility. The safety of your patients is of the utmost importance, and a sense of security may make a difference in their experience. The Affinity® Three Birthing Bed includes features to address safety that both your patients and their caregivers will appreciate. Specifications:

- Length: 90" (2290 mm)
- Width: 39" (990 mm)
- Bed Height: 38" (970 mm)
- Mattress Length: 78" (1980 mm)
- Mattress Width: 35" (890 mm)
- Maximum Head Elevation: 63°
- Maximum Lift Capacity: 500 lbs. (227 kg)
- Caster Size: 6" (150 mm) / 8" (200 mm)
- Weight: 450-460 lbs. (204-208 kg)



Condition	Refurbished
-----------	-------------

Hill-rom Burke Triflex II Bariatric Hospital Bed

The Tri-Flex II Bariatric bed features a frame designed to provide dependable care. This bed can accommodate most users needs and is available with optional foam and low air loss surfaces.

Features

Expandable width frame and surface adjust for greater user comfortPatient controls located on the siderail (head and knee)Trendelenburg - allows the user to place the bed in a tilted position which enables their head/neck to be lower than their lower extremitiesReverse Trendelenburg - allows the user to place the bed in a tilted position which enables their lower extremities to be lower than their head/neckChair Positioning - allows the user to be in a fully seated position in the bedModular design for easy assemblyBattery Back up - if a power outage occurs the battery back up can be used to operate the bed



Condition	Refurbished
-----------	-------------

Hill-rom CareAssist ES Medical Surgical Bed

The CareAssist® ES Medical Surgical Bed is a fully electric hospital bed providing essential and dependable technology in your home. Bring your hospital bed home and be assured of comfortable and safety. The CareAssist® ES Med Surg Bed provides easy to use solutions to help caregivers manage patient safety, positioning, easy maintenance and comfort. With or without using a patient transfer sheet, the CareAssist expands to meet your needs, you can use either an 80 inch or 84 inch long mattress.



Key Product Information:

- Compatible with most powered and non-powered mattresses
- Trendelenburg (foot of bed raised) positioning up to 16 degrees
- Top and bottom half rails integrated into the bed
- SafeView® alerts provide status of bed
- Hand held pendant clips into the rail
- Nightlight and battery backup
- Audible "brake not set" alarm
- SlideGuard® bed articulation technology with auto contour
- Electric vascular foot position
- One-button dining chair position
- Four clocking rails allows you to independently raise or lower any of the side rail sections
- One-hand bed extender, Bed extends from 80 inches to 86

Hill-Rom® CareAssist® ES Medical Surgical Bed Specifications:

Bed Height:	Low position floor to top of deck: 15.75" (42.55 cm)
High position floor to top of deck:	32.5" (82.55 cm)
Trendelenburg:	Both trend and reverse trendelenburg
Sleep surface Width:	36.00" (91.44 cm)
Sleep surface normal length:	80.00" (203.20 cm)
Sleep surface extended length:	86.00" (218.44 cm)
Overall Width:	40" (102 cm)
Overall normal length:	100" (254 cm)
Overall extended length:	94.04" (238.86 cm)
Max head of bed elevation:	65°
Max knee elevation:	20°
Max Trend and reverse Trend:	16°
Weight Limit:	500 lbs
Technical and Quality Assurance Standards:	UL 60601-1, CSA® a C22.2 No. 601.1, IEC 60601-2-38, IEC 60601-1, IEC 60601-1, IEC 60601-1-2, EN ISO 9001 and EN 13485

Unit comes with Standard Comfort Mattress

Condition	Refurbished
-----------	-------------

Hospital Beds

Hill-Rom Centra Series 850 & 852

The Hill-Rom® Centra® Series 850 & 852 hospital beds are perfect for hospitals & surgery centers. These beds offers a variety of features that are perfect for your patient care needs. Specifications:
Three motors• High-low travel time: 21 seconds• Sleep surface: 36" x 82"• Bed frame clearance with 5" casters: 7.5"• Low position: 17"• High position: 30"• Central pedal for braking and steering• Night light• CPR release• Hill-Rom Century Series 835 & 837



Condition	Refurbished
-----------	-------------

Hill-Rom TotalCare Bed

The TotalCare® P500 therapy bed delivers the versatility to also respond to the needs of the immobile patients with a moderate to high risk for pressure ulcer development. The best strategy to avoid immobility related complications is to get patients mobile as early as possible. Aside from the Progressive Mobility™ support coming with the TotalCare® bed functionality.



The surface technology addresses pressure, shear, friction and heat—the four main extrinsic causes of pressure ulcers. Its weight based pressure redistribution system permanently communicates with the bed's scale system to generate the optimum pressure redistribution settings in four different zones and in any bed position.

Advanced Microclimate® Technology removes excess heat and moisture from the surface for cooler, dryer skin, helping to prevent moisture lesions and increase patient comfort. A rapid response turn assist with a reminder function helps achieving turn protocol compliance while relieving strain and discomfort for caregivers and patients.

Features:

- Weight based pressure redistribution in any bed position
- Advanced Microclimate® technology
- 3-layer Microclimate Management™ surface cover with in-built X-ray sleeve
- Point-Of-Care® siderail controls with Graphical Caregiver Interface feature (GCI).
- Digital Head of bed angle indicator and alarm
- 3 mode bed exit sensor with alarm silence and alarm suspend modes
- FullChair® and chair egress with stand assist
- One-button Boost® patient repositioning support
- FlexaFoot™ electrical bed retraction and extension
- HandsFree® emergency CPR and Trendelenburg
- History data on patient weight, head of bed elevation and chair position
- Turn Assist with reminder function
- Optional IntelliDrive® powered transport
- Emergency Trendelenburg/Reverse Trendelenburg: 20°/20°(Siderail Trendelenburg/reverse Trendelenburg is 15°)

Condition	Refurbished
-----------	-------------

Hill-Rom VersaCare

The VersaCare® bed is a full-featured patient care bed system that offers the lowest height in acute care, an integrated air surface that helps automate the delivery of wound care therapy, and many features to reduce the risk of caregiver injury. In addition, this latest innovation from Hill-Rom® offers superior patient comfort and intuitive controls to help caregivers effectively and easily operate the system to its fullest capacity.



Med-Surg units include patients with a wide variety of needs. Some are awake and mobile, sitting up and getting out of bed with little or no assistance. Others need more help. The VersaCare® bed is designed to help prevent or reduce the impact of falls for patients at a range of acuity levels.

For more active patients, the VersaCare® bed's lower position – the lowest height of any acute care bed – means it's safer to get in and out of bed. For higher acuity patients, the VersaCare® bed makes it easy to position patients up in chair without leaving the bed. And with the LowChair® feature, caregivers can keep the bed at the same low height while placing patients in a chair position, without compromising on the patient's freedom of movement or the caregiver's time.

The VersaCare® bed's unique surface shape supports the patient as he or she sits on the edge of the surface, providing firm support for egress. Nightlights at the side of the bed are focused on where the patient's feet will land. Comfortable, easy-to-grip siderails protect without restraining.

Specifications:

- Length: 94.5" (2400 mm)
- Width: 40" (1020 mm)
- Minimum Bed Height: 18" (457 mm)
- Maximum Bed Height: 37" (940 mm)
- Mattress Length: 86" (2180 mm)
- Mattress Width: 35.5" (902 mm)
- Maximum Head Elevation: 65°
- Maximum Lift Capacity: 500 lbs. (227 kg)
- Caster Size: 5" (127 mm)

Condition	Refurbished
-----------	-------------

Did you know that we provide maintenance and service for every piece of equipment listed in this catalog?

Stryker GoBed II Med/Surg Bed w. Scale

Stryker's GoBed II Med/Surg bed was developed with extensive input and feedback from nursing professionals to add efficiency and mobility to patient care practices. The GoBed II is a motorized, fully-electric equipped bed with numerous one-button or one-handed controls allowing caregivers to spend their time helping patients recover quickly in a secure and comforting environment. This bed features full length siderail coverage to aid in patient safety, a low bed height of 14.5 inches promoting proper patient positioning for safe entry and exit, integrated pump holder conveniently located on footboard, and 120-volt outlet to help maximize space at the foot end of the bed. This bed is ideal for healthcare facilities for Med/Surg purposes, educational and training facilities and residential settings.



14.5-inch (37 cm) low-bed height• Retractable fifth-wheel steering• StayPut™ retractable bed frame design• 6-inch (15 cm) casters• Trend./reverse Trend. • Four independent electric motors• Trend. display on footboard (only available with scale option)• Centrally-located steer and four-wheel brake mechanism• Four drainage bag hooks• Eight IV pole/traction equipment sockets• Roller bumpers• Manual back-up for head and knee control• Integrated pump holder• Nurse controls on footboard and siderails• Fixed patient controls on siderails• Degree indicator for head elevation• Patient restraint locations• Electronic function lockout controls• One-handed dampened siderail release• Auto contour• Photo-sensitive night light

Specs:Overall Length - 94.25"Overall Width - 40"Weight Capacity - 500lbs.Height Range - High (29") Low (14.5")Litter posititoning Backrest 0°– 60° (siderails down)Knee Gatch 0°– 28°Trendelenburg/Rev. - ±14°Patient Surface - 36" x 84"Unit has Scale

Condition	Refurbished
-----------	-------------

Stryker LD304 Birthing Bed

The Stryker LD304 Maternity Birthing Bed provides comfort and convenience for both patients and clinicians. With convenient controls, a superior steering system and automatic pelvic tilt, caregivers can easily guide patients through the birthing process without back strain. The ultra-plus mattress keeps patients comfortable. Includes a wooden headboard and an under-bed night light.



DimensionsOverall Length - 93" (236 cm)Overall WidthSiderails Up: 41" (104 cm)Siderails Down: 37" (94 cm)Lift CapacityBed - 500 lb. (228 kg)Head Section - 250 lb. (114 kg)Foot Section - 300 lb. (136 kg)Height Range (to litter top)High - 35.5" (90 cm)Low - 17.5" (44 cm)Litter PositioningBackrest - 0° – 70°Pelvic Tilt - 4°Maximum Trend. - 8°Patient Surface - 33" x 81" (84 x 206 cm)Mattress Thickness - 3" (8 cm)Seat Depth - 12.5" (32 cm)Caster Diameter - 6" (15 cm)

Condition	Refurbished
-----------	-------------

Stryker Secure 2 bed w/ Mattress

The Secure II is a full-featured medical and surgical bed that combines practical engineering with smart bed technology to help provide a patient platform that not only meets the patients' needs, but also the needs of those who care for them.



Features such as a retractable bed frame, centrally-located brake, intermediate side rail positioning and open architecture design help deliver ease, flexibility and an unparalleled level of safety for both patients and caregivers.

Additional Features:Single-button, in bed scale provides accurate, repeatable weights with the patient in any positionBed Awareness option monitors brake, siderail position, bed exit and low bed height, and provides proactive alert if parameters changeFour foley bag hooks, and four IV receptaclesPatient restraint locationsCompletely reconditioned and works like new

Stryker Scure II Med /Surg Bed - include:MattressHead and FootboardsFull Trendelenburg and reverse Trendelenburg

Condition	Refurbished
-----------	-------------

Hospital Beds

Stryker Secure 3 Med/Surg Bed

Stryker's current S3 med/surg bed replaces the popular Secure II (Model 3005) bed being one-generation newer incorporating the latest technological improvements to the bed design and features. The S3 is a versatile, fully electric, adjustable Med/Surg hospital bed that features an open architecture with one-hand intermediate and foot side rails that have both caregiver and patient bed position controls, low bed height, and retractable bed frame that keeps patients comfortable regardless of bed articulation, removable head and footboards that include caregiver bed position controls and patient lockouts, an integrated one-touch scale to easily weigh the patient, Stryker's patented Chaperone® bed exit system and Zone Control® technology and downward motion interrupt for patient safety.

Another optional safety feature is iBed Awareness which monitors local bed status information, alerting caregivers visually, audibly or remotely if preset parameter are compromised. By providing better bedside data and greater flexibility to connect, caregiver's can improve care efficiencies to help prevent adverse events. The S3® is an ideal choice for healthcare facilities, educational and training facilities, and home-use due to its durability, functionality, versatility and ease of use.



Features

Centrally-located brake pedal

4-wheel, steel brake with swivel lock system

6" (15 cm) casters • One-touch 30-degree head of bed angle and lock

One-hand, intermediate siderails

StayPut bed frame design

Open architecture – flat litter deck

Fowler angle gauge

Four IV pole and traction equipment socks

Four foley bag hooks

Steer control

Integrated pump holder

Foot prop rod

Downward motion interrupt

CPR release and removable headboard

Dual pedestal lift design

Three-sided bed motion controls

One-button cardiac chair position

Foot end lock-out controls

Foot end digital trend. and fowler angle display

Trend./reverse trend. controls

Neutral accent panels

Night light

iBed upgradable

DC motor

Specification
Overall length: 95" (241.3 cm)
Overall width:
Siderails Up: 41.5" (105.4 cm)
Siderails Down: 39.5" (100.3 cm)
Weight capacity: 500 lb. (227 kg)
Height range: (to little top)
High: 30" (76.2 cm)
Low: (standard with 6" casters) 16" (40.6 cm)
Litter positioning:
Backrest: 0°-60°
Knee gatch: 0°-40°
Trend./reverse trend: +12°/-10°
Patient surface: 35" x 84" (89 x 213 cm)
Caster diameter: 6" (15 cm)

Condition	Refurbished
-----------	-------------

Did you know that we provide maintenance and service for every piece of equipment listed in this catalog?



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers

Washer Disinfectors

Surgical & Exam Lighting

Surgical & Exam Tables Stretchers
& Patient Transportation

Vital Sign & Patient Monitors

Blanket/Fluid Warmers

Scrub Sinks & Accessories

Defibrillators

Electrosurgical Units

Surgical Endoscopy

Imaging

Tourniquet Systems

Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxo MEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product Index

Product	Page
Alterra 1232 Hi-Low Bed.....	1
Hill-ROM Affinity 3 Birthing Bed.....	2
Hill-Rom Advance Series.....	1
Hill-Rom Advanta 2 Bed.....	1
Hill-Rom Advanta Series.....	1
Hill-Rom Affinity 4 Birthing Bed.....	2
Hill-Rom Affinity II & III.....	3
Hill-Rom Centra Series 850 & 852.....	4
Hill-Rom TotalCare Bed.....	4
Hill-Rom VersaCare.....	4
Hill-rom Affinity 4 Birthing Bed.....	3
Hill-rom Burke Triflex II Bariatric Hospital Bed.....	3
Hill-rom CareAssist ES Medical Surgical Bed.....	3
Stryker GoBed II Med/Surg Bed w. Scale.....	5
Stryker LD304 Birthing Bed.....	5
Stryker Secure 2 bed w/ Mattress.....	5
Stryker Secure 3 Med/Surg Bed.....	6

Did you know that we provide maintenance and service
for every piece of equipment listed in this catalog?



Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's Daughters
Virginia Eye Institute



OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTRASTEAM
SUSSMAN

Washer Disinfector

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBEROPTICS
WELCH ALLYN
EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE

WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424

www.auxomedical.com
sales@auxomedical.com

 auxo MEDICAL