



auxoMEDICAL



**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

-

2020

Sterilizers

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.



ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: support@auxomedical.com

Table of Contents

Sterilizers..... 1-4



auxoMEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

New & Refurbished
Medical Equipment

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

Skytron Integrity 175 Sterilizer

Building on Skytron's history of quality and innovation, the 175 sterilizer is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

Features:

- Spacious 171.72" x 171.72" x 261.2" sterilization chamber
- Processes up to two 25-pound trays
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered sliding door
- Recessed foot-pedal and touchscreen door control
- Stainless steel and copper piping



Condition	New
-----------	-----

Skytron Integrity 215 Sterilizer

Building on Skytron's history of quality and innovation, the 215 sterilizer is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

Features:

- Spacious 211.72" x 211.72" x 38" sterilization chamber
- Processes up to six 25-pound trays (150 pounds per cycle)
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered sliding door
- Recessed foot-pedal and touchscreen door control
- Stainless steel and copper piping
- Dual-walled, fully jacketed stainless steel pressure vessel



Condition	New
-----------	-----

Skytron Integrity 270 Sterilizer

Building on Skytron's history of quality and innovation, the Integrity 270 sterilizer is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

Features:

- Processes up to sixteen 25-pound trays (400 pounds per cycle)
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered door locking/unlocking mechanism
- Stainless steel and copper piping
- Dual-walled, fully jacketed stainless steel pressure vessel



Condition	New
-----------	-----

Statim 2000 Sterilizer

The STATIM 2000 Classic boasts a 6-minute cycle, but has the added advantage of employing a fully removable chamber so you can aseptically transport a set of freshly sterilized instruments directly to the procedure area.

The STATIM 2000 is fully automatic. Simply select one of the three sterilization programs and with a single touch of the keypad, begin the cycle. An advanced microprocessor controls all functions and maintains optimal sterilization conditions, from start to finish.

The STATIM 2000 is a reliable cassette sterilizer that's compact enough to fit in procedure areas and fast enough to sterilize jewelry, expensive Lasik, PRK, endoscopic and cataract instruments, including phaco hand pieces, between procedures - it can complete an unwrapped sterilization cycle in just six minutes and a wrapped cycle in 14 minutes.

- No plumbing required. Integrated reservoir.
- Front panel controls.
- Low water cut-off. This prevents burnout of boiler elements.
- Safety switch automatically activates when pressurized
- Chamber Volume: 1/2 Gallon (1.9 Liters)
- Reservoir: 1 Gallon (3.8 Liters)
- Cycle Time: 6 minutes



Steris 400 Series 16" Sterilizer

The completely refurbished V116 steam sterilizer has an interior chamber size of 16x16x26". Includes a rack & two shelf loading system. Single vertical power door (foot pedal actuated) with touch screen operation and printer.

Touch-screen with 30-line x 40-character display area• Ink-on-paper impact printer• Help screens for programming and troubleshooting alarm conditions• Automatic check of control program and cycle data maintains process integrity• Service reprogrammable flash ROM memory

16 x 16" sterilizers – chamber clearance is 12" (304 mm) for top shelf and 14" (357 mm) for bottom shelf.

20 x 20" sterilizers – chamber clearance is 18" (457 mm) for intermediate shelves, and 15" (381 mm) for bottom shelf.

Optional Stainless steel side panels available.Optional 30kW steam generator for electrical units.



Condition	Refurbished
-----------	-------------

Steris 400 Series 20" Sterilizer

The completely refurbished V120 steam sterilizer has an interior chamber size of 20x20x38". Includes a rack & two shelf loading system. Single vertical power door (foot pedal actuated) with touch screen operation and printer.Touch-screen with 30-line x 40-character display area

- Ink-on-paper impact printer
- Help screens for programming and troubleshooting alarm conditions
- Automatic check of control program and cycle data maintains process integrity
- Service reprogrammable flash ROM memory

20 x 20" sterilizers – chamber clearance is 18" (457 mm) for intermediate shelves, and 15" (381 mm) for bottom shelf.



Condition	Refurbished
-----------	-------------

Steris Amsco Century 1262

26 x 26 x 39" interior chamber dimensions Door is 26"x26", interior is circular thus 34" W Vertical Sliding automatic door

Electrical: Control: 120 VAC, 2A, single-phase, 60 Hz Vacuum Pump: 208/240 VAC, 6A per phase, 3-phase, 60 Hz, or 480 VAC, 3A per phase, 3-phase, 60 Hz

Plumbing: Drain: 2" ODT drain terminal (Floor drain capacity must handle peak water consumption.)

Water: Pressure – 20 to 50 psig (1.38 to 3.45 bar), dynamic Consumption – 130 gal/hr (495 L/hour)

Peak Flow – 15 gal/min (57 L/min) Size – 1" NPT

Steam: Pressure – 50 to 80 psig (3.45 to 5.52 bar), dynamic, condensate free, and 97% to 100% vapor quality Consumption – 148 lb/hr (67 kg/hour) Peak Flow – 255 lb/hr (116 kg/hour) Size – 1" NPT



Condition	Refurbished
-----------	-------------

Steris Amsco Century Medium Series

Fully Refurbished V148, V160 and V136 Century Series Medium steam sterilizers: designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The pre-vacuum sterilizers are equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 26x37.5x48" chamber, the sterilizers provide enough space fit all of your equipment and sterilize with speed and reputable quality in mind.

Features:

- Florescent display for easy reading
- Ink on paper impact printer
- 26x37.5x48" chamber
- Pre-vacuum Steam Sterilizer
- Fast & efficient
- Emergency manual exhaust valve
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Communication interface compatibility w/ most PC devices
- Electronic water conservation valve



Condition	New
-----------	-----

Steris Amsco Century Small Series

The Century® Sterilizer meets the needs of today's O.R. and other areas of the hospital. The Century® Sterilizer features foot pedal operation of the automatic sliding door for hands-free loading and unloading. The vacuum fluorescent tough panel permits direct operator cycle selection. The large display of time remaining during the cycle allows the operator to perform other tasks while awaiting cycle completion.



Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Steris Amsco Century V148

The Steris Amsco Century V148 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V148 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 24" x 36" x 48" chamber, the V148 sterilizer provides enough space fit all of your equipment and sterilize with speed and reputable quality in mind. Features:

- Florescent display for easy reading
- Ink on paper impact printer
- 24" x 36" x 48" chamber
- Pre-vacuum Steam Sterilizer
- Fast & efficient
- Emergency manual exhaust valve
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Communication interface compatibility w/ most PC devices
- Electronic water conservation valve



Condition	Refurbished
-----------	-------------

Steris Amsco Century V160

Fully Refurbished V160 Century Series Medium steam sterilizer: designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V148 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 26" x 37.5" x 60" chamber, the V160 sterilizer provides enough space fit all of your equipment and sterilize with speed and reputable quality in mind.

- Florescent display for easy reading
- Ink on paper impact printer
- 26" x 37.5" x 60" chamber
- Pre-vacuum Steam Sterilizer
- Fast & efficient
- Emergency manual exhaust valve
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Communication interface compatibility w/ most PC devices
- Electronic water conservation valve
- 1 Year Parts & Labor Warranty, 6 months on PM Parts, 5 warranty years on pressure vessel.
- Optional Stainless steel side panels available.
- Standard Preventative Maintenance must be completed while unit is under warranty



Condition	Refurbished
-----------	-------------

Steris Amsco Eagle 3033

Steris Amsco 3033- Completely Refurbished Amsco Eagle 3000 Stage 3 series, 24"x 36" x 36", single door sterilizer Recess mounted, or optional cabinet package. 3033 Pre-vac sterilizer, (3031 variant for gravity operation) with operating temperature of 250-270F.

- Stage 3 controller with thermal printer
- Single Radial arm manual door (optional hinge)
- Recess Mount
- Interior chamber size of 24x36x36"
- Pre-vacuum configuration or Gravity configuration
- Thermal display with touch pad operation
- Hydrostatically tested pressure chamber

All Auxo medical reconditioned steam sterilizers are fully disassembled and rebuilt using all OEM parts and procedures.

Warranty: 1 Year Parts & Labor, 6 months on PM Parts, 5 years on pressure vessel.

Optional Stainless steel side panels and Cart available.

Standard Preventative Maintenance must be completed while unit is under warranty



Condition	Refurbished
-----------	-------------

Steris Amsco Eagle 3043

The Amsco Steris Eagle 3043 offers an economical option to meet the needs of today's OR, bulk processing in Central Services, and Ambulatory Surgery Centers. When quality is paramount but budgets are constrained, you can count on STERIS. Quality versus cost is no longer a choice that has to be made; enjoy both with STERIS Certified, Pre-Owned sterilizers featuring Amsco® Eagle® steam sterilizer technology. Designed for fast, efficient sterilization of heat- and moisture-stable materials, prevacuum sterilizers are equipped with prevacuum, gravity, flash, express, leak test, and daily air removal test cycles.

Features:

- Stage 3 controller with thermal printer
- Single Radial arm manual door (optional hinge)
- Recess Mount
- Interior chamber size of 24x36x48"
- Pre-vacuum configuration or Gravity configuration
- Thermal display with touch pad operation
- Hydrostatically tested pressure chamber
- All Auxo medical reconditioned steam sterilizers are fully disassembled and rebuilt using all OEM parts and procedures.



Condition	Refurbished
-----------	-------------

Steris Amsco Eagle 3053

Steris Amsco 3053- Completely Remanufactured Amsco Eagle 3000 Stage 3 series, 24"x 36" x 48", single door sterilizer Recess mounted, or optional cabinet package. Model 3053. Pre-vac sterilizer, (3051 variant for gravity operation) with operating temperature of 250-270F

- Stage 3 controller with thermal printer
- Dual Radial arms manual door (optional hinge) , Pass through
- Recess Mount with side panels
- Interior chamber size of 24x36x60"
- Pre-vacuum configuration or Gravity configuration
- Thermal display with touch pad operation
- Hydrostatically tested pressure chamber



Condition	Refurbished
-----------	-------------

Steris Amsco V136

Completely refurbished , 24"x 36" x 36", single door sterilizer recess mounted, or optional cabinet package (G136 variant for gravity operation) with operating temperature of 250-270F. Designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V136 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. The Amsco V136 sterilizer comes complete with a one year parts and labor warranty (6 months on preventive maintenance part), five year warranty on the pressure vessel.



- Florescent display for easy reading
- Ink on paper impact printer
- 24" x 36" x 36" chamber
- Rack and 2 shelves
- Recessed
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Electronic water conservation valve

Condition	Refurbished
-----------	-------------

Steris Eagle 3000 Medium Series

Steris Eagle 3000 Stage 3 Series Sterilizers can be used in a hospital or laboratory and offer a choice of two configurations. Gravity sterilizers operate by the gravity downward displacement air removal principle. Vacumatic sterilizers have a mechanical air evacuation system for high-speed, high-volume sterilization. Both types of sterilizers use steam-under-pressure as the sterilizing agent. Sterilizers are furnished ready for connection to building steam supply.



Condition	Refurbished
-----------	-------------

Steris Eagle 3000 Small Series

Eagle® 3000 Stage 3™ Series Sterilizers feature an advanced microcomputer control system, providing latest standard for cycle setup, selection and monitoring. Once settings are made and the cycle is started, microcomputers accurately monitor and control system operations. The control column on this sterilizer includes an easy-to-read vacuum fluorescent display and an integral thermal printer.



Condition	Refurbished
-----------	-------------

Steris V116 Steam Sterilizer

Steris V116 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V116 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles.



Features:

- Florescent display for easy reading
- Ink on paper impact printer
- 16" x 16" x 26" chamber
- Rack and 2 shelves
- Freestanding unit
- Vertical Door
- Pre-vacuum Steam Sterilizer
- Fast & efficient
- Emergency manual exhaust valve
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions

Condition	Refurbished
-----------	-------------

Steris V120 Steam Sterilizer

The Steris Amsco / Century V120 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V120 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 20" x 20" x 38" chamber, the V120 sterilizer provides enough space to fit all of your equipment and sterilize with speed and reputable quality in mind.



Features:

- Maximum reliability, performance and value.
- Prevacuum.
- Easy to use.
- Touchscreen-operated control system.
- Process control system provides continuous cycle monitoring and control, as well as service diagnostic assuring consistent cycle performance.

Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers
Washer Disinfectors
Surgical & Exam Lighting
Surgical & Exam Tables Stretchers
& Patient Transportation
Vital Sign & Patient Monitors
Blanket/Fluid Warmers
Scrub Sinks & Accessories
Defibrillators
Electrosurgical Units
Surgical Endoscopy
Imaging
Tourniquet Systems
Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product	Page
Skytron Integrity 175 Sterilizer.....	1
Skytron Integrity 215 Sterilizer.....	1
Skytron Integrity 270 Sterilizer.....	1
Statim 2000 Sterilizer.....	1
Steris 400 Series 16" Sterilizer.....	2
Steris 400 Series 20" Sterilizer.....	2
Steris Amsco Century 1262.....	2
Steris Amsco Century Medium Series.....	2
Steris Amsco Century Small Series.....	2
Steris Amsco Century V148.....	3
Steris Amsco Century V160.....	3
Steris Amsco Eagle 3033.....	3
Steris Amsco Eagle 3043.....	3
Steris Amsco Eagle 3053.....	4
Steris Amsco V136.....	4
Steris Eagle 3000 Medium Series.....	4
Steris Eagle 3000 Small Series.....	4
Steris V116 Steam Sterilizer.....	4
Steris V120 Steam Sterilizer.....	4

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER

Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's
Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions

Surgery Center Services of America
Chesapeake Healthcare Planning
Sodexo USA



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's
Daughters
Virginia Eye Institute

OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTROSTEAM
SUSSMAN

Washer Disinfectant

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPERTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBROPTICS
WELCH ALLYN

EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424
www.auxomedical.com
sales@auxomedical.com

 **auxo**MEDICAL