



auxoMEDICAL

**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

–
2020



Surgical Lighting

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.

ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.

Table of Contents

Surgical Lighting.....	1-4
------------------------	-----



auxo MEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

**New & Refurbished
Medical Equipment**

Dual VistOR PRO LED Ceiling Mount

The VistOR PRO Procedural Light features a large 16-inch, highly-polished faceted reflector which keeps shadowing to a minimum and maximizes the light output. This pre-focused light also combines a great depth of field, smooth, drift-free positioning and an easy-to-change bulb that doesn't require the use of tools.

Features:

- 64,500 Lux light intensity from a single quartz halogen bulb
- Color correct, extremely cool light output
- Large pre-focused pattern
- 36-inch depth of field
- Large faceted reflector provides shadow free illumination
- Color temperature: 4,000 K
- Depth of Field: 36 inches
- Bulb Life: 2,000 hours
- Bulb Voltage/Wattage: 12VAC/90W
- Available as a Ceiling Mount, Wall Mount, Freestanding Mobile Unit or Track Mounted on the Shuttle Track



Condition	New
-----------	-----

Nuvo VisTOR EX LED Examination Light Wall Mount

Significantly brighter, whiter and cooler than any previous technology—at the greatest value. Large flood style pattern delivering a 6-inch pattern of light. The light can be paired with the VistOR PR and VistOR PRO models to provide a combination of lighting.

Product SpecificationsColor Temperature 4300° K pure white illuminationLight Intensity 55,000 lux at 24" (60 cm)Intensity Level 2-stage dimmingLED Life Efficiency 50,000 hour-rated LED lifeRadiated Heat Near ZeroDetachable Sterile Handle OptionalPattern Size 6.0"Reflector Diameter 5"Depth of Field (in) InfiniteVoltage 120 - 240 VAC 50/60 HzAgency Approvals UL 2601; CUL; CE Mark



Condition	New
-----------	-----

Nuvo VisTOR EX Portable LED Examination Light Floor Mount

Significantly brighter, whiter and cooler than any previous technology—at the greatest value. Large flood style pattern delivering a 6-inch pattern of light. The light can be paired with the VistOR PR and VistOR PRO models to provide a combination of lighting.



Product SpecificationsColor Temperature 4300° K pure white illuminationLight Intensity 55,000 lux at 24" (60 cm)Intensity Level 2-stage dimmingLED Life Efficiency 50,000 hour-rated LED lifeRadiated Heat Near ZeroDetachable Sterile Handle OptionalPattern Size 6.0"Reflector Diameter 5"Depth of Field (in) InfiniteVoltage 120 - 240 VAC 50/60 HzAgency Approvals UL 2601; CUL; CE Mar

Condition	New
-----------	-----

NUVO VistOR MSLED Dual Arm Lighting System



The NUVO Surgical Light provides superior performance for all aspects of surgery. The light head provides excellent shadow reduction, a high-intensity light pattern and optimal color rendition. Its unique, patented interface allows the user to fully control the light's intensity, pattern, position and orientation from the light head. Appropriate for the vast majority of general and specialized surgical procedures including Micro Surgery, Orthopedics, OB/GYN, Ophthalmology, and Cosmetic Surgery.

Features:

- Light Intensity (LUX): 120,000
- Color Temperature Adjustment: 4300°K
- Pattern Intensity: 5 Settings, 50 -100 %
- Pattern Size (Inches): 9-13
- Red Rendering Index: 95
- Controls in the Sterile Field: Yes
- LED Projected Life (hrs): 50,000

Condition	New
-----------	-----

Surgical Lighting

NUVO VistOR PRO LED Single Head Ceiling Procedure Light

Adding unsurpassed brilliance in the ORideal within the surgical environment, the VistOR Pro offers unsurpassed longevity and flexibility in a package that's incredibly affordable. The large diameter LED array provides excellent shadow control – even with just a single light head. Pressing the center button on the sterilizable handle allows 5-stage dimming plus ambient light mode (for minimally-invasive procedures) for easy management of brightness within the illumination field.



The VistOR Series offers a selection of three different-sized lights and provides unparalleled optics for ambulatory surgery centers, procedure rooms, trauma rooms, clinics and physician office.

Whether you're performing surgery, procedures or exams, it's absolutely critical that the lighting conditions be optimal so you can visualize your work with complete accuracy. Nothing less is acceptable. At Nuvo, our sole focus is making that happen.

Features Of The VistOR PRO LED Exam And Procedure Light:

Cross Focus Technology™ provides superior shadow control. On/off and intensity settings are easily controlled from the sterilizable handle. 36 high efficiency ultracool LED's provide superior illumination and unmatched shadow control. Precision-balanced and articulated arm to allow the greatest reach and range of positioning with surprising ease. Optional Battery Back-Up available for portable and ceiling/wall mounts. Complete range of mounting options: Single Ceiling, Dual Ceiling, Wall, and Portable Mounting options. We also offer our versatile ShuttleTrak™ ceiling-mounted lights for procedures that require extended overhead coverage. It comes in Single or Dual Trolley models. Proudly made in the USA. Product Specifications Color Temperature 4300° K pure white illumination Light Intensity 100,000 lux at 1 meter Intensity Settings 5 (from 50% – 100%) LED Life Efficiency 50,000 hour-rated LED life Radiated Heat Near-zero Controls on Sterile Handle Yes Voltage 100 - 240 VAC 50/60 Hz

Condition	New
-----------	-----

NUVO Vu LED Surgical Light

When you're performing surgery, it's critical that you can see what you're doing down to the most minute detail. That's why the Nuvo Vu LED Surgical Light should be at the top of the list for your hospital operating rooms. The Nuvo Vu LED Surgical Light is brighter, whiter, cooler and more capable in the operating room than any previous technology. So you get the most accurate colors with unparalleled shadow control and virtually no heat emission. All this for a value that's hard to match.



Features:

- Adjustable color temperature (3,800° to 4,800° K)
- Brilliant 160,000 lux output
- 50,000 hour-rated LED life
- Cross Focus Technology provides superior shadow control
- Near zero radiated heat
- Complete range of mounting options
- Green Mode for low-light, minimally-invasive surgery

Condition	New
-----------	-----

Nuvo Vu OR Lighting Touch Screen Tablet Control (Optional)



Touchscreen wall control enhances efficiency in the OR

Nuvo gives you the option of a touchscreen wall control. The wall control panel communicates wirelessly, managing one, two, or three lights as well as the HD Video Camera. While you focus on operating, your nurse or assistant can control light intensity, color temperature, Green mode, plus camera functions — all remotely.

Condition	New
-----------	-----

Ritter Midmark 354 Light



The Ritter 354 Procedure Spotlight has everything healthcare professionals need – but without the heat. It's an ideal task light for all medical environments where illumination is required for external examinations and procedures. And with its pure, crisp and defined light pattern, it can be adjusted to suit individual needs.

Features:

- Illumination at 24" (60.9 cm)
- Focal length: 7,000 fc (75,300 lux)
- Illumination at 36" (91.4 cm)
- Focal length: 4,000 fc (43,000 lux)
- Maximum pattern diameter at 36" (91.4 cm)
- Focal length: 10" (25.4 cm)
- Minimum pattern diameter at 36" (91.4 cm)
- Focal length: 3" (7.6 cm)

Condition	Refurbished
-----------	-------------

Did you know that we provide maintenance and service for every piece of equipment listed in this catalog?

Ritter Midmark 355 Procedure Light

The Ritter 355 Procedure Light is the perfect balance between quality and affordability. It's counterbalanced to provide precise positioning without drifting. A 4,200 K color temperature, large pattern size with an even distribution of light, and unsurpassed shadow control provides the most accurate color rendition for improved diagnosis.

Features:

- Requires less than 2 ½ pounds of force to move the light
- Fully autoclavable handle
- Truer color temperature with less heat

This price includes:

- Dual head
- Ceiling mount



Condition	Refurbished
-----------	-------------

Skytron Aurora Series

Skytron's Aurora® LED Surgical Lighting provides advanced bright, cool and color correct illumination to the surgical field and a selectable color temperature of 4,000K (soft white) or 4,500 (bright white) based upon surgeon preference for every procedure. Focusable LEDs extend illumination range while retaining maximum intensity at any distance. Typical applications include General, Orthopedic, Plastics, Vascular, Cardiovascular, OB/GYN, Trauma, Neurosurgery, and Urology.



Condition	Refurbished
-----------	-------------

Skytron ST9 Exam Lights

ST9 Exam Light provides high quality, high performance illumination within an extremely durable light weight polymer light head, creating the perfect balance between performance, maneuverability, and durability. The ST9 features advanced high intensity, shadow free optics without glare. Applications include LDR, LDRP, OB/GYN, Exam, Emergency Room & Special Procedures.



Features:

- High Intensity Illumination with up to 45,200 lux (4,200 fc) at a working distance of 42"
- Advanced VSRD (Vertically Segmented Reflector Design) Optics produce a myriad of rich and vibrant lighting segments blending into a high quality, soft white column of light
- Radically Durable, Featherlight Polymer Light head Construction provides for optimum maneuverability and durability
- Cool, Color Corrected Light provides natural Tissue rendering with a Color Rendering Index of 93
- Non-Sterile Positioning Handle is provided on the light head for positioning control

Condition	Refurbished
-----------	-------------

Skytron Stellar Series

Stellar® provides advanced lighting performance with high intensity, shadow-free and focusable optics, offering precise surgeon focus control for every procedure while retaining maximum intensity. Extended working distances, beyond one meter, are not uncommon for many surgical procedures such as OB-GYN, Orthopedics, Vascular, and Urology. Stellar's advanced VSRD Optics provide surgeons with sterile focus control. Just a 1/4 turn of the sterile center focus handle captures maximum intensity at any distance without having to physically reposition lightheads.



Condition	Refurbished
-----------	-------------

Surgical Lighting

Steris Harmony LC 500 OR Lights

The cost-effective Harmony® LC Surgical Lighting System provides bright halogen lighting and a flexible, easy-to-maneuver suspension system. The Harmony® LC Surgical Lighting system is a contemporary product recently emerging on the aftermarket and is offered refurbished on a limited basis.

Features

- Easier assessment of tissues and vessels with true color; intense, white light; and excellent shadow control
- Surgeons have the power to adjust pattern size and intensity
- Helps to keep the room cool with less heat-causing light
- Supports integration needs with in-field monitors or in-light cameras
- Wide range of configurations available to meet your needs



Steris SQ240

The Steris® SQ240™ Surgical Lighting Fixture offers a two-arm unit mounted on a central hub. The rotary suspension shaft of the surgical light head moves on tapered roller bearings at one end of the horizontal arm. The unit produces up to 10,000-foot candles and provides continuous lighting by automatically switching to an auxiliary bulb when needed. The SQ240™ is also available in Single Mount Lights.



Condition	Refurbished
-----------	-------------

Condition	Refurbished
-----------	-------------

Steris Harmony Series

Steris® offers a wide variety of high-quality surgical lighting products to meet the needs of any operating room – from large teaching hospitals to smaller Ambulatory Surgery Centers – because every patient deserves the best possible light. Harmony lights maintain the truest color for quick assessment of tissues and vessels to help surgeons provide the best possible patient outcomes. The optical coatings filter out 98% of heat-carrying infrared rays, which creates a cooler, more comfortable surgical environment.



Steris VLED OR Lights

The Harmony vLED Surgical Lighting System is a configurable, modular lighting system, providing high quality illumination for surgical and diagnostic applications.



The system can be configured to meet the needs of a wide variety of applications ranging from simple diagnostic procedures performed in outpatient surgery centers, to more complicated procedures performed in major acute care suites. Locations of use include operating rooms, exam and treatment rooms, hospital emergency departments and hospital intensive care units.

A wall-mounted panel provides a user interface for controlling light heads and camera. The wall-mounted panel monitors LED module status and can be used to control light intensity, as well as camera features such as focus, camera rotation and zoom. In addition, the light heads can be turned on and off, and their intensity can be controlled using conveniently located surgeon's controls on the light handle of each lighthead. The patented surgeon in-light controls include a one-touch command to turn all light heads on simultaneously to their previous settings.

Light heads provide cool, shadow controlled lighting. The system can be ordered with one, two or three light heads mounted to a suspension system capable of continuous 360° rotational positioning. The light heads are available in both adjustable and fixed pattern models. The illumination pattern size is adjusted by rotating the handle. Most Harmony vLED Surgical Lighting Systems can be configured (at time of order) to include one or two video monitor support arms.

Condition	Refurbished
-----------	-------------



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers

Washer Disinfectors

Surgical & Exam Lighting

Surgical & Exam Tables Stretchers
& Patient Transportation

Vital Sign & Patient Monitors

Blanket/Fluid Warmers

Scrub Sinks & Accessories

Defibrillators

Electrosurgical Units

Surgical Endoscopy

Imaging

Tourniquet Systems

Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxo MEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product Index

Product	Page
Dual VistOR PRO LED Ceiling Mount.....	1
NUVO VistOR MSLED Dual Arm Lighting System.....	1
NUVO VistOR PRO LED Single Head Ceiling Light.....	2
NUVO Vu LED Surgical Light.....	2
Nuvo VisTOR EX LED Examination Light Wall Mount.....	1
Nuvo VisTOR EX Portable LED Examination Light Floor Mount.....	1
Nuvo Vu OR Lighting Touch Screen Tablet Control (Optional).....	2
Ritter Midmark 354 Light.....	2
Ritter Midmark 355 Procedure Light.....	3
Skytron Aurora Series.....	3
Skytron ST9 Exam Lights.....	3
Skytron Stellar Series.....	3
Steris Harmony LC 500 OR Lights.....	4
Steris Harmony Series.....	4
Steris SQ240.....	4
Steris VLED OR Lights.....	4



Did you know that we provide maintenance and service
for every piece of equipment listed in this catalog?



Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's Daughters
Virginia Eye Institute



OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTRASTEAM
SUSSMAN

Washer Disinfector

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBEROPTICS
WELCH ALLYN
EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS

GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424

www.auxomedical.com
sales@auxomedical.com

 auxo MEDICAL