



auxoMEDICAL

**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

-

2020



Ventilators

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.



ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: support@auxomedical.com

Table of Contents

Ventilators..... 1-3



auxoMEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

New & Refurbished
Medical Equipment

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

ALLIED Healthcare AHP300 Ventilator

The AHP300 is a new broad-featured transport ventilator designed to meet the specific needs of modern EMS response. It features the broad set of ventilation modes of critical care ventilators and a single panel user interface that is easy to use, easy to clean rugged and visible even in bright light conditions. It has quick start programmable presets which let the institution pre-program up to five ventilation protocols into memory and recall them instantly with the press of separate rapid recall buttons. This feature allows the responder to typically begin ventilating a patient with their desired protocol in less than 30 seconds. It also features innovative ventilation characteristics like Dynamic Flow Optimization and Seamless Mode Transition. Includes Ventilator Power Cord Oxygen Hose Instruction Manual Strap



Broad Featured Capability Invasive and non-invasive support for infant (5 kg and up) to adult patients, with a flexible array of ventilation options including Pressure Control/Volume Control ventilation, SIMV/Assist Control spontaneous breath support, and stand-alone CPAP/BiLevel support.

The AHP300 offers value-added features

Condition	New
-----------	-----

Autovent 2000 Automatic Transport Ventilator, 1.5-Second Version

In some instances, product photos may vary from actual product models, sizes and/or colors listed on the product detail page. The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. 1.5 Second-Version



Condition	New
-----------	-----

Autovent 2000 Automatic Transport Ventilator, 2-Second Version

The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. second-version 2000



Control Modules: Weight: 24 oz. (680 grams). Dimensions: 6" (W) x 1.75" (D) x 3.5" (H) (15.24 cm x 4.4 cm x 8.9 cm). Supply Pressure Range: 40 - 60 psig. Operating Temperature: 0° F to 125° F. Storage Temperature: -40° F to 160° F. Frequency - AV2000: 8 to 20 BPM. Tidal Volume - AV2000: 400 to 1200 ml. Flow Rate: 16 to 48 lpm. Inspiratory Time - AV2000: approx. 1.5 or 2.0 seconds. Expiratory Time: 1.5 to 6.0 seconds I-E Ratio: 1:1 to 1:4. Case Material: Polyester. Input connection: Plated brass. Output Connection: Plated brass minute volume: 0 to 24 lpm. Gas Consumption Driving Gas: 0.5 lpm maximum. Driving Gas: 0.5 lpm maximum. Patient Valve Assembly:

Weight: 16 oz. (450 g). Flow: 0 - 48 lpm at 50 psig. Peak Inspiratory Flow: 48 lpm at an airway pressure drop of less than 2.5 H2O. Delivery Pressure: 60 ± 5 cm H2O. (Insp.) Crack Pressure: 0 to -2 cm (-0.8") H2O. Exhalation Resistance: 0 - 10 lpm: 1.5 cm H2O (max). 11 - 70 lpm: 3.8 cm H2O (ma

Condition	New
-----------	-----

Autovent 3000 Automatic Transport Ventilator, 1.5-Second Version

In some instances, product photos may vary from actual product models, sizes and/or colors listed on the product detail page. The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. 1.5



Condition	New
-----------	-----

Autovent 3000 Automatic Transport Ventilator, 2-Second Version

The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. second-version 3000



AutoVent 3000 ventilator for adult and pediatric use. Includes control module, patient valve, supply tubing, 6 foot oxygen hose and non-rebreathing valve with hose. 2 second inspiratory time. Control Modules Supply Pressure Range: 40-60 psig Operating Temperature: 0° F to 125° F Storage Temperature: -40° F to 160° F Frequency - AV3000: 8 to 28 BPM Tidal Volume - AV3000: 200 to 1200 ml Flow Rate: 16 to 48 lpm Inspiratory Time - AV3000: approx. 2.0 seconds Expiratory Time: 1.5 to 6.0 seconds I-E Ratio: 1:1 to 1:4 Weight: 24 oz. (680 grams) Dimensions: 6" W x 1-3/4" D x 3-1/2" H Case Material: polyester Input connection: plated brass Output Connection: plated brass Minute Volume: 0 to 24 lpm Gas Consumption Driving Gas: 0.5 lpm maximum Driving Gas: 0.5 lpm maximum

Condition	New
-----------	-----

Bipap Vision Ventilator System

Integrated Display Screen displays easy-to-read real-time graphics in waveform or bar scale format, providing the following measured and calculated numerical values to assist with ventilator management:

IPAP Total Breath Rate Exhaled Minute Ventilation % Patient Triggered Breaths Exhaled Tidal Volume EPAP Ti/Tot Peak Inspiratory Pressure Patient or Total Leak

Modes CPAP and S/T are available for maintaining airway patency and providing pressure support ventilation. PAV * /T (Proportional Assist Ventilation/ Timed mode) helps to promote improved patient comfort and enhanced patient/ventilator synchrony through independent adjustment of Volume and Flow assist based on the patient's disease state. Oxygen Control up to 100% Allows you to blend oxygen within the unit in concentrations ranging between 21-100% without auxiliary oxygen tubing attached to the mask.

Pressure regulation at the mask allows you to regulate pressure therapy to assure delivery where it counts (not inside the machine), and automatically adjusts for variations.

Upgradeable features provide flexibility in your purchase decision, now and in the future.



Condition	Refurbished
-----------	-------------

Carefusion Pulmonetic LTV 950 Portable Ventilator

An intuitive interface makes the LTV 950 easy to understand and operate. The unit features extensive patient monitoring capabilities and a comprehensive alarms package for optimum patient care. With a footprint and size at a fraction of comparably equipped ventilators (13.4 lbs), the LTV 950 allows for greater patient mobility.

Patient Comfort Designed to provide optimum patient comfort and care, the LTV 950 features:

- Variable Rise Time
- Variable Sensitivity Flow Trigger
- Variable Flow Termination
- Variable Time Termination
- Leak Compensation
- NPPV Mode
- High-Pressure Alarm Delay

Oxygen Bleed In Low Flow / Low Pressure Versatile Power Options For increased patient mobility and ease of transport, the LTV can be used with a variety of power sources:

- AC power
- LTV Transport Battery System
- LTV Universal Power Supply
- 3-Hour External Battery
- 9-Hour External Battery
- Automobile Lighter Adapter



Condition	Refurbished
-----------	-------------

Covidien Puritan Bennett 840 Respiratory Ventilator

The refurbished, Covidien Puritan Bennett® 840™ Ventilator System, for purchase or rental, is the flagship product in our line of critical care ventilators. It is responsive to patients and offers superior comfort, delivering sensitive, precise breaths to critically ill neonatal through adult patients.

Product Features Technological Sophistication – With its high-performance pneumatics, dual-microprocessor electronics, and DualView™ touchscreens. Upgradeable – The 840 Ventilator can be upgraded and customized with various software options to meet your clinical needs, today and in the future. Enhanced noninvasive ventilation (NIV) that improves comfort for patients and performance for clinicians. Designed for use on patients with a stable respiratory drive. Seamless Communications – The 840 Ventilator supports communication with all major patient monitoring and hospital information systems. It provides seamless electronic data transfer into a patient's medical record. When used with Clinivision® Mobile Patient Charting software, the package provides clinicians the opportunity to practice evidence-based medicine.



Condition	Refurbished
-----------	-------------

Pulmonetics LTV 1000 Ventilator

The Pulmonetics LTV 1000 is a positive pressure, invasive/noninvasive, conventional ventilator. The LTV 1000 can be used for a variety of settings such as ICU, transport, CPAP, sub-acute, and home care. The vent can be used on pediatric and adult patients, and can be used for acute to long-term care. The LTV 1000 is small and lightweight and offers continuous flow, so it is ideal for portable ventilation.

Features:

- Variable Rise Time
- Variable Sensitivity Flow Trigger
- Variable Flow Termination
- Variable Time Termination
- Leak Compensation
- NPPV Mode
- High Pressure Alarm Delay



Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Puritan Bennett 7200AE Ventilator

The Puritan Bennett 7200 respiratory ventilator has an advanced micro-processor design that offers CMV, assist/control, SIMV, CPAP, PEEP, and sigh controls. This critical care ventilator has a built-in nebulizer, air/oxygen blender, and variable pressure waveforms. Data is clear on the digital display and easy to input with the keyboard control.



"Product Features

- State-of-the-art gas delivery system.
- Demonstrated track record of reliability
- Full-featured microprocessor-controlled ventilator with 7202 display screen.
- Gas delivery features - Flow-by, Pressure Support, Pressure Control Ventilation
- Patient monitoring features - Graphics 2.0, Respiratory Mechanics, Auto-PEEP
- Communications features - RS-232, Digital Communications Interface (DCI 2.0)
- Optional 7250 Metabolic Monitor
- Integrated safety systems

Condition	Refurbished
-----------	-------------

Respironics PLV-102 Ventilator

Respironics LifeCare PLV-102 Ventilator (with internal oxygen blender) - The LIFE CARE PLV-102 is a microprocessor controlled, electrically powered portable volume ventilator that may also be pressure limited. Specifically designed for long-term use in the home, it is equally suitable for transport, short-term and institutional use. The ventilator is small and lightweight and will operate using a 12 volt DC external battery as well as 120 volt AC power. Reliability, versatility, and ease of operation are all hallmarks of this machine, making it the ideal choice for this setting.



The ability of the PLV-102 to ventilate a wide range of patients makes it the clinicians' choice for portable ventilation. Low flow rates (10 L/min) and low tidal volume (50 ml) allow for ventilation of pediatric patients, too. Additionally, adult patients with increased minute ventilation demands can be effectively ventilated with peak flow availability up to 120 L/min.

- Digital Displays:
 - Tidal Volume
 - I:E Ratio
 - Respiratory Rate
 - Inspiratory Flow Rate

Condition	Refurbished
-----------	-------------



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers
Washer Disinfectors
Surgical & Exam Lighting
Surgical & Exam Tables Stretchers
& Patient Transportation
Vital Sign & Patient Monitors
Blanket/Fluid Warmers
Scrub Sinks & Accessories
Defibrillators
Electrosurgical Units
Surgical Endoscopy
Imaging
Tourniquet Systems
Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product	Page
ALLIED Healthcare AHP300 Ventilator.....	1
Autovent 2000 Automatic Transport Ventilator, 1.5-Second Version.....	1
Autovent 2000 Automatic Transport Ventilator, 2-Second Version.....	1
Autovent 3000 Automatic Transport Ventilator, 1.5-Second Version.....	1
Autovent 3000 Automatic Transport Ventilator, 2-Second Version.....	1
Bipap Vision Ventilator System.....	2
Carefusion Pulmonetic LTV 950 Portable Ventilator.....	2
Covidien Puritan Bennett 840 Respiratory Ventilator.....	2
Pulmonetics LTV 1000 Ventilator.....	2
Puritan Bennett 7200AE Ventilator.....	3
Respironics PLV-102 Ventilator.....	3

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER

Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's
Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions

Surgery Center Services of America
Chesapeake Healthcare Planning
Sodexo USA



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's
Daughters
Virginia Eye Institute

OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTROSTEAM
SUSSMAN

Washer Disinfectant

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPERTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBROPTICS
WELCH ALLYN

EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424
www.auxomedical.com
sales@auxomedical.com

 **auxo**MEDICAL