



auxoMEDICAL



**New &
Refurbished**

Medical
Equipment

Product
Catalog

2019

-

2020

Warming Cabinets

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

We are pleased to announce that we can now service, repair and maintain your **anesthesia equipment!**

We are happy to provide you with solutions to best support your needs regardless of if they are biomedical, sterilization and now anesthesia.



ANESTHESIA

Biomedical & Sterile Processing

At Auxo Medical, we're always looking to expand our service and preventative maintenance capabilities to better serve you and keep your patients safe.

Biomedical Equipment

We provide a full line of Biomedical preventative maintenance, repairs, and calibrations onsite to keep you operating under NFPA99 standards and AAMI recommendations.

Sterile Processing

Our Specialty is sterilization and keeping your facility operating under joint commission standards.

We offer a variety of preventative maintenance programs that keeps your capital equipment running with little to zero down time.

We have the expertise to provide the equipment you need, the reliability to build a lasting relationship, and the responsiveness and services you need and expect.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: support@auxomedical.com

Table of Contents

Warming Cabinets..... 1-3



auxoMEDICAL

2019-2020 CATALOG

**This catalog is
a list of our
most popular
equipment.**

**Please contact
us for a
comprehensive
list of all
products.**

New & Refurbished
Medical Equipment

Order Now! Call 1-888-728-8448 Online www.auxomedical.com

Amsco QDJ01 Tabletop Warming cabinet

18 Inches Deep x 30 Inches Wide / Countertop Single Compartment Warming Cabinet

STERIS's Amsco Warming Cabinets are designed to be more than just blanket warmers; they also store and warm sterile intravenous fluids (IV fluids) and surgical irrigation solutions, as well as other linens, to programmed digital temperature settings for patient applications.



Condition	Refurbished
-----------	-------------

Blickman 7927TG Warming Cabinet

Models 7927TG (glass door) warming cabinets (blanket warmers) are a compact size for when space is at a premium. These single compartment countertop blanket warmers come standard with self-closing magnetic gasket, adjustable tension doors. An adjustable shelf is available as an option. It measures 20 - 5/8" deep, allowing for placement on most counter or tabletop surfaces. Unit is fully self contained with no heat vents, as are all Blickman blanket warmers. Temperature control from 90°F - 160°F. To set up a solutions only warming cabinet, the unit can be easily programmed by the user to heat only between 90°F - 110°F



Additional Information Model Number 7927TG Part Number 14B7927243 Exterior Dims: Width 24" (60.96cm) Depth 20-5/8" (52.39cm) Height 24-1/2" (62.23cm) Interior: CUBIC CAPACITY 2.87 CU/FT Shelves 0 Amps 3.5 Total Wattage 350 Ship Wt. (lb) est. 185

Condition	Refurbished
-----------	-------------

Blickman Warming Cabinet models 7921TG and 7921T

Models 7921TG and 7921TS single door warming cabinets are 74-1/2" tall and ideal for warming blankets or solutions. An open and adjustable shelf design adds to flexibility and expanded use. Easy-to-use controls for setting internal temperature are conveniently located above the door. Each unit comes standard with integrated casters for ease of mobility.



Features:

- Proportional-Integral-Derivative (PID) digital controllers feature tactile membrane (stainless-steel dome), full numeric keypads, over-temp audible and visual warnings as well as automatic cut-off
- Built-in temperature data recorder with real-time clock
- 160 x 180 pixel graphical display with dual color (red/green)
- Provides at a glance trend view, showing a digital representation of recent temperature history
- Self-closing, magnetic-gasket, adjustable tension doors
- Temperature control range allows for setting from 90° to 160°.

Condition	Refurbished
-----------	-------------

Blickmann Digital Full-size Warmer

Model 7924D (solid door) double door warming cabinets (blanket warmers) are 74-1/2" tall and are the perfect answer for warming blankets and solutions in two separately heated compartments. Set temperature control via key and pass-code lockout. Self-closing, magnetic gasket, adjustable tension doors. No set temperature restriction between compartments. Each warming cabinet unit comes standard with integrated casters for ease of mobility.



Temp Range 90-160 deg. F

Width 30" (76.20cm) Depth 26-5/8" (67.63cm) Height 74-1/2" (189.23cm) Shelves 2 Amps 15 Total Wattage 1500

Condition	Refurbished
-----------	-------------

Blickmann Digital Tabletop Warmer

The Blickman Warmer Cabinet Model 7922 is a Single-door warming cabinet and provide ideal flexibility for warming solutions and blankets in heated Shelves

- Open, adjustable shelf design adds convenience and flexibility.
- External Dimensions : 24" W X 36" H X 30" L - Internal Dimensions : 22H x 27L x 22 W inches - Weight : 140 Lbs Voltage : 110/120V , 50-60 Hz, 500W - 90-160 degrees F - Digital Read-out and control



Condition	Refurbished
-----------	-------------

CMP Blanket and Solution Warming Cabinet

CMP has been building high quality warming cabinets for over 65 years and has become the industry standard for hospital and emergency room blanket and solution storage. CMP stainless steel Warming Cabinets can be found in almost every hospital in the U.S. and can be configured for your unique requirements.

Models: DSW1AE-18, DSW1AE-18 (74H), DSW1AE-24, DSW1AE-24 (74H), DSW2AE-1-18, DSW2AE-1-24, DSW2AE-2-18, DSW2AE-2-24, DSW2AE-2-24, DSW3AE-18 (18W), DSW3AE-24

Features:

- Stainless Steel Construction
- Electric Heated Compartments
- Over Temperature Alarm
- Separately Controlled Chambers
- Temperature Adjusting Lockout
- Optional Quick-Temp USB temperature monitoring system
- Many Standard Sizes Available
- Microprocessor Control System
- Customizable Temperature Settings/Limits
- Wire Shelves/Baskets to maximize airflow



Condition	New
-----------	-----

Warming Cabinets

MAC Medical Warming Cabinets

These warmers are equipped with independent, digitally controlled heating chambers that offer actual temperature and set point displays. Simple plug and play data download requires no additional software.

Features:

- Automatically records data every 30 minutes
- Internal storage up to 2 years
- Displays actual temperature and set points at all times
- Independent digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Fully insulated to provide uniform heating
- Audible & visual over-temp alarm
- Remote alarm output



Condition	New
-----------	-----

Olympic Warmette 56920 Blanket Warmer

The Olympic Warmette® Model 20 Warming Cabinet features a compact design which allows medical facilities to have warm blankets and solutions closer to where they are needed. The Warming Cabinet can be placed in or near rooms that are in need of warm blankets and solutions thereby eliminating the process of having to retrieve them from an inconvenient, centrally-located warmer. Convenient placement of the Blanket Warming Cabinet allows professionals to provide solutions and blankets to patients that are closer to the desired storage temperature.



The Hospital Blanket Warmer features simple controls which allow staff members to easily set their desired temperature. The temperature setting knob can be locked at five set temperature points: 100°F, 110°F, 120°F, 130°F, and 140°F. The Medical Blanket Warmer has a lockable temperature setting which allows professionals to keep the cabinet at their desired temperature without the worry of having the internal temperature tampered with. A signal light notifies staff members when the interior temperature has exceeded the user's desired temperature. The radiant heating element on the Warmer Cabinet allows fluids and blankets to be heated uniformly.

Thanks to the transparent doors on the Nursery Warming Cabinets, staff members can keep track of inventory without having to open the cabinet which would result in unwanted temperature drops. The Blanket Warmer can be placed on a counter, mounted to a wall, or placed on top of a floor stand and features large comfortable handles that swing fully open for easy loading and unloading. The Blanket and Fluid Warming Cabinet operates silently which makes it ideal for nurseries, birthing rooms, ER, OR, and recovery rooms.

Condition	Refurbished
-----------	-------------

Pedigo 2120 Fluid Warmer Refurbished with Caster Base

Deluxe fluid warmers feature 4 digit LED display, window doors, control lock programming and auto time programming

Overall Size: 20.75 inW x 28.5 in D x 33.2 in H Interior - 3.4 cubic feet - 14.125 in W x 21.625 in D x 20 in H

Two Epoxy coated Baskets with capacity of 14L per basket includes 5" Caster Base Designed specifically for the warming of irrigation and intravenous fluids. Irrigation temperature adjusted from: 98°F to 150°F (37°C to 65.5°C) Intravenous temperature range: 98°F to 104°F (37°C to 40°C) Safety override feature prevents overheating Lockout feature Individual temperature controls for irrigation fluids and intravenous fluids allow for precise temp settings for the proper fluid type



Condition	Refurbished
-----------	-------------

Single Chamber Blanket/Fluid Warmer

Blanket/Fluid Warmer -Single Chamber -Overall size: 30W x 26.6D x 36H -Chamber size: 26W x 23D x 25H -2 adjustable shelves -8.63 cubic feet of storage -Temp range 90-160 F.



Condition	Refurbished
-----------	-------------

Skytron J2 Series Warming Cabinets

Skytron's J2 Series Warming Cabinets are designed to maintain clinically ideal temperatures for solution bottles, solution bags and blankets. Six models are available to provide optimum clinical storage requirements.

Counter Top Models: SS2201-J2, SS2201-J2G Free Standing Models: SS2205-J2, SS2205-J2G, SS2207-J2, SS2207J2G

Features:

- Double wall stainless steel construction with number 4 finish
- Field reversible left or right hinged doors
- Optional double paned, tempered glass doors
- Heating compartments fully insulated with circulating fan
- Two or three adjustable stainless shelves
- Digital temperature controls for individual compartments
- Recirculating fan in each compartment
- Adjustable temperature range



Condition	Refurbished
-----------	-------------

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Steris AMSCO QDJ05 Single Compartment Warming Cabinet

The Steris Amsco Warming Cabinet is the solution for all warming and storage needs within the operating room. This Steris Cabinet is designed as a blanket warmer as well as a warmer and storage container for sterile IV fluids and irrigation supplies. The newly designed system keeps linens and fluids warm, secure and compliant at the proper temperature settings. The Steris QDJ-05 Single Compartment Warming Cabinet is equipped with digital temperature settings to ensure proper the proper warming temperature.



Features:

- Door locks (electronic and mechanical) for security
- Double wall door construction with 1.5" insulation
- LED display shows the temperature instantly
- Temperature selection range is 90-160 degrees F
- Temperature lockout for control and heat consistency
- Convenient Celsius and Fahrenheit temperature selection
- Ability to warm and maintain temperature of linens, blankets, IV fluid and irrigation solutions
- High quality stainless steel to withstand heavy and continuous use.
- Internal storage shelves are provided for organization.
- Outside dimensions (H x W x D) at 24" Depth: 24 1/4 x 30 x 26 1/2" (616 x 762 x 673 mm)
- Inside Dimensions (H x W x D) at 24" Depth: 13 1/2 x 24 x 22 7/8" (343 x 610 x 581 mm)

Condition	Refurbished
-----------	-------------

Steris AMSCO QDJ06 Dual Compartment Warming Cabinets

The Steris Amsco Warming Cabinet is the solution for all warming and storage needs within the operating room. This Steris Cabinet is designed as a blanket warmer as well as a warmer and storage container for sterile IV fluids and irrigation supplies. The newly designed system keeps linens and fluids warm, secure and compliant at the proper temperature settings.



Features:

- LED readout for instant temperature confirmation
- Adjustable shelving for flexible organization and storage
- Reliable heating with uniform and efficient heat transfer
- Visual access to inventory with dual-pane tempered glass or traditional solid steel doors
- Temperature lock-out for consistency and control
- Celsius and Fahrenheit selection for convenience and accuracy
- Free-standing or recessed installation
- Steel or dual-paned tempered glass doors
- 18" (457mm) or 24" (610 mm) depth

Condition	Refurbished
-----------	-------------

Steris Analog Tabletop Warmer

The Steris Amsco M70 Series line of blanket warmers was one of the most popular lines ever manufactured. This series features a tabletop version with a rotary knob that controls temperature. The Amsco® M70® also features over-temperature alarms (visual and audible) and a heating indicator (visible). On full size models, the upper compartment of the Amsco® M70® is the same as the tabletop version, but with a separate lower cabinet that is larger in size than the upper. The lower cabinet of the full size model comes as either a non-heated storage compartment or a heated version, which would maintain a constant 120° F temperature for warming of blankets only.



Condition	Refurbished
-----------	-------------

Steris QDJ03

The Steris / Amsco QDJ-03 series of tabletop blanket warmers is a perfect solution to save space in your department or surgery center. The QDJ-03 is a cost-saving solution for blanket warmer needs in your medical facility versus the cost of a comparable new unit. Comes with a rack and is completely refurbished to excellent operational and cosmetic condition.



Steris Model QDJ03-Easy grip door handle-Temperature set dial, °F-LCD temperature window-Overtemp warning light-4.5 foot heavy duty power cord-21 x 23 in. shelf included-power on/off rocker switch-Fully Refurbished and Biomedically Inspected

Inner Dimensions: 13-1/2 x 24 x 22-7/8" Outer Dimensions: 29.25"W x 26.5"Deep

Condition	Refurbished
-----------	-------------

Steris QDJ04 Dual Compartment Blanket Warmer

Steris Model QDJ04-2 Compartments-Independent temperature controls for top & bottom-Contemporary design-Stainless steel construction-New circulation fan-Digital displays-Overheat indicators-Completely Remanufactured and Biomedically Inspected



Inner dimensions of the upper compartment:• width: 23.5 in. • depth: 22 • height: 13 Inner dimensions of the bottom compartment:• width: 23.5 in. • depth: 22 • height: 38 Overall dimensions of the unit:• width: 30 in. • depth: 28 • height: 75.5

Condition	Refurbished
-----------	-------------



EXPAND

Your Medical Sales Relationship

PRODUCTS

Autoclaves & Sterilizers
Washer Disinfectors
Surgical & Exam Lighting
Surgical & Exam Tables Stretchers
& Patient Transportation
Vital Sign & Patient Monitors
Blanket/Fluid Warmers
Scrub Sinks & Accessories
Defibrillators
Electrosurgical Units
Surgical Endoscopy
Imaging
Tourniquet Systems
Infusion & Syringe Pumps

Auxo Medical is focused on building quality relationships with our clients to ensure customer satisfaction and long term success.

Every customer need is based on a unique situation. You can expect Auxo Medical to tailor a plan to support your situation every time — with fast answers, and a solution designed with equipment, installation and maintenance services that match your needs.

Our focus is on delivering the highest quality refurbished and new equipment. We provide consistent value, often saving as much as 50% off the price of new systems.

The end result you can expect will be same every time, quality equipment and service, coupled with market leading value, provided by a knowledgeable sales staff.



auxoMEDICAL

5401 Distributor Drive, Richmond VA 23225

Toll Free: 888-728-8448 | Local: 804-233-4424 | Fax: 804-233-4204

Sales: sales@auxomedical.com | Support: ops@auxomedical.com

Product	Page
Amsco QDJ01 Tabletop Warming cabinet.....	1
Blickman 7927TG Warming Cabinet.....	1
Blickman Warming Cabinet models 7921TG and 7921T.....	1
Blickmann Digital Full-size Warmer.....	1
Blickmann Digital Tabletop Warmer.....	1
CMP Blanket and Solution Warming Cabinet.....	1
MAC Medical Warming Cabinets.....	2
Olympic Warmette 56920 Blanket Warmer.....	2
Pedigo 2120 Fluid Warmer Refurbished with Caster Base.....	2
Single Chamber Blanket/Fluid Warmer.....	2
Skytron J2 Series Warming Cabinets.....	2
Steris AMSCO QDJ05 Single Compartment Warming Cabinet.....	3
Steris AMSCO QDJ06 Dual Compartment Warming Cabinets.....	3
Steris Analog Tabletop Warmer.....	3
Steris QDJ03.....	3
Steris QDJ04 Dual Compartment Blanket Warmer.....	3

Did you know that we provide **maintenance and service** for every piece of equipment listed in this catalog?

Auxo Medical is providing the highest quality of service and refurbished equipment and ensuring that our customers meet the highest operational standards of modern medical technology.

Here are a few of our existing regional and national customers!

CUSTOMER

Snapshot

DIRECT EQUIPMENT CUSTOMERS

VCU Health System
HCA Healthcare
National Spine and Pain Centers
Sentara Healthcare

National Institute of Health (Bethesda) King's
Daughters Children's Hospital INNOVA Hospital System
Northside Hospital (Atlanta)



EQUIPMENT DISTRIBUTORS

Mckesson Medical Surgical
Moore Medical
SpectraMed Solutions

Surgery Center Services of America
Chesapeake Healthcare Planning
Sodexo USA



STERILIZATION/BIOMEDICAL SERVICE CUSTOMERS

HCA Healthcare
VCU Health System
Bon Secours Health System
Old Dominion University Biology Dept
VCU Massey Cancer Center
Virginia Eye Institute.
College of William and Mary
Children's Hospital of King's
Daughters
Virginia Eye Institute

OrthoVirginia
Novozymes Biological
Aramark
Blue Ridge Orthopedics
Virginia Biotechnology Research Park
Hess Orthopedics
Tidewater Physicians Multispecialty Group
SurgCenter Development
McGuire VA Hospital/Dept. of Vet. Affairs



EQUIPMENT WE SERVICE

Autoclaves

AMSCO
COX
MIDMARK
PELTON & CRANE
SCICAN
TUTTNAUER

Sterilizer

PRIMUS
SKYTRON
STERIS
GETINGE
BETASTAR
CONTINENTAL

Boilers

REIMERS ELECTRASTEAM
STERIS
CHROMALOX
PACIFIC STEAM
ELECTROSTEAM
SUSSMAN

Washer Disinfectant

STERIS
GETINGE
STEELCO
BELIMED
MIELE
SCICAN
SCIENTEK
SKYTRON

Surgical Tables

SURGICAL TABLES INC.
MIDMARK
SKYTRON
STERIS
MAQUET
STRYKER
OSI

Exam Tables

UMF MEDICAL
MIDMARK
BREWER
CLINTON INDUSTRIES

Imaging Tables

SURGICAL TABLES INC.
OAKWORKS
MEDSTONE ARCOMA
BIODEX
STERIS
OSI

Stretcher

HILL-ROM
MAC MEDICAL
MIDMARK
STRYKER
HAUSTED
GENDRON

Hospital Beds

HILL-ROM
STRYKER

Infusion Pumps

BAXTER
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Syringe Pumps

BAXTER
MEDFUSION
B. BRAUN
ALARIS
ABBOTT
HOSPIRA

Patient Monitors

BIOLIGHT
CRITICARE
MINDRAY/DATASCOPE
GE/DATEX-OHMEDA
EDAN
PHILIPS
WELCH ALLYN

Anesthesia Monitors

BIOLIGHT
MINDRAY/DATASCOPE
GE/DATEX OHMEDA
EDAN
PHILIPS

Ultrasonic Cleaners

STEELCO
STERIS
CREST ULTRASONICS
DURASONIC
SHARPERTEK
SONIC SYSTEMS INC

Surgical lights

BURTON
NUVO
SKYTRON
STERIS
MAQUET
PHILIPS
BOVIE
MEDICAL ILLUMINATION

Exam Lights

BOVIE MEDICAL CORP
SKYTRON
STERIS
MIDMARK
NUVO
PHILIPS

Anesthesia Machines

GE

Ventilators

Surgical Headlights

LUXTEC
GULF FIBROPTICS
WELCH ALLYN

EKG/ECG Machines

BURDICK
EDAN DIAGNOSTICS
GE
WELCH ALLYN
ELI

Scrub Sinks

MAC MEDICAL
STERIS AMSCO

Defibrillators

PHILIPS
ZOLL MEDICAL
MEDTRONIC

AEDs

CARDIAC SCIENCE
DEFIBTECH
PHILIPS
ZOLL MEDICAL
MEDTRONIC

Incubators

HILL-ROM
STRYKER

Fluid & Blanket

Warmer

CONTINENTAL METAL
MAC MEDICAL
SKYTRON
STERIS
BLICKMANN
OLYMPIC
PEDIGO

Electrosurgical Units

AARON BOVIE
BUFFALO
CONMED
COVIDIAN VALLEYLAB
ERBE
ETHICON
MEDTRONIC

Light Sources

LUXTEC
STRYKER
CONNED
WOLF
OLYMPUS
SMITH & NEPHEW

Freezers

Centrifuges

Warming Baths

BASIC SERVICES PROVIDED

- Equipment
- Inventories
- Equipment
- Acquisition &
- Liquidation
- Electrical Safety
- Inspections
- Equipment
- Repair
- Incoming
- Inspections
- 24 Hour support
- Periodic
- Inspections &
- Maintenance
- Equipment
- Calibrations

CONTACT

Tanner Seese
Director of Sales

804-233-4424
www.auxomedical.com
sales@auxomedical.com

 **auxo**MEDICAL