


# Fixed and Adjustable Height Processing Sinks

Ergonomically  
*Friendly*



Electronic adjustable height sink shown  
with optional lower shelves, pegboard  
backframe, and accessories

## NEW FEATURES

- 12" Angled Backsplash
- Backsplash Mounted Faucets
- 30" Deep Standard Sink Top
- 20" Deep Standard Basin
- 2 1/2" Front of Sink to Basin
- Individual Plumbing for Adjustable Height Sinks
- Universal Workflow
- Larger Seismic Brackets
- Pegboard Backframe and Accessories Options
- Bolt-on Side Splash Option
- H-Brace with Lower Shelf Option



Note: ETL compliance  
for adjustable height  
models only

[View more information on the back >>](#)



**auxo**MEDICAL®

**MADE  
IN THE USA**

## STANDARD FEATURES

- Open pipe base
- Single, dual, or triple basins
- Fully welded construction with reinforced tops
- 37" standard height for stationary sinks
- Adjustable feet
- Marine edge
- Flat working surface
- Etched gallon markings
- Basin size: 24"W x 20"L x 10"D

## OPTIONAL FEATURES

- **Ergonomically friendly** electronic adjustable height with four position programming
  - Electronic adjustable height range: 35"-42"
- Overhead LED Light
- Mixing valve faucet
- Pre-rinse faucet
- Deionized (DI) water faucet
- Air gun kit (with 5 tips)
- Water gun kit (with 8 tips)
- Dosing soap pumps
- Eye wash
- Larger basin: 30"W x 20"L x 10"D
- Hand lever drains
- Foot operated drains (stationary height models only)
- Copper drain system (not compatible with deionized water)

## Fixed Height Processing Sinks

Base Model #	Dimensions
PS307237-1000	30"D x 72"W x 37"H, Single Basin
PS309637-1000	30"D x 96"W x 37"H, Dual Basin
PS3012037-1000	30"D x 120"W x 37"H, Triple Basin

*Note: Model # will change based on features added*

## Electronic Adjustable Height Processing Sinks

Base Model #	Dimensions
PS307235-1000	30"D x 72"W x 35"H, Single Basin
PS309635-1000	30"D x 96"W x 35"H, Dual Basin
PS3012035-1000	30"D x 120"W x 35"H, Triple Basin

### **Power Requirements (Electric Adjustable Height Sinks)**

120 Volt, 60 Hz, single phase, 15 amp GFCI protected electrical outlet, installed per local building codes, and protective earth.

*Note: Model # will change based on features added*

To place an order, contact us at **(804) 233-4424** or **(888) 728-8448**.



5401 Distributor Drive, Richmond, VA 23225  
p. 804.233.4424 / 888.728.8448 f. 804.233.4202 w. [www.auxomedical.com](http://www.auxomedical.com)

Information regarding this product is subject to change without prior notice.