

INTRODUCING THE HSS SERIES SURGICAL STOOL



HSS500-S SADDLE SEAT STYLE



HSS500-WD WATER DROP SEAT STYLE



HSS500-W WEDGE SEAT STYLE

BUILT FOR PERFORMANCE
AND COMFORT

HAUSTED SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES



Adjustable Contoured Backrest

300 POUNDS
300 lb (136 kg)
Weight Capacity
EVENLY DISTRIBUTED

360° Rotating
Urethane Arm Rest
for Comfortable
Arm Positioning

Dual-density Foam Seat
with **DURA-GEL**[®] Insert
for Enhanced Comfort



180° Arm Rest Support for more
Stable Arm Positioning

Adjustable Height Range
21 - 28 inches (53 - 71 cm)

OPTIONAL ACCESSORY:
Riser Kit Increases Height Range
an additional 2 inches (5 cm)
Item No: H081025

ABS Thermoformed
Hood for Easy Cleaning



Complete the HSS Series Surgical Stool
Experience with the new Hausted
ESC² Ophthalmology Platform
Item No: ESC2EYEST



FOOT-CONTROLLED HYDRAULICS
with Three-Spindle Base and
Four 3 inch (8 cm) Dual Tente[®]
Casters for Greater Stability



Front
Locking
System

ENGINEERED TO PROVIDE
ENHANCED
SURGEON COMFORT



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



EC REP

CE Partner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu

graham
field