



# auxoMEDICAL<sup>®</sup>

## Surgical Tables Inc. **SpineXR Table**



Easy C-arm access

Motorized rotation

Unmatched stability

Lightweight and easy to move

Motorized table length

Made in the USA

An Affordable and  
Versatile Option



## Standard Accessories

- Hand Control
- Patient Straps

## Optional Accessories

- Carbon Fiber Top
- Electric Floor Locks
- Cervical Traction



Adjustable Prone Patient Supports



Adjustable Foot Rest



Adjustable Headrest



6" Side Rail



Adjustable Arm Board



Disposable Headrest Comfort Liners



Heel Supports



Shin Supports



Disposable Pad Covers  
(standard or comfort)



Shown with optional accessories.

## Features

- **Motorized adjustable table length 88"-102"** (fits 4'8" to 6'6" patient)
- **Carbon fiber patient support rails** (compatible with other brand table accessories)
- **4 Swivel, total-lock castors, 5"**
- **Audible message on power-up** indicating table is ready to use and to lock castors
- **Hand control for Height, Trendelenburg and Lateral Rotation**
- **Table length adjustment on safety-covered switch** located on table base
- **Emergency stop**
- **Removable power cord**
- **Illuminated power switch**

## Specifications

- **500 lb Weight Limit**
- **Table Dimensions:** 88"-102"L x 29.5"W
- **Power:** 120V/15A/60Hz
- **Height Adjustment:** 30.5"-50.5"
- **Trendelenburg:** +/- 12°
- **Table Weight:** 880 lbs

