TABLE OF CONTENTS

AEDs 2
- Cardiac Science Powerheart G5, Semi-Auto, Dual Language 2
- Philips Heartstart FRX AED w/ Wall Mount Box 3
- Philips OnSite AED 4
- AED Wall Cabinet Small with Standard Audible Alarm 4
- Cardiac Science Powerheart G5 New, ICPR, Semi-Auto, Dual Language 5
- Fully Automatic Zoll AED Pro 6
- Powerheart AED G3 Plus 6
- Philips Heart Start FRx AED 7
- Heartsine Samaritan PAD 350P AED 8
- Zoll AED Plus 8
- Cardiac Science Powerheart AED G3 8
- Philips HeartStart Onsite 9



AEDS

When automated external defibrillators are reliable and easy to operate, they can help save victims of sudden cardiac arrest. Auxo Medical knows that during a crisis, dependable and patient-friendly AEDs need to be rescue ready.

We carry the top brands of automated external defibrillators, including HeartSine and Philips. All of Auxo Medical's refurbished AED units come with easy-to-follow instructions. They are ready for years of reliable service.

Easy-to-understand visual and voice prompts guide users during chaotic situations. Auxo Medical provides top of the line units that prompt CPR after an initial shock, to make sure life saving measures are efficient and timely. We have variety of AED units for you to choose from.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED





CARDIAC SCIENCE POWERHEART G5, SEMI-AUTO, DUAL LANGUAGE

The **Cardiac Science Powerheart G5 Semi-Auto AED with Dual Language capability** is a state-of-the-art emergency device designed to maximize survival chances during sudden cardiac arrest. Easily switch between English and Spanish at the push of a button, this AED features RESCUE READY® technology, conducting daily self-checks of all critical components and regular high-voltage system charges to ensure optimal readiness. The device provides clear visual and voice prompts, along with RescueCoach® CPR prompting and a metronome to guide rescuers through effective chest compressions. Its non-polarized electrode pads simplify placement, reducing steps during emergencies. Equipped with an IntelliSense battery that monitors performance and storage conditions, the G5 delivers a shock in less than 10 seconds after CPR, making it a reliable, user-friendly solution for both professional and lay responders seeking to provide immediate, life-saving assistance.
Read More
SKU: AM-CSG5S-80C-S-R
Price: Request Quote



PHILIPS HEARTSTART FRX AED W/ WALL MOUNT BOX

The **Philips HeartStart FRx AED** is a state-of-the-art defibrillator, notable for its durability and user-friendly design. Backed by a 1-year manufacturer's warranty, this model uniquely offers an optional infant/child key, enabling the use of adult pads on pediatric patients. Lightweight and rugged, the FRx is designed to be both portable and resilient. This package includes the FRx Biphasic unit, a new long-life lithium battery, Smart Pads II electrodes, a user's manual, and a convenient carry case. Also included are biocertification and a physician's prescription, as required by the FDA. **A wall mount box is provided for easy and accessible storage**.

Read More
SKU: AM-PHSAED-FRX
Price: Request Quote



PHILIPS ONSITE AED

The **Philips OnSite AED** is designed for straightforward operation, prioritizing ease of use in emergency situations. Activating the device is as simple as pulling the handle, which initiates guided audio instructions to assist users through the process. Utilizing SMART Biphasic Therapy, the AED delivers a safe and effective shock. A built-in self-test ensures the device is always ready for use, and the Quick Shock feature minimizes interruptions to chest compressions. Integrated SMART Pads sense when they are properly applied to the patient's skin. The AED also features clear, intuitive visual cues with icon-driven operation. Each unit includes a new battery, adult pads, a carrying case, and a user manual to ensure preparedness.

Read More
SKU: AM-PHONSITE-AED
Price: Request Quote



AED WALL CABINET SMALL WITH STANDARD AUDIBLE ALARM

The **AED Wall Cabinet Small with Standard Audible Alarm**, designed to ensure immediate accessibility and security for your AED device. Measuring 16 inches in length, 8 3/8 inches in depth, and 14 5/8 inches in height, our compact cabinet features a robust 20-gauge cold rolled steel construction with a durable powder coat paint finish. The clear acrylic window and magnetic door latch provide easy visibility and access, while the prominent graphics indicate the presence of a defibrillator and a warning for the alarmed system. Our standard audible alarm, which can be set with a keyed system, emits a high decibel sound ranging from 80 to 120 dB when the cabinet door is opened, effectively deterring theft and alerting those nearby to a potential emergency. Importantly, **the cabinet accommodates a variety of AED models, including Philips, Heartsine, Cardiac Science, and Zoll**, ensuring versatility and readiness in critical situations. As an industry standard, our cabinets are designed to remain unlockable, ensuring that life-saving AEDs are always accessible when needed.

Read More
SKU: AM-HST-CAB-A
Price: Request Quote



CARDIAC SCIENCE POWERHEART G5 NEW, ICPR, SEMI-AUTO, DUAL LANGUAGE

The **Powerheart G5** is Cardiac Science's latest addition to their line of automated external defibrillators (AED or AED Machine). The G5 is loaded with capabilities to allow professional and lay responders to give a patient their best chance of surviving a sudden cardiac arrest event.

Read More
SKU: AM-CSG5S-80C-S
Price: Request Quote





FULLY AUTOMATIC ZOLL AED PRO

ZOLL Fully Automatic AED Pro with AED Cover, CPR-D-padz® Electrode, ZOL8000-0860-01 AED Pro Non-Rechargeable Lithium Battery, and Carry Case Operator's Guide Included

The ZOLL Fully Automatic AED Pro is an advanced medical device that provides life-saving assistance in emergency situations. This automated external defibrillator (AED) is designed to be user-friendly and can be used by both trained professionals and bystanders. The fully automatic feature allows the device to analyze the patient's heart rhythm and deliver a shock if necessary without any manual intervention. This makes the AED Pro an ideal tool for first responders and healthcare professionals who need a quick and reliable device to treat sudden cardiac arrests. The AED Pro also offers advanced features such as Real CPR Help, which provides real-time feedback on the quality of cardiopulmonary resuscitation (CPR) being performed. With its durability, ease of use, and advanced features, the ZOLL Fully Automatic AED Pro ensures that timely and effective treatment is provided to those in need, ultimately increasing the chances of survival for patients experiencing a cardiac emergency.

Read More
SKU: ZOL93010340499991010
Price: Request Quote






POWERHEART AED G3 PLUS

With its RHYTHMx analysis algorithm, the G3 Pro lets first responders proactively monitor for life-threatening arrhythmias. The earlier a responder recognizes cardiac arrest, the shorter the time to defibrillation, and the better the chance for survival.

Read More
SKU: N/A
Price: Request Quote





PHILIPS HEART START FRX AED

The **Philips HeartStart FRx AED** is a cutting-edge defibrillator designed for ease of use and durability. It stands out with its optional infant/child key, allowing the use of adult pads on pediatric patients, a feature unique to this model. Lightweight and rugged, the FRx is built to withstand demanding environments without compromising portability. It comes equipped with a carrying case, SMART Pads II, a lithium battery pack, and a Heart Smart CPR/AED Rescue Kit. Furthermore, the device offers lifetime tracking and notification technology, a lighted display, and an interactive user interface to guide users through the rescue process. A 1-year manufacturer's warranty provides additional peace of mind.

Read More
SKU: AM-HEARTSTART-FRX
Price: Request Quote



HEARTSINE SAMARITAN PAD 350P AED

The HeartSine Samaritan PAD is an automated external defibrillator (AED) that analyzes the heart rhythm and delivers an electrical shock to victims of Sudden Cardiac Arrest in order to restore the heart to normal rhythm. This user-friendly AED provides easy to follow visual and audio prompts to guide the rescuer from pad application to shock delivery, if required.

Download Brochure
Read More
SKU: AM-HS350P-AED
Price: Request Quote



ZOLL AED PLUS

The **Zoll AED Plus** is an advanced automated external defibrillator that stands out with its Real CPR Help feature, providing real-time feedback to guide users through effective CPR. Equipped with Biphasic Technology, the device ensures precise and effective defibrillation during life-saving interventions. The AED Plus comes complete with a durable carry case, making it easy to transport and protect. It includes one set of long-lasting lithium batteries and one set of Stat-Padz, ensuring the device is always ready for use in emergencies. To maintain optimal functionality, the AED Plus also features an AED Check Tag for regular maintenance checks. The package is rounded out with an operator’s manual, offering comprehensive instructions for both novice and experienced users.

Download Brochure
Read More
SKU: AM-ZAEDP
Price: Request Quote



CARDIAC SCIENCE POWERHEART AED G3

In the chaos that follows sudden cardiac arrest (SCA), concerned but untrained bystanders are scared. Do they know what to do? The Powerheart™ AED G3 Plus with RescueCoach™ voice prompts to talk rescuers through the steps. When the rescuer applies the defibrillation pads, the G3 Plus Automatic analyzes the heart rhythm and “knows” when to deliver (or not deliver) the shock. The shock is delivered automatically, with no button to push, and no human intervention (semi-automatic version also available). After the shock, the unit prompts for CPR with a built-in metronome that sets the pace for proper chest compressions.

Download Brochure
Read More
SKU: AM-PHG3
Price: Request Quote



PHILIPS HEARTSTART ONSITE

Philips® designed the HeartStart® OnSite Defibrillator for the ordinary person in the extraordinary moment. The HeartStart® is designed to be the easiest to set up and use and the most reliable defibrillator available. Philip’s innovative technology, based on extensive research and user feedback, has produced a defibrillator so easy to use that you can potentially save the life of a coworker, friend, or anyone else stricken with sudden cardiac arrest. Using clear, clam voice instructions, it guides you through each step of defibrillation, including CPR coaching. Integrated SMART Pads placed on the victim’s bare skin sense and adapt the defibrillator’s instruction to your actions every step of the way. The HeartStart® can be used to treat infants and children as well as adults.

Download Brochure

Read More
SKU: AM-PHONSITE-AED-N-CAB
Price: Request Quote

INDEX

A
AED Wall Cabinet Small with Standard Audible Alarm 4

C
Cardiac Science Powerheart AED G3 8
Cardiac Science Powerheart G5 New, ICPR, Semi-Auto, Dual Language 4
Cardiac Science Powerheart G5, Semi-Auto, Dual Language 2

F
Fully Automatic Zoll AED Pro 5

H
Heartsine Samaritan PAD 350P AED 7

P
Philips Heart Start FRx AED 6
Philips Heartstart FRX AED w/ Wall Mount Box 3
Philips HeartStart Onsite 8
Philips OnSite AED 3
Powerheart AED G3 Plus 6

Z
Zoll AED Plus 7