TABLE OF CONTENTS

Biomedical Parts & Accessories 3
Vetland V1400 Series Optional EtCO2 Module 3
Respironics CO2 Monitoring Module 4
Mindray V12 V21 Nellcor Superpod Module 4
Philips InvelliVue MP5 Fixed Wall Mount Kit 5
Mindray MPM Masimo Set Mortara Module SpO2/ ECG/ NIBP/ 2-IBP 2-Temp V 2.0 115-022715-00 5
GE Carescape E-MiniC - C02 Module 6
Mindray Microstream CO2 Module Passport 8 and 12 6
Mindray Sidestream CO2 Module Passport 8 and 12 7
Mindray Add-on WiFi module material package FRU epm10M / epm12M 7
Stryker 940 Cast Cutter - Stainless Steel Blade (#STB043) 8
Stryker 940 Cast Cutter - Titanium Nitrided Blade (#STB045) 8
Mindray Passport and ePM Monitor Paper Thermal - 50mmx20m - Box of 12 (OEM: 0683-00-0505-02) 9
Mindray Passport and ePM Monitor Paper Thermal - 50mmx20m - Single Roll 9
Criticare / Datascope Patient Monitor Printer Paper (5PK) 10
Edan Patient Monitor Thermal Chart Paper (10 Rolls per Box) - Edan # 01.57.78035 10
Respironics CO2 Monitoring Module Only M2741A 11
Mindray IBP Module (pn 115-047286-00) includes module 12
Dinamap Mobile Stand 12
Patient Monitor Mobile Stand - Universal 13
Adjustable Swivel and Tilting Monitor Wall Mount - Universal 14
GE E-PRESTN Carescape Multi Parameter Module 14
Biolight Monitor Printer Paper - Pack of 4 Rolls 14
Mindray ePM AG Module 15
Mindray Sidestream CO2 Module DPM6, DPM7 15
Zimmer ATS Dual Blue Hose 16
ATS Dual Red Hose 16
Cuff Cylindrical 34 Inch ATS Zimmer 17
Cuff Large 42 Inch ATS Zimmer 17



BIOMEDICAL PARTS & ACCESSORIES

The Biomedical Parts & Accessories category offers a wide selection of replacement components and add-ons that help healthcare facilities and biomedical technicians keep essential equipment running safely and efficiently. Shoppers can expect to find replacement batteries and chargers, power supplies and cords, patient monitoring accessories such as ECG leads, SpO2 sensors, and blood pressure cuffs and hoses, and consumables like electrodes, filters, tubing, and recorder paper.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED









VETLAND V1400 SERIES OPTIONAL ETCO2 MODULE

The **Vetland V1400 Series Optional $EtCO_2$** Module offers versatile capnography options, supporting either mainstream or sidestream measurement to suit various veterinary monitoring needs. Designed for seamless integration with Vetland V1400 series monitors, this module provides accurate end-tidal CO_2 readings, enhancing respiratory and anesthesia monitoring during veterinary procedures.
*Please note this listing is for 1 (one) EtCO2 Module - Mainstream OR Side Stream
Read More
SKU: V-OPT-CO2
Price: Request Quote



RESPIRONICS CO2 MONITORING MODULE

The **Respironics CO_2 Monitoring Module** is an ideal solution for non-intubated patients, featuring a proprietary sample cell that protects the internal NDIR components from contamination, ensuring long-term durability and reliable measurements. Its common connector allows easy interchangeability between mainstream and sidestream monitoring, supported by a robust, maintenance-reducing pump that requires no calibration. The system's small, lightweight LoFlo sensor can be shared across multiple rooms or during patient transport, and it can be mounted on bed rails, IV poles, or integrated into existing monitoring systems. Included with this package are one disposable adult nasal cannula, an adult/pediatric airway adapter set with dehumidification tubing, and a sampling line kit with dehumidification tubing, facilitating versatile and effective respiratory monitoring across diverse clinical settings.

Read More
SKU: AM-12.08.078137-R
Price: $1,295.00



MINDRAY V12 V21 NELLCOR SUPERPOD MODULE

The **Mindray V12 V21 Nellcor Superpod Module** is a refurbished, highly versatile monitoring accessory designed to enhance the capabilities of Mindray V Series monitors. This module provides comprehensive patient monitoring with ECG and respiration assessment, temperature measurement, dual invasive blood pressure (IBP), Nellcor SpO_2, and non-invasive blood pressure (NIBP). Compatible with V12 and V21 models, it offers a flexible, integrated solution for adult, pediatric, and neonatal care. Please note that this listing includes only the module; accessories are available separately for complete setup. OEM Part Number: 0998-00-1802-0102A.

Read More
SKU: 0998-00-1802-0102A
Price: $795.00



PHILIPS INVELLIVUE MP5 FIXED WALL MOUNT KIT

The **Philips IntelliVue MP5 Fixed Wall Mount Kit** features a sturdy 6-inch (15.2 cm) fixed arm designed for secure wall mounting of the patient monitor. It allows for tilt and rotation adjustments to optimize viewing angles and patient access. The kit includes cable hooks positioned beneath the arm for organized hanging of patient leads, enhancing clutter-free workspace and safety. Please note that the 19-inch GCX channel for mounting is sold separately, providing flexibility to customize the mounting setup according to clinical needs.

Read More
SKU: AM-PH-0061-01
Price: $195.00



MINDRAY MPM MASIMO SET MORTARA MODULE SPO2/ ECG/ NIBP/ 2-IBP 2-TEMP V 2.0 115-022715-00

The **Mindray MPM Masimo Set Mortara Module (115-022715-00)** is a versatile multi-parameter patient monitoring module that integrates SpO_2, ECG, NIBP, two invasive blood pressures (2-IBP), and two temperature channels into compatible Mindray patient monitors such as Passport 12m, Passport 17m, DPM 6+, and DPM 7+. Featuring advanced software version 2.0, this lightweight module weighs approximately 2 pounds and measures 5-1/2 inches by 3 inches by 4 inches, making it highly portable and easy to incorporate into clinical workflows. While it does not include any accessories, this module provides comprehensive vital sign monitoring to enhance patient assessment and care in various medical settings.

Read More
SKU: AM-11-022715-00-R
Price: Request Quote



GE CARESCAPE E-MINIC - C02 MODULE

The **GE Carescape E-MiniC - C02 Module** is a compact, reliable solution for precise carbon dioxide monitoring during anesthesia and critical care procedures. Designed to seamlessly integrate with GE's Carescape anesthesia systems, this module provides real-time, high-accuracy CO_2 measurements essential for patient safety. It includes a **Mini D-Fend Water Trap**, which effectively captures moisture and condensation from the breathing circuit, helping to maintain sensor accuracy and extend its lifespan. Perfect for surgical, intensive care, and respiratory therapy environments, the E-MiniC offers streamlined performance in a space-efficient design, ensuring clinicians have dependable data to optimize ventilation and patient outcomes.

Read More
SKU: AM-GE-E-MiniC
Price: Request Quote



MINDRAY MICROSTREAM CO2 MODULE PASSPORT 8 AND 12

The **Mindray Microstream CO2 Module** (Part Number 115-030771-00, Oridian) delivers accurate CO2 monitoring via the Microstream sampling technology. Compatible with a variety of anesthesia and patient monitoring systems, it provides fast CO2 response, reliable measurements, and straightforward setup. The Passport 8 and Passport 12 variants offer flexible mounting and seamless integration for diverse clinical environments. This module includes an IBP (invasive blood pressure) connection option for integrated physiology monitoring, supporting efficient workflows in perioperative and critical care settings.


Read More
SKU: AM-115-030771-00
Price: Request Quote



MINDRAY SIDESTREAM CO2 MODULE PASSPORT 8 AND 12

The **Mindray Sidestream CO2 Module** provides reliable, non-invasive CO2 monitoring by sampling exhaled breath from the patient's airway. Designed for integration with compatible anesthesia and patient monitoring systems, it offers accurate CO2 measurements, fast response times, and easy setup. Suitable for a range of clinical settings, the **Passport 8 and Passport 12 variants** are compact, robust modules that support seamless workflow, flexible mounting, and dependable performance for both routine and complex patient care.

Read More
SKU: AM-115-011185-00
Price: Request Quote



MINDRAY ADD-ON WIFI MODULE MATERIAL PACKAGE FRU EPM10M / EPM12M

WiFi module material package FRU
epm10M / epm12M

Read More
SKU: 115-059744-00
Price: Request Quote



STRYKER 940 CAST CUTTER - STAINLESS STEEL BLADE (#STB043)

The **Stryker 940 Cast Cutter Stainless Steel Blade (Replacement #STB043, OEM #940-25**) is designed for use with the Stryker 940 cast saw. This blade is specifically **optimized for plaster casts**, offering dependable performance for that material. It features a 2.36-inch oval blade and is made from stainless steel, providing solid durability and reliable cutting for plaster applications. Each replacement package includes 1 blade.

QTY: 1 Blade

Read More
SKU: 0940-025-000
Price: $114.95



STRYKER 940 CAST CUTTER - TITANIUM NITRIDED BLADE (#STB045)

The Stryker 940 Cast Cutter Titanium Nitrided Blade (Replacement #STB045, OEM# 940-035) is a precision blade designed for cutting synthetic casts. The STB045 blade fits the Stryker 940 cast saw and features a durable titanium-nitrided coating for enhanced longevity and smooth performance. It has a 2.36-inch oval blade and is **specifically optimized for synthetic cast materials**, delivering reliable cuts with improved strength and corrosion resistance.

QTY: 1 Blade

Read More
SKU: 0940-035-000
Price: $114.95





MINDRAY PASSPORT AND EPM MONITOR PAPER THERMAL - 50MMX20M - BOX OF 12 (OEM: 0683-00-0505-02)

The **Mindray Passport and ePM Monitor Paper Thermal**, measuring **50mmx20m**, is an essential accessory designed specifically for use with all Mindray monitors. This high-quality thermal paper ensures accurate and reliable recording of vital signs and other important patient data. **Each box contains 12 rolls**, providing a convenient supply for busy healthcare environments where consistent monitoring is crucial. The paper's dimensions are optimally suited for seamless compatibility with Mindray devices, ensuring smooth operation and clear, legible printouts that are vital for patient record-keeping and analysis. Whether in a hospital, clinic, or other medical settings, this thermal paper is an indispensable tool for maintaining efficient and effective patient care.

Mindray OEM Part Number: 0683-00-0505-02

Read More
SKU: 0683-00-0505-02
Price: $39.95



MINDRAY PASSPORT AND EPM MONITOR PAPER THERMAL - 50MMX20M - SINGLE ROLL

The **Mindray Passport and ePM Monitor Paper Thermal**, measuring **50mm by 20m**, is an essential accessory for healthcare professionals utilizing Mindray monitors. Designed as a **single roll**, this thermal paper ensures accurate and reliable recording of vital patient data across all Mindray monitoring systems. Its compatibility with various Mindray models makes it a versatile choice for medical facilities, ensuring seamless integration and uninterrupted functionality. The 50mm width provides ample space for detailed readings, while the 20-meter length ensures longevity and reduced frequency of replacements, streamlining operations in fast-paced healthcare environments. This thermal paper is crucial in maintaining precise patient documentation, ultimately supporting high-quality patient care.

Read More
SKU: AM-0683-00-0505-02-EA
Price: $3.95



CRITICARE / DATASCOPE PATIENT MONITOR PRINTER PAPER (5PK)

The **Criticare/Datascope Patient Monitor Printer Paper**, available in a convenient **5-pack**, is an essential accessory for medical facilities utilizing patient monitoring systems. Each roll of blank white printer paper measures 57mm in width and extends to 40 meters in length, ensuring ample space for printing vital patient data and monitoring results. This high-quality paper is designed to seamlessly integrate with Criticare and Datascope monitors, providing clear and accurate printouts that are crucial for patient record-keeping and analysis. Its durable construction ensures smooth operation without jamming, allowing healthcare professionals to focus on delivering excellent care without disruptions. Whether for routine check-ups or critical monitoring, this printer paper is a reliable choice for busy medical environments.

Read More
SKU: AM-553
Price: $54.28



EDAN PATIENT MONITOR THERMAL CHART PAPER (10 ROLLS PER BOX) - EDAN # 01.57.78035

The **Edan Patient Monitor Thermal Chart Paper**, identified by Edan #01.57.78035, is an essential accessory for healthcare facilities utilizing the EDAN IM50 and M80 patient monitors. This high-quality thermal chart paper is specifically designed to accurately capture and document vital patient data, ensuring precise and reliable record-keeping. **Each box contains 10 rolls**, making it a cost-effective solution for medical environments that require consistent and uninterrupted monitoring. The **50x20** dimensions of the chart paper are crafted to fit seamlessly into the EDAN patient monitors, allowing healthcare professionals to focus on patient care without concerns about compatibility or installation. This product is a vital component in maintaining the efficiency and accuracy of patient monitoring systems, supporting healthcare providers in delivering high-quality care.

Read More
SKU: AM-01.57.78035
Price: $60.03



RESPIRONICS CO2 MONITORING MODULE ONLY M2741A

The **Respironics CO2 Monitoring Module Only M2741A** is a state-of-the-art device that utilizes advanced infrared absorption spectroscopy to accurately measure end-tidal CO2 (EtCO2) and respiration rates. This sophisticated sensor is capable of delivering a precise capnogram for respiratory rates of up to 150 breaths per minute, making it highly effective for a wide range of patients, including infants and smaller individuals, thanks to its efficient sampling rate of 50 ml/min. Within just 20 seconds, the module provides detailed waveform data, ensuring timely and reliable monitoring. Its user-friendly design features simple snap-in sensors, and the device boasts a rugged construction with no moving parts, making it both shock-resistant and IPX4 splash-proof. Additionally, the module is engineered to require no regular calibration, enhancing its convenience and reliability in a variety of clinical settings.

Read More
SKU: AM-152.001.0001
Price: Request Quote



MINDRAY IBP MODULE (PN 115-047286-00) INCLUDES MODULE

IBP Module (pn 115-047286-00) includes module
Can be used with ePM 12M only

The IBP (Invasive Blood Pressure) Module, identified by part number 115-047286-00, is a specialized component designed for the specific application with the Mindray ePM 12M monitoring systems. This module enables the ePM 12M to monitor invasive blood pressure, offering a critical functionality for patients who require continuous blood pressure monitoring with high accuracy. The compatibility of the IBP Module is exclusive to the ePM 12M, ensuring that it seamlessly integrates and functions within this system only. This design ensures that healthcare professionals can rely on the precision and reliability of the invasive blood pressure readings, crucial for managing patients in critical or intensive care settings where monitoring blood pressure with great precision is necessary for informed clinical decisions. The integration of the IBP Module with the ePM 12M enhances the monitor's capabilities, making it a powerful tool in the hands of healthcare providers for patient care.

Read More
SKU: 121-001580-00
Price: Request Quote



DINAMAP MOBILE STAND

5 caster Height Adjustable, clamp mount, with Basket and Handle

Read More
SKU: AM-MOBSTAND-DM
Price: Request Quote



PATIENT MONITOR MOBILE STAND - UNIVERSAL

This new **Universal Patient Monitor Mobile Stand** is an essential piece of equipment designed to facilitate seamless patient monitoring. This stand comes with a Universal Plunge Plate Slide Mount, which can be specified to pair with any monitor. The installation includes a mounting plate that can be tilted to provide the optimal viewing angle for healthcare professionals.

• Basic Mobile Stand • Universal Plunge Plate Slide Mount • Specify Monitor that will be paired with it • Mounting plate tilt for optimal viewing angle • Monitor Height range 33.5-44.5" • Utility basket 13"wx8"dx5"h • Net Weight 10lbs • Wheel Base diameter 27”

Read More
SKU: AM-MOBSTAND-6001
Price: Request Quote



ADJUSTABLE SWIVEL AND TILTING MONITOR WALL MOUNT - UNIVERSAL

Universal Adjustable Swivel and Tilting Monitor Mount

Long fixed arm plus a short swivel arm. Swivel and tilting head, basket included. Standard 13 inch rail, aluminum matte finish on rail. White paint on arm. Channel underneath arm.

Includes: Basket

Read More
SKU: AM-JW6002FW
Price: Request Quote



GE E-PRESTN CARESCAPE MULTI PARAMETER MODULE

GE E-PRESTN Carescape Module Refurbished. This Multi Parameter Module comes with ECG up to 12 leads, 2 IBP Invasive Blood Pressure, SPO2, NIBP, RESP and Temperature.

OEM Part Number:
E-PRESTN

Read More
SKU: AM-E-PRESTN
Price: Request Quote



BIOLIGHT MONITOR PRINTER PAPER - PACK OF 4 ROLLS

Read More
SKU: TPP-50
Price: $16.95



MINDRAY EPM AG MODULE

The ePM 12M paired with the ePM multi-gas module, automatically identifies and measures the inspired and expired values of O2, CO2, N2O and 5 anesthetic agents providing a compact, yet powerful OR monitoring solution.

Read More
SKU: 121-001936-00
Price: Request Quote



MINDRAY SIDESTREAM CO2 MODULE DPM6, DPM7

The **Mindray Sidestream CO_2 Module for DPM6 and DPM7 Monitors** provides reliable, accurate monitoring of exhaled carbon dioxide levels using sidestream technology. Designed specifically for compatibility with Mindray DPM6 and DPM7 patient monitors, this module enables continuous ventilation assessment crucial for anesthesia, critical care, and emergency applications. Its compact, easy-to-install design ensures seamless integration into existing monitoring setups, supporting clinical decision-making with real-time capnography data.

Warranty: 1 Year

Read More
SKU: AM-6800-30-50137
Price: $1,995.00



ZIMMER ATS DUAL BLUE HOSE

Dual Channel Blue Hose

Read More
SKU: 60-4017-002QB
Price: Request Quote



ATS DUAL RED HOSE

Zimmer Dual Channel Red Hose

Read More
SKU: 60-4017-002QR
Price: Request Quote



CUFF CYLINDRICAL 34 INCH ATS ZIMMER

34 INCH DPSB ATS CYLINDRICAL CUFF AUTOMATIC TOURNIQUET SYSTEMS (A.T.S.) REUSABLE CYLINDRICAL CUFFS (NON- STERILE)

Read More
SKU: 60-7500-006-00
Price: Request Quote



CUFF LARGE 42 INCH ATS ZIMMER

DPSB CONTOUR CUFF - XL (30 IN-39 IN) AUTOMATIC TOURNIQUET SYSTEMS (A.T.S.) XL CONTOUR CUFFS (REUSABLE, NON- STERILE)

Read More
SKU: 60-7500-008-00
Price: Request Quote

INDEX

A
Adjustable Swivel and Tilting Monitor Wall Mount - Universal 12
ATS Dual Red Hose 15

B
Biolight Monitor Printer Paper - Pack of 4 Rolls 13

C
Criticare / Datascope Patient Monitor Printer Paper (5PK) 8
Cuff Cylindrical 34 Inch ATS Zimmer 15
Cuff Large 42 Inch ATS Zimmer 16

D
Dinamap Mobile Stand 10

E
Edan Patient Monitor Thermal Chart Paper (10 Rolls per Box) - Edan # 01.57.78035 9

G
GE Carescape E-MiniC - C02 Module 4
GE E-PRESTN Carescape Multi Parameter Module 12

M
Mindray Add-on WiFi module material package FRU epm10M / epm12M 6
Mindray ePM AG Module 13
Mindray IBP Module (pn 115-047286-00) includes module 10
Mindray Microstream CO2 Module Passport 8 and 12 5
Mindray MPM Masimo Set Mortara Module SpO2/ ECG/ NIBP/ 2-IBP 2-Temp V 2.0 115-022715-00 4
Mindray Passport and ePM Monitor Paper Thermal - 50mmx20m - Box of 12 (OEM: 0683-00-0505-02) 7
Mindray Passport and ePM Monitor Paper Thermal - 50mmx20m - Single Roll 8
Mindray Sidestream CO2 Module DPM6, DPM7 14
Mindray Sidestream CO2 Module Passport 8 and 12 5
Mindray V12 V21 Nellcor Superpod Module 3

P
Patient Monitor Mobile Stand - Universal 11
Philips InvelliVue MP5 Fixed Wall Mount Kit 3

R
Respironics CO2 Monitoring Module 2
Respironics CO2 Monitoring Module Only M2741A 9

S
Stryker 940 Cast Cutter - Stainless Steel Blade (#STB043) 6
Stryker 940 Cast Cutter - Titanium Nitrided Blade (#STB045) 7

V
Vetland V1400 Series Optional EtCO2 Module 2

Z

Zimmer ATS Dual Blue Hose 14