TABLE OF CONTENTS

Centrifuges 2
StatSpin Express 4 Centrifuge 3
Serofuge 2002 Centrifuge 3



CENTRIFUGES

Auxo Medical offers a curated selection of centrifuges designed for various laboratory and medical applications. Our inventory features the StatSpin Express 4 and the Clay Adams Sero-Fuge 2002. These centrifuges are essential tools for separating components of liquids, such as blood, urine, or other biological samples, enabling further analysis and diagnosis.

The StatSpin Express 4 is a high-speed horizontal centrifuge known for its rapid processing capabilities, enhancing laboratory efficiency. It is ideal for applications requiring quick separation of samples, such as blood cell separation for hematology analysis. The Clay Adams Sero-Fuge 2002, on the other hand, is a versatile centrifuge designed for a broader range of applications, including serological testing, blood banking, and general laboratory procedures.

Centrifuges offered through Auxo Medical are known for their reliability, ease of use, and safety features. They are designed to meet the demands of busy laboratories, providing accurate and consistent results. The benefits of using these centrifuges include time-saving efficiency, enhanced sample quality, and improved diagnostic accuracy. Order your centrifuges through Auxo Medical today!

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED



STATSPIN EXPRESS 4 CENTRIFUGE

The **StatSpin Express 4 high-speed horizontal centrifuge** revolutionizes laboratory efficiency by providing rapid turnaround times and enhancing productivity, particularly in Lean laboratory settings. Positioned adjacent to any STAT analyzer, it ensures that serum or plasma gel tubes are ready in mere minutes, thanks to its horizontal spin that creates a flat gel barrier. The device is user-friendly, featuring only three settings—3, 5, or 10 minutes—eliminating the need for speed adjustments or complicated settings. Safety is paramount, with a lid that remains locked until the rotor stops completely and an alarm that signals the end of the cycle. This allows for immediate processing of samples as they arrive, removing the inefficiencies of batching and waiting. Capable of handling up to eight samples simultaneously, the Express 4 significantly boosts throughput. In just three minutes, it delivers high-quality serum or plasma separation, making it ideal for rapid separation of gel tubes, creating flat gel barriers, preparing samples for STAT chemistry, and facilitating single-piece flow for swift sample processing. Key features include its compact, fast, and powerful design for improved productivity, ease of use, universal power supply for global compatibility, maintenance-free brushless motor, and patented rotor design.

Read More
SKU: AM-SS-EX4
Price: $2,890.00



SEROFUGE 2002 CENTRIFUGE

The **Clay Adams Sero-Fuge, Model 2002 Centrifuge**, is a versatile and highly efficient system designed to meet the rigorous demands of blood grouping, cross matching, typing, and other cell washing procedures. This advanced centrifuge from Clay Adams boasts a 12-place head rotor and operates at two speeds, making it adaptable for various laboratory needs. Its superior features include a locking lid for enhanced safety, digital controls and displays for precise operation, and a gasket seal to ensure optimal performance. The dual-speed operation allows for settings of either 2400rpm/500xG or 3400rpm/1049xG, providing flexibility depending on the specific requirements of the procedure. Additionally, the Sero-Fuge 2002 includes an unbreakable 12-place plastic rotor designed to accommodate blood tubes up to 12.7 mm in diameter and 75 mm in length, ensuring durability and reliability. This centrifuge stands out as a top choice for laboratories seeking high performance and versatility in their cell washing and blood testing processes.

LEAD TIME: 6 Weeks

Read More
SKU: AM-SEROFUGE-2002
Price: $2,995.00

INDEX

S

Serofuge 2002 Centrifuge 3
StatSpin Express 4 Centrifuge 2