TABLE OF CONTENTS

PPE 2
Non-Woven 3 Ply Disposable Face Mask with Adjustable Nose Wire Box of 50 3
K95 Disposable fold face masks 3
Reusable Face Masks 3
Branded Ear Savers for Use with Surgical Mask 4
Non-Branded Ear Savers for Use with Surgical Mask 4
Disposable Anti-Fog Lens Medical Face Shield 5
FR1 Face Shield, Velcro style Strap 6



PPE

Auxo Medical provides several **face mask and face shield** options to protect the wearer's body from injury or infection. The hazards addressed by PPE include physical, electrical, heat, chemicals, biohazards, and airborne particulate matter.

PPE protects patients who are at high risk for contracting infections through a surgical procedure or who have a medical condition. This barrier reduces the chance of touching, being exposed to, and spreading germs. In addition, protects those who are exposed to substances or potentially infectious material brought in by visitors and healthcare workers. Removing and disposing of contaminated PPE prevents exposing both the wearer and other people to infection.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED



NON-WOVEN 3 PLY DISPOSABLE FACE MASK WITH ADJUSTABLE NOSE WIRE BOX OF 50

Read More
SKU: AM-XT1
Price: $14.25



K95 DISPOSABLE FOLD FACE MASKS

Qty: 1 box

Lightweight, non-woven polypropylene material
Elastic loops fit comfortably around ears
Nose clip molds mask securely to face
Fluid resistant and disposable
FDA Certified Facility
20 Per Box

Read More
SKU: AM-62280-KN95
Price: $14.95



REUSABLE FACE MASKS

This face mask is one way to combat the virus. Social distancing, hand washing, and other CDC recommendations should also be followed.
This face covering does not replace the N95 used in healthcare settings and should not be used when caring for those with COVID-19.
Surgical masks provide protection against large droplets, splashes or sprays of bodily or hazardous fluids. They do not provide the wearer with reliable protection from inhaling smaller airborne particles and are not considered respiratory protection.

Features

- Soft 3-ply cotton fabric
- Antimicrobial finish that contains silver and copper
- Machine washability tested up to 15 times before discarding. Tumble dry on high.
- Sold in quantities of 50

Read More
SKU: AM-TRFM-100
Price: $241.30



BRANDED EAR SAVERS FOR USE WITH SURGICAL MASK

Wearing a mask isn’t fun to begin with but wearing it all day multiple days in a row can quickly create a lot of pain from rubbing on ears. These ear savers provide an alternative to attach the mask to behind the head, saving the ears of our frontline workers!

Fully Customizable! Design with YOUR company's branding!
Sold in quantities of 500, 1000, 2000, 5000 or 10000

Read More
SKU: AM-ITBES-2000
Price: Request Quote



NON-BRANDED EAR SAVERS FOR USE WITH SURGICAL MASK

Wearing a mask isn't fun to begin with but wearing it all day multiple days in a row can quickly create a lot of pain from rubbing on ears. These ear savers provide an alternative to attach the mask to behind the head, saving the ears of our frontline workers!

Lightweight, durable, comfortable plastic ear saver 5.5"L x 1"W (loop points 2" W, 1.5" apart) White

Sold in quantities of 500, 1000, 2000, 5000 or 10000

Read More
SKU: AM-ITES-500, AM-ITES-1000, AM-ITES-10000, AM-ITES-5000, AM-ITES-2000
Price: Request Quote



DISPOSABLE ANTI-FOG LENS MEDICAL FACE SHIELD

The Face Shield is designed for single-use and as a personal protective shield for medical use only. The anti-fog lens makes for uncompromised peripheral, upward, and downward vision. The lightweight design fastens with hook-and-loop and consists of a 100% nylon strap and 1" polyurethane foam head cushion. This face shield is to be used in accordance with CDC guidelines for emergency use only. Compatible with the N95 type mask.

Features:

Single-Use
Anti-Fog Lens
Authorized by FDA under an EUA
Lightweight design
The face shield is to be used in accordance with CDC guidelines for emergency use only.
Compatible with N95 type masks.

Do not use in the presence of high-intensity heat source or flammable gas
Shield is not impact resistant and should not be used as eye protection

Materials:

100% nylon strap
1" polyurethane foam head cushion

MASK AND GOGGLES NOT INCLUDED - SHIELDS ONLY

Sold in cases of 25 or 100

Read More
SKU: AM-FS-1000-AF
Price: $560.84



FR1 FACE SHIELD, VELCRO STYLE STRAP

Qty 25 disposable Shields per case

Sold in cases of 25

Introducing the FR1 Face Shield – made and assembled entirely in the USA, in partnership with a premium US manufacturer with extensive health care experience.

FDA approved wide lens material with uncompromised peripheral, upward and downward vision.

Lightweight design that fastens securely to your head ensuring the face shield moves with you. Industrial foam padding for increased comfort. The FR1 shield is to be used in accordance with CDC guidelines for emergency use only.

Compatible with N95 type face masks (not included). Ensure you are using respiratory protection, this device is NOT a respirator.

Read More
SKU: AM-FR1-V-25
Price: $197.15

INDEX

B
Branded Ear Savers for Use with Surgical Mask 3

D
Disposable Anti-Fog Lens Medical Face Shield 4

F
FR1 Face Shield, Velcro style Strap 5

K
K95 Disposable fold face masks 2

N
Non-Branded Ear Savers for Use with Surgical Mask 4
Non-Woven 3 Ply Disposable Face Mask with Adjustable Nose Wire Box of 50 2

R
Reusable Face Masks 3