
TABLE OF CONTENTS

Rentals 3
GE Dash 3000 Patient Monitor - RENTAL 4
Dinamap Pro 400 V2 Vital Signs Monitor - RENTAL 4
GE DINAMAP Procare 400 Vital Signs Monitor - RENTAL 4
Mindray Accutorr V Vital Signs Monitor - RENTAL 5
Mindray Spectrum Patient Monitor - RENTAL 5
Welch Allyn Spot Vital Signs Model 420 NTB - RENTAL 6
Mindray DPM 4 - RENTAL 6
Mindray Passport V with CO2 - RENTAL 7
Mindray Passport V - RENTAL 7
Datascope Passport 2 - RENTAL 8
Steris Amsco Eagle 3043 - RENTAL 8
Steris Amsco Eagle 3053 - RENTAL 9
Steris Amsco Century V160 - RENTAL 9
Steris Amsco Century V148 - RENTAL 10
Steris V116 Steam Sterilizer - RENTAL 10
Steris Eagle 3000 Small Series - RENTAL 11
Hill-Rom P8000 - RENTAL 11
Midmark 510 Stretcher - RENTAL 12
Hill-Rom P8050 OB-GYN Stretcher - RENTAL 12
Steris Harmony LC 500 OR Lights - RENTAL 13
Zimmer ATS 2000 - RENTAL 13
Zimmer ATS1200 Single Channel Tourniquet System - RENTAL 14
Zimmer ATS 3000 - RENTAL 14
Steris Caviwave Single Chamber Ultrasonic Cleaning System - RENTAL 15
Durasonic DS20L Ultrasonic Cleaner - RENTAL 15
Steris Analog Tabletop Warmer - RENTAL 16
Steris QDJ04 Dual Compartment Blanket Warmer - RENTAL 17
Amsco QDJ01 Tabletop Warming cabinet - RENTAL 17
Steris Synergy Washer - RENTAL 17
Steris Reliance 444 - RENTAL 18
Datex-Ohmeda S5 Compact - RENTAL 19
Skytron 3500 Elite - RENTAL 19
Skytron 6500 - RENTAL 19

Steris Amsco 3080 - RENTAL 20
Steris Amsco Model 3085 - RENTAL 20
Medfusion 3500 V4 Syringe Infusion Pump - RENTAL 21
Braun Vista Basic - RENTAL 21
Baxter Infus O.R. - RENTAL 22
Baxter Flo-Gard 6201 - RENTAL 22
Baxter AS50 Syringe Pump - RENTAL 23
Alaris Medley 8110 Syringe Pump Module - RENTAL 24
Alaris Medley 8100 - RENTAL 24
Alaris Medley 8015 - RENTAL 25
Hill-Rom Advance Series - RENTAL 25
Hill-Rom Advanta Series - RENTAL 26
Hill-Rom TotalCare Bed - RENTAL 26
Hill-ROM Affinity 3 Birthing Bed - RENTAL 27
GE Corometrics 120 Fetal Monitor - RENTAL 27
GE/Corometrics 172 Twins Fetal Monitor - RENTAL 28
Midmark 647 Power Podiatry Chair - RENTAL 28
Midmark Ritter 75E - RENTAL 29
Midmark Ritter 75L - RENTAL 30
Midmark Ritter 222 Barrier-Free Exam Table - RENTAL 30
Midmark 204 Exam Table - RENTAL 31
Valleylab Force 2 - RENTAL 31
Valleylab Force FX & Force FXC - RENTAL 31
Burdick 8500 ECG Machine - RENTAL 32
GE MAC 1200 - RENTAL 32
Mortara ELI 250 ECG Machine - RENTAL 33
Lifepak 20E Defibrillator - RENTAL 33
LIFEPAK 20 - RENTAL 34
Zoll M Series 3 lead Biphasic, Pacing AED - RENTAL 34
Steris Powerpack Chomalox Boiler - RENTAL 35
Ritter M11 Autoclave, M11-001 - RENTAL 35
Ritter Midmark M9, M9-001 - RENTAL 36
Midmark Ritter M11 - RENTAL 36
Tuttnauer 3870M - RENTAL 37
Datascope Passport V w/ Gas Module III - RENTAL 37
CardioCap 5 with 5 Agent Gas with Printer - RENTAL 38
Datascope Passport 2 With 5 Gas Module - RENTAL 38
Zoll AED Plus - RENTAL 39
GE Datex-Ohmeda S/5 Aespire - RENTAL 39






RENTALS

Auxo Medical understands that medical equipment can come at an expensive price which is why we offer our customers the opportunity to test out our high quality equipment rentals. We also understand that medical equipment rental isn't always just about finances. Sometimes, customers wish to try out a piece of equipment before investing in it so rentals seem ideal. Our rentals consist of name brands like Steris, Mindray, Zimmer, Skytron, and so much more.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED



GE DASH 3000 PATIENT MONITOR - RENTAL

The GE Dash 3000 monitor comes with a large color screen and was designed for any hospital or clinical setting. The GE Dash 3000 is available with the following options

Read More
SKU: AM-GE-DASH3000-RT
Price: Request Quote



DINAMAP PRO 400 V2 VITAL SIGNS MONITOR - RENTAL

Redesigned for quality-driven healthcare professionals who demand proven world-class technology in every parameter, the GE Medical DINAMAP Pro 400 V2 brings together proven DINAMAP algorithms in a monitor with enhanced reliability features, and the durability you expect from a GE Dinamap Vital Signs Monitor.

Read More
SKU: N/A
Price: Request Quote



GE DINAMAP PROCARE 400 VITAL SIGNS MONITOR - RENTAL

The DINAMAP ProCare 400 monitor enhances productivity with easy to use features and an intuitive menu structure. With a properly-sized BP cuff in place, just press a button to take a BP reading. It’s that easy. All DINAMAP ProCare Monitors are fully upgradeable so they can grow with your changing needs. With GE Healthcare you get more than a monitor. You get a strategic partner and a total solution.

Read More
SKU: AM-GE-DPC400-U
Price: Request Quote



MINDRAY ACCUTORR V VITAL SIGNS MONITOR - RENTAL

The new Accutorr V vital signs monitor is an easy to use, non-invasive blood pressure monitor with more features and capability. Every unit comes equipped with a recorder, lithium ion battery and your choice of SpO2 technology including Masimo SET®, or Nellcor® OxiMax® SpO2. The Accutorr V is capable of storing up to 1200 measurements easily accessible by patient ID through bar code technology. Plus using the Accutorr CS Charting Solution you can conveniently connect the Accutorr V to your facility’s EMR. If maximum flexibility, high value, and reliability match your vital signs needs, then look no further than the Accutorr V monitoring solution.

Read More
SKU: N/A
Price: Request Quote



MINDRAY SPECTRUM PATIENT MONITOR - RENTAL

With 8 selectable waveforms and 3 or 5-lead ECG , the Datascope Spectrum monitor supplies the features you need for critical patient care. A vivid 12.1" display comes with auto adjustable numerics and waveforms for optimal visibility. The standard graph and list trends with optional memory expansion of the Spectrum ensure comprehensive and convenient record keeping. Built for performance and function, the Datascope Spectrum is bundled into the compact and sophisticated monitor you need.

Read More
SKU: N/A
Price: Request Quote



WELCH ALLYN SPOT VITAL SIGNS MODEL 420 NTB - RENTAL

The Welch Allyn Spot Vital Signs Model 420 NTB monitor provides vital signs in seconds with hospital-grade technology that’s easy to use, accurate and affordable. Measures noninvasive blood pressure (NIBP), pulse rate, SpO2 & temperature

Read More
SKU: N/A
Price: Request Quote



MINDRAY DPM 4 - RENTAL

The DPM 4 monitor is a compact, mobile anesthesia companion made for the versatile clinic. Continuous patient monitoring on a high-resolution color display allows for uninterrupted anesthesia care, no matter how many places it's needed. High-efficiency is of paramount importance to anesthesia clinicians and the DPM 4 was created for that purpose. Clear, precise data transmission and easy integration with current anesthesia equipment means more patient time and less mechanical configuration.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



MINDRAY PASSPORT V WITH CO2 - RENTAL

Passport V is built on a solid foundation of experience and knowledge. It follows a proven approach to acquisition and presentation of patient data, thereby delivering impressive reliability and performance. Enhanced capability, more optional features and extraordinary value make Passport V a compelling alternative for any healthcare environment requiring continuous physiological monitoring.

Download Brochure
Read More
SKU: AM-PASSPORT5-CO2-SS-RT
Price: Request Quote



MINDRAY PASSPORT V - RENTAL

The Passport V offers a comprehensive range of options that make it well-suited for most hospital departments. The flexible feature-set includes your choice of CO 2 and SpO2 technologies. Additional offerings such as multiple invasive pressures, external device interfacing and a secondary display solution are also available for increased monitoring of your more critical patients.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



DATASCOPE PASSPORT 2 - RENTAL

The Passport 2 features an incredible intelligent display interface which allows the user to evaluate patient status at a glance. Data is presented clearly, concisely and logically. Waveforms are aligned to digital data for immediate interpretation. Graph and list displays are easily configured to personal preference. The 4 Auto-Set Alarms program is quickly configured and extremely easy to use.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



STERIS AMSCO EAGLE 3043 - RENTAL

The Amsco Steris Eagle 3043 offers an economical option to meet the needs of today’s OR, bulk processing in Central Services, and Ambulatory Surgery Centers. When quality is paramount but budgets are constrained, you can count on STERIS. Quality versus cost is no longer a choice that has to be made; enjoy both with STERIS Certified, Pre-Owned sterilizers featuring Amsco® Eagle® steam sterilizer technology. Designed for fast, efficient sterilization of heat- and moisture-stable materials, prevacuum sterilizers are equipped with prevacuum, gravity, flash, express, leak test, and daily air removal test cycles.

Read More
SKU: N/A
Price: Request Quote



STERIS AMSCO EAGLE 3053 - RENTAL

Steris Amsco 3053- Completely Remanufactured Amsco Eagle 3000 Stage 3 series, 24"x 36" x 48", single door sterilizer Recess mounted, or optional cabinet package. Model 3053. Pre-vac sterilizer, (3051 variant for gravity operation) with operating temperature of 250-270F

Read More
SKU: AM-3053-DD-RT
Price: Request Quote



STERIS AMSCO CENTURY V160 - RENTAL

Fully Refurbished V160 Century Series Medium steam sterilizer: designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V148 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 26" x 37.5" x 60" chamber, the V160 sterilizer provides enough space fit all of your equipment and sterilize with speed and reputable quality in mind.

Read More
SKU: N/A
Price: Request Quote



STERIS AMSCO CENTURY V148 - RENTAL

The Steris Amsco Century V148 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V148 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. With a 24" x 36" x 48" chamber, the V148 sterilizer provides enough space fit all of your equipment and sterilize with speed and reputable quality in mind.

Read More
SKU: N/A
Price: Request Quote



STERIS V116 STEAM STERILIZER - RENTAL

Steris V116 Steam Sterilizers are designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V116 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



STERIS EAGLE 3000 SMALL SERIES - RENTAL

Eagle® 3000 Stage 3™ Series Sterilizers feature an advanced microcomputer control system, providing latest standard for cycle setup, selection and monitoring. Once settings are made and the cycle is started, microcomputers accurately monitor and control system operations. The control column on this sterilizer includes an easy-to-read vacuum fluorescent display and an integral thermal printer.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

HILL-ROM P8000 - RENTAL



Hill-Rom’s rugged procedural stretcher is the ideal choice for Emergency, OR, PACU, or ASC settings. This stretcher offers unsurpassed mobility, comfort, and quality. Hill-Rom’s innovative Steering Plus™ system allows crisp cornering and control. In addition, the system self adjusts when passing over uneven floor surfaces, preventing vibrations from reaching the patient, as experienced with larger, fixed-wheel systems.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

MIDMARK 510 STRETCHER - RENTAL



Refurbished Midmark 510 Stretcher

Read More
SKU: N/A
Price: Request Quote



HILL-ROM P8050 OB-GYN STRETCHER - RENTAL

With its fast and easy patient positioning, Hill-Rom’s OB-GYN stretcher (P8050) is ideal for use in both maternity and emergency. Built to last year after year, this stretcher offers outstanding mobility, convenience, and quality.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



STERIS HARMONY LC 500 OR LIGHTS - RENTAL

The cost-effective Harmony® LC Surgical Lighting System provides bright halogen lighting and a flexible, easy-to-maneuver suspension system. The Harmony® LC Surgical Lighting system is a contemporary product recently emerging on the aftermarket and is offered refurbished on a limited basis.

Read More
SKU: N/A
Price: Request Quote



ZIMMER ATS 2000 - RENTAL

The built-in safety, reliability, and convenience of the A.T.S.® 2000 Tourniquet frees you to focus on surgery. Dual displays that monitor actual cuff pressure, precise microprocessor controls, and one-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.
The dual-cuff display is easy to read, even from a distance, and touch pad controls are simple to read and use. Fast cuff inflation and deflation save time. A large selection of cuffs, disposable and reusable, fits virtually all patients. The lightweight unit can be conveniently positioned, and positive-locking connectors have angled ports to keep lines clear of the surgical site.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



ZIMMER ATS1200 SINGLE CHANNEL TOURNIQUET SYSTEM - RENTAL

The Zimmer A.T.S. 1200 Tourniquet System provides accurate pressure monitoring for patient safety and confidence. Microprocessor controls and built in safety features with self-diagnostic tests make monitoring reliable in a single port, dual cuff surgical tourniquet system.

Read More
SKU: N/A
Price: Request Quote



ZIMMER ATS 3000 - RENTAL

The A.T.S.® 3000 Automatic Tourniquet System is a dual-port, dual-cuff system with microprocessor controls and dedicated ports for supplying and measuring pressure independently. With the innovative Limb Occlusion Pressure (LOP) feature, the A.T.S.® 3000 Tourniquet combines the latest in advanced tourniquet technology with the well-established Zimmer®tradition of safety, reliability, and convenience.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



STERIS CAVIWAVE SINGLE CHAMBER ULTRASONIC CLEANING SYSTEM - RENTAL

The Caviwave Ultrasonic Cleaning Systems are freestanding units that automatically clean, flush, lubricate, and dry instruments (depending on unit configuration) prior to final wash and disinfection cycles. It provides powerful micromechanical scrubbing frequency of 132kHz which more than triples the commonly used 40kHz.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



DURASONIC DS20L ULTRASONIC CLEANER - RENTAL

These DuraSonic ultrasonic cleaners include four models with sizes ranging from 11 to 25 liters. Each model is constructed using durable industrial style 40kHz transducers. These provide increased cleaning power along with built-in sweep frequency to ensure uniform cleaning activity throughout the bath. The five larger sizes have built-in drains. Each model has a digital timer and mechanical thermostat, with a power switch/fuse in the rear of the machine.

Read More
SKU: N/A
Price: Request Quote



STERIS ANALOG TABLETOP WARMER - RENTAL

The Steris Amsco M70 Series line of blanket warmers was one of the most popular lines ever manufactured. This series features a tabletop version with a rotary knob that controls temperature. The Amsco® M70® also features over-temperature alarms (visual and audible) and a heating indicator (visible). On full size models, the upper compartment of the Amsco® M70® is the same as the tabletop version, but with a separate lower cabinet that is larger in size than the upper. The lower cabinet of the full size model comes as either a non-heated storage compartment or a heated version, which would maintain a constant 120° F temperature for warming of blankets only.



Download Brochure
Read More
SKU: N/A
Price: Request Quote



STERIS QDJ04 DUAL COMPARTMENT BLANKET WARMER - RENTAL

Steris Model QDJ04

Read More
SKU: N/A
Price: Request Quote



AMSCO QDJ01 TABLETOP WARMING CABINET - RENTAL

STERIS's Amsco Warming Cabinets are designed to be more than just blanket warmers; they also store and warm sterile intravenous fluids (IV fluids) and surgical irrigation solutions, as well as other linens, to programmed digital temperature settings for patient applications.
18 Inches Deep x 30 Inches Wide / Countertop
Single Compartment Warming Cabinet

Read More
SKU: N/A
Price: Request Quote



STERIS SYNERGY WASHER - RENTAL

The **Steris Synergy Washer, available as a rental**, is a versatile and efficient solution designed to meet the highest standards of cleaning and disinfection. Featuring double power doors, this washer offers the flexibility of a single-door operation through a specialized programming function, allowing for customized usage without compromising its dual-door capability. Equipped with the advanced Stage 3000 control system, it boasts nine adjustable cycles to cater to various cleaning needs. Users can choose between steam heat options (208, 240, or 480 volts) or an electric heating configuration (480 volts) to suit their specific requirements. The system includes one enzyme and one detergent injection pump, facilitating automatic chemical dosing ranging from 1/8 to 2.0 oz/gal (1.0 to 16 mL/L). Additionally, an injection pump is provided for the automatic addition of a lubricant during the Thermal Rinse phase, ensuring optimal maintenance of instruments. With a remarkably good rack-to-chamber ratio of 96%, this washer ensures efficient space utilization, compliant with EN/ISO 15883 standards. While the stainless steel trim and wall flashing package are not included, essential accessories such as the Transfer Cart and Manifold rack can be purchased separately to enhance functionality.

Read More
SKU: N/A
Price: Request Quote



STERIS RELIANCE 444 - RENTAL

The Steris Reliance® 444 Single-Chamber Washer/Disinfector is a mechanical washer equipped with an Eagle® 3000 Stage 3™ control system. It meets the highest standards of cleaning performance and maximum throughput. The Reliance® 444 is a washer that adapts to a wide variety of loads and operations, and will process more than the capacity of an ordinary washer its size.

Download Brochure
Technical Specs
Read More
SKU: AM-444-RRM-RT
Price: Request Quote



DATEX-OHMEDA S5 COMPACT - RENTAL

A complete monitor for anesthesia, the Datex-Ohmeda S/5 anesthesia monitor is a multi-parameter monitor. The Datex-Ohmeda S/5 anesthesia monitor features a large screen for easy viewing in the operating room. It is an ideal monitor for any setting, including the hospital, surgery center, small OR and recovery. The Datex-Ohmeda S/5 anesthesia monitor is a top of the line unit from the leader in anesthesia monitoring.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



SKYTRON 3500 ELITE - RENTAL

The Skytron 3500 Elite features an incredible 21" top slide, 500 lb lift/articulation, return-to-level and return-to-center. An ultra modern operating table, able to satisfy all of your C-Arm and general surgery needs.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

SKYTRON 6500 - RENTAL



The Skytron® 6500 Surgical Table is a general surgery table filled with features. The 500 lb. (Elite) and 850 lb. (Hercules) weight capacity and the C-Arm compatibility allow this OR table to be put into any operating room. This table can also adjust to the trendelenburg, back, side, foot, and flex positions.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

STERIS AMSCO 3080 - RENTAL



The Amsco® 3080 Surgical Table features a radiolucent top and can be adjusted to the same variety of positions as the 2080 model. It comes standard with a remote control for easy and efficient patient positioning. This is one of the most durable, flexible, and reliable surgical tables available.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

STERIS AMSCO MODEL 3085 - RENTAL



The Amsco 3085 General Surgical Table is one of the most reliable and popular surgical tables. Recognized as the industry standard due to its sheer versatility and reliability, the Steris Amsco 3085 has been trusted for years by healthcare facilities across the world. The familiar, easy to use hand controls quickly perform table articulations for any patient posture needed.

Download Brochure
Technical Specs
Read More
SKU: N/A
Price: Request Quote

MEDFUSION 3500 V4 SYRINGE INFUSION PUMP - RENTAL



The Medfusion 3500 Syringe Pump is an innovative pump designed for any department within the hospital or physicians practice. The refurbished Medfusion 3500 Syringe Pump can intuitively detect different types of syringes to add greater safety and accuracy to drug administration. This "smart pump" accepts seven different manufacturers syringes, is ergonomically friendly with single hand loading, easily mounts to IV pole or bassinet, and uses FlowSentry technology for patient safety and drug accuracy. This syringe pump is sold refurbished and include a 6-month warranty.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



BRAUN VISTA BASIC - RENTAL

Braun® understands that alternate site clinicians and patients need a pump that is cost-effective, simple to operate, and has the ability to do more than just continuous infusions. The Vista® Basic was designed to ensure simplicity in all aspects of its use: from set installation, to programming, to special mode use and alarm correction. Based on a design that has sols over 80,000 units worldwide, the Vista® Basic provides high reliability packaged in one of the easiest pumps to use.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



BAXTER INFUS O.R. - RENTAL

The Infus O.R.™ syringe infusion device aids in administering many intravenous agents given during anesthetic procedures. It conveniently delivers narcotics, muscle relaxants and vasoactive drugs via SMART LABELS – a system that displays appropriate infusion parameters and rates for specific medication, helping decrease the likelihood of calculation errors in fast-paced surgical environments.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



BAXTER FLO-GARD 6201 - RENTAL

The Flo-Gard® 6201 Volumetric Infusion Pump from Baxter® can deliver a wide variety of fluids over a broad range of infusion rates. Uses only standard Baxter®solution administration sets. Configurable settings allow the pump to be custom-tailored to best suit the hospital’s needs. Pumps a wide variety of fluids, including blood and fat emulsions. Occlusion sensors detect both upstream and downstream restrictions.

The Programmed Delivery Profile (PDP) enables manual programming or ramped calculation of up to 10 sequential infusion programs for situations where multiple flow raters are indicated. Slide clamp loading option provides flow shutoff when I.V. set is removed from the pump. Locking control panel helps prevent patient tampering. Volume-Time programming automatically calculates flow rate. An automatic self-test routine checks for proper function before use. Five-hour memory retains infusion data after power-off.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



BAXTER AS50 SYRINGE PUMP - RENTAL

The Auto Syringe AS50 Infusion Pump is designed to meet the fluid and drug delivery requirements of today's changing clinical environment. It provides accurate, continuous or intermittent infusion of intravenous solutions, drug solutions, whole blood, and packed red blood cells. The pump is indicated for infusion via intravenous (IV), intra-arterial (IA), epidural, or subcutaneous routes of administration. It can be piggybacked into an ongoing infusion line to deliver secondary solutions automatically, precisely, and economically.

Infusion rates are programmable from 0.01 to 438 mL/hr. The pump accepts standard disposable syringes from 1 mL to 60 mL in size. A numeric keypad simplifies programming and makes the pump easy to use. During programming, the status display prompts the user to enter data about the syringe and the infusion. While the pump is running, the status display shows information about the infusion. Safety and effectiveness are reinforced by pre-programmable bolus operation, titration of a dose without interruption of fluid flow, and easily understood alarm and alert messages. The pump includes an optional serial communication feature, which provides remote computer control and monitoring capabilities. A standard 9-pin connector allows an RS-232 communication link between the pump and the host computer.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



ALARIS MEDLEY 8110 SYRINGE PUMP MODULE - RENTAL

The Alaris Syringe module is a syringe pump designed to work with the Alaris PC unit. It combines precision instrumentation with the safeguards of the Guardrails Suite MX software, helping ensure that the right dose or rate is delivered to your patient.The syringe module is intended for facilities that use infusion for the delivery of fluids, medications, blood and blood products using continuous or intermittent delivery through clinically acceptable routes of administration (e.g. intravenous, intra-arterial or enteral).

The syringe module uses a unique pressure sensing disc which shortens the time to alarm, accurately monitors in line pressure, reduces delays upon infusion start up and decreases inadvertent bolus following occlusion release.

Read More
SKU: N/A
Price: Request Quote



ALARIS MEDLEY 8100 - RENTAL

The Alaris® Pump module is designed for the delivery of fluids, medications, blood and blood products using continuous or intermittent delivery for adult, pediatric or neonatal patients. Infusion delivery can occur via clinically acceptable routes of administration (e.g., intravenous, intra-arterial, subcutaneous, epidural, enteral or irrigation of fluid spaces). It is the only large volume pump available today that can help protect every infusion type, including continuous and intermittent deliveries. Clinicians can attach up to four pump modules to a single Point of Care (PC) unit allowing four different infusions. The pump modules can be attached to either side of the PC unit to create the desired configuration. The Alaris® Pump module supports the Guardrails® Suite MX software to reduce the risk of medication errors.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



ALARIS MEDLEY 8015 - RENTAL

The Alaris® System is a modular system intended for adult, pediatric and neonatal care in today’s growing professional healthcare environment. It consists of the PC Unit, the Guardrails® Suite MX, and up to 4 detachable infusion and/or monitoring modules (channels). The Auto-ID Module can be included as a fifth module. The Alaris® System supported by this DFU uses a new generation PC Unit (Model 8015) which provides wireless connectivity right out of the box, and an enhanced display (including color) to clearly communicate critical programming, infusion, monitoring and hospital-defined policy information. Alaris® System wireless communication makes it easier than ever to increase the safety of iv medication and continuously improve clinical best practices, regardless of existing wireless infrastructure.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

HILL-ROM ADVANCE SERIES - RENTAL



The Hill-Rom® Advance® Series hospital bed is perfect for hospitals & surgery centers. This bed offers a variety of features that are perfect for your patient care needs.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

HILL-ROM ADVANTA SERIES - RENTAL



The Hill-Rom® Advanta® Series hospital bed is perfect for hospitals & surgery centers. This bed offers a variety of features that are perfect for your patient care needs.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

HILL-ROM TOTALCARE BED - RENTAL



The TotalCare® P1900 therapy bed delivers the versatility to also respond to the needs of the immobile patients with a moderate to high risk for pressure ulcer development. The best strategy to avoid immobility related complications is to get patients mobile as early as possible. Aside from the Progressive Mobility™ support coming with the TotalCare® bed functionality.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

HILL-ROM AFFINITY 3 BIRTHING BED - RENTAL



In labor and delivery, time and safety are top of mind and every second counts. We have designed every detail of the Affinity 3 birthing bed to help make the birthing process safe for the patient and caregiver. Additional back saving features like the no lift Stow and Go™ foot section and ergonomic push handles make the Affinity 3 a smart choice for caregivers.

Read More
SKU: N/A
Price: Request Quote



GE COROMETRICS 120 FETAL MONITOR - RENTAL

Corometrics 120 maternal/fetal monitors help you manage the labor, delivery and recovery process for every patient. Every 120 series monitor features the clinical excellence that has set the standard in M/F monitoring for more than 30 years. Patented Smart BP, Heart Beat Coincidence, Spectra Alerts, DINAMAP non-invasive blood pressure algorithm and 9-crystal pulse Doppler ultrasound transducers... You can trust the data as an accurate assessment of mom and baby's progress. Labor can be unpredictable. Be prepared with a monitor that can help you manage the entire process - even in the most difficult cases. Smart BP feature - delays auto determination during contractions Vital Signs History feature saves maternal information in memory for quick review or print summary Recovery mode sets the recorder to standby allowing delivery to print only the vital signs and save paper

Read More
SKU: N/A
Price: Request Quote



GE/COROMETRICS 172 TWINS FETAL MONITOR - RENTAL

Whether you perform antepartum fetal monitoring in an office setting or in the hospital, Corometrics* 170 Series Fetal Monitor can help. You can trust our fetal monitors to help you improve the standard of care for your expecting mothers. We designed these compact, lightweight monitors to help you effectively and efficiently care for mom and baby by accurately and noninvasively tracking baby’s heart rate for non-stress tests.

The 170 series is used with our Nautilus* watertight and lightweight Toco and 9-crystal Ultrasound transducers for comfortable and reliable monitoring.

Read More
SKU: N/A
Price: Request Quote



MIDMARK 647 POWER PODIATRY CHAIR - RENTAL

The 647's new, ergonomic foot section design gives physicians what they need at the point of care to perform even the lengthiest procedures. Touch controls and sealed foot control are standard on this model, as are covered glides on the foot section and debris tray. You'll be surprised at how easy it is to keep your chair clean and operating smoothly.

Download Brochure
Read More
SKU: N/A
Price: Request Quote





MIDMARK RITTER 75E - RENTAL

This surgical chair comes standard with power adjustment, drain pan, articulating head rest, adjustable stirrups, accessory rails, foot section lockable in four procto positions, paper roll holder, paper straps, and quick release upholstery tops. Colors may vary.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



MIDMARK RITTER 75L - RENTAL

To help keep you and your staff fresh and alert all day, the Midmark® 75L does the heavy work. Full eight-way power, essential low height, unprecedented patient comfort and a multitude of efficiency features significantly improve your practice, making the Midmark® 75L the table of choice. With a low height of 19 inches and an exclusive Premium Comfort System™, the 75L provides easy access for patients who are pregnant, elderly, or have disabilities. This drives efficiency into your practice by enabling patients to access the table themselves, while helping you and your staff position patients quickly and efficiently. Colors may vary.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



MIDMARK RITTER 222 BARRIER-FREE EXAM TABLE - RENTAL

The days of limited exam rooms for disabled patients are gone. With Ritter® Barrier-Free® tables, more of your exam rooms can be accessible to your special needs patients. This not only puts your patients at ease, but it also makes your practice more efficient.

Download Brochure
Read More
SKU: N/A
Price: Request Quote





MIDMARK 204 EXAM TABLE - RENTAL

The Ritter 204 Manual Exam Table features the latest in design, comfort and functionality in a manual exam table. It provides an enlarged footstep for easier access to make patients feel secure, three side drawers for optimal supply management, modern styling to coordinate with your office, an ergonomic top design to help patients feel comfortable, and additional weight capacity to handle heavy patients.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



VALLEYLAB FORCE 2 - RENTAL

The Force™ 2 ESU electrosurgical cautery unit provides cutting and coagulation in monopolar or bipolar outputs. Provides flexibility through varying degrees of hemostasis by allowing you to choose one of three preset modes. Permits two surgeons to fulgurate from a single generator, for added convenience and efficiency. Accessories are activated only when keyed, reducing the potential for injury caused by inadvertent activation. Ensures controlled, precise desiccation with less destruction of peripheral tissue, making the Force™ 2 generator ideal for laparoscopic procedures.

Download Brochure
Read More
SKU: N/A
Price: Request Quote





VALLEYLAB FORCE FX & FORCE FXC - RENTAL

The Force FX™ generator provides surgeons with the advantage of accurately delivering the selected power setting, even through resistant tissue. Instant Response technology provides surgeons with improved performance at lower power settings, minimizing the risk of tissue damage and neuromuscular stimulation, and decreasing the need to “turn up the generator.” The Force FX™ generator adjusts automatically, responding to tissue changes, maintaining power delivery, and minimizing drag. Capacitive coupling is reduced by 30-50% when using Instant Response™ technology. This reduction is achieved by limiting the RMS voltage and the high frequency harmonics.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



BURDICK 8500 ECG MACHINE - RENTAL

The Burdick 8500 ECG / EKG shows you what you need, comfortably. Burdick 8500 ECG equipment has a friendly user interface, multi-angle 7-inch widescreen, and color display that makes reading waveform data easy, and protocol choices streamline your workflow.

Read More
SKU: AM-B8500-U-RT
Price: Request Quote



GE MAC 1200 - RENTAL

The MAC* 1200 digital, 12-lead electrocardiograph system offers comprehensive ECG solutions with practical features that help meet the needs of hospitals, clinics, office based practices and clinical trials. Its advanced algorithm capabilities, seamless connectivity to the CS and MUSE* Cardiology Information Systems and easy-to-use features provide the quality and convenience users require.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



MORTARA ELI 250 ECG MACHINE - RENTAL

Compact and lightweight, the ELI 250c Electrocardiograph provides comprehensive functionality in an easy-to-use, portable device. Multiple connectivity options ensure patient data transfers easily from the ECG to the location of your choice. Experience the benefits of bidirectional communication via internal modem, LAN or wireless LAN with Mortara's E-Scribe™ and Athena products, as well as third-party EMR, HIS, and PACS systems via XML, PDF and DICOM®. Best 10 technology helps you capture the best 10 seconds of ECG data from system memory, reducing the need for repeat ECGs. Wireless data acquisition helps increase workspace flexibility by gathering 12-lead data without the need for a tethered cable.

Read More
SKU: N/A
Price: Request Quote



LIFEPAK 20E DEFIBRILLATOR - RENTAL

Clinically advanced and packed with power, the LIFEPAK® 20e defibrillator/monitor is highly intuitive for first responders, and also skillfully combines AED function with manual capability so that ACLS-trained clinicians can quickly and easily deliver advanced therapeutic care.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



LIFEPAK 20 - RENTAL

The LIFEPAK® 20 defibrillator/monitor is the ideal crashcart defibrillator putting early, effective defibrillation into the hands of BLS responders while transitioning easily to full therapy and monitoring capabilities for the ALS-trained code team. The 20 is highly intuitive, making it easy for infrequent AED-trained responders to quickly understand and use. With an easy push of a latch, the 20 automatically converts to a manual defibrillator, including waveform displays. In manual mode, all of the advanced cardiovascular life support tools become available. The 20offers noninvasive pacing, Masimo SET® pulse oximetry (SpO2), ECG monitoring (3- or 5-wire) and synchronized cardioversion.

Download Brochure
Read More
SKU: N/A
Price: Request Quote


ZOLL M SERIES 3 LEAD BIPHASIC, PACING AED - RENTAL



The smallest, lightest, most advanced defibrillator monitor ever made, with all the capabilities you'll ever need! The Zoll M Series is Only 11 pounds (5kg). Half the weight and a third the size of any unit in its class, the Zoll M Series Defibrillator has everything you need for cardiac resuscitation.

Download Brochure
Read More
SKU: N/A
Price: Request Quote

STERIS POWERPACK CHOMALOX BOILER - RENTAL



Operating Power: 30KW x
Available Voltage: 208V/480V 3 phase

Read More
SKU: N/A
Price: Request Quote



RITTER M11 AUTOCLAVE, M11-001 - RENTAL

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



RITTER MIDMARK M9, M9-001 - RENTAL

The M9 autoclave couldn't be easier and is time tested for reliability. You simply load it - set it - and forget it. Select your cycle and the time is set for you. Press "Start" and the M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you fast drying time for your instruments.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



MIDMARK RITTER M11 - RENTAL

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



TUTTNAUER 3870M - RENTAL

This is Tuttnauer's largest manual table-top autoclave. Believe it or not there is still a demand for the 3870m analog unit among many healthcare facilities and veterinary practices. This beefy, no frills autoclave covers all the basics needed for proper sterilization. It can always be counted on even under harsh conditions. The large 15" x 30" chamber accommodates longer instruments or a large number of smaller instruments, making it a perfect match for spay and nueter clinics, equine practices, busy animal hospitals and multi-practice medical offices.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



DATASCOPE PASSPORT V W/ GAS MODULE III - RENTAL

The Passport V offers a comprehensive range of options that make it well-suited for most hospital departments. The flexible feature-set includes your choice of CO 2 and SpO2 technologies. Additional offerings such as multiple invasive pressures, external device interfacing and a secondary display solution are also available for increased monitoring of your more critical patients.

Read More
SKU: N/A
Price: Request Quote



CARDIOCAP 5 WITH 5 AGENT GAS WITH PRINTER - RENTAL

GE™ CardioCap 5 monitor is a compact, all-in-one monitor designed for use in the operating room, ambulatory surgery unit, induction room and PACU. The Cardiocap/5 offers full vital signs monitoring capabilities from the monitoring of oxygenation and circulation to advanced airway gas analysis and ventilation measurement with Patient Spirometry."

Download Brochure
Read More
SKU: N/A
Price: Request Quote



DATASCOPE PASSPORT 2 WITH 5 GAS MODULE - RENTAL

The Gas Module II adds state-of-the-art gas monitoring and analysis capabilities to the Datascope vital signs monitors. The Gas Module II with optional Spirometry measures airway pressures and lung volumes directly at the patient’s airway. This data, along the loops and numeric data that you need can be displayed on the Spectrum OR™. The Gas Module II is a breath-by-breath gas analyzer, designed to meet the comprehensive anesthesia monitoring requirements of virtually every hospital and freestanding surgical center - whatever its size, specialty, or patient base.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



ZOLL AED PLUS - RENTAL

When a cardiac arrest occurs, the fact is that only half of the victims will need a shock. But nearly all will require CPR. You deserve an AED that helps you all the time. And, only one AED can actually see when you are doing CPR and help you do it well. You need more than just commands, without assistance. That’s not smart, and it’s certainly not help. Zoll’s AED Plus® features Real CPR Help, a CPR feedback tool that is able to actually see what you are doing and provide feedback to help you do it well. Audio and visual prompts help you rescue with confidence and clarity unmatched by any other AED.

Download Brochure
Read More
SKU: N/A
Price: Request Quote



GE DATEX-OHMEDA S/5 AESPIRE - RENTAL

The S/5 Aespire® was developed using GE’s unique approach to anesthesia workstation development – close and continuous collaboration with clinicians. It combines the proven performance of the 7100 Ventilator or 7900 Ventilator with Datex-Ohmeda® monitoring technology into an affordable system solution designed for clinicians throughout the world. The S/5 Aespire® utilizes the Advanced Breathing System, the ABS™. It’s design, based on your needs and expertise, is fully integrated into the workstation. Integrated electronics detect complete engagement or disengagement of the ABS™ while fewer parts and connections reduce the potential for leaks and misconnects, providing greater patient safety.

S5 Aespire 7100 Brochure
S5 Aespire 7900 Brochure
Read More
SKU: N/A
Price: Request Quote

INDEX

A

Alaris Medley 8015 - RENTAL 23
Alaris Medley 8100 - RENTAL 23
Alaris Medley 8110 Syringe Pump Module - RENTAL 22
Amsco QDJ01 Tabletop Warming cabinet - RENTAL 15

B

Baxter AS50 Syringe Pump - RENTAL 21
Baxter Flo-Gard 6201 - RENTAL 21
Baxter Infus O.R. - RENTAL 20
Braun Vista Basic - RENTAL 20
Burdick 8500 ECG Machine - RENTAL 30

C

CardioCap 5 with 5 Agent Gas with Printer - RENTAL 36

D

Datascope Passport 2 - RENTAL 6
Datascope Passport 2 With 5 Gas Module - RENTAL 37
Datascope Passport V w/ Gas Module III - RENTAL 36
Datex-Ohmeda S5 Compact - RENTAL 17
Dinamap Pro 400 V2 Vital Signs Monitor - RENTAL 2
Durasonic DS20L Ultrasonic Cleaner - RENTAL 14

G

GE Corometrics 120 Fetal Monitor - RENTAL 26
GE Dash 3000 Patient Monitor - RENTAL 2
GE Datex-Ohmeda S/5 Aespire - RENTAL 38
GE DINAMAP Procare 400 Vital Signs Monitor - RENTAL 3
GE MAC 1200 - RENTAL 31
GE/Corometrics 172 Twins Fetal Monitor - RENTAL 26

H

Hill-Rom Advance Series - RENTAL 24
Hill-Rom Advanta Series - RENTAL 24
Hill-ROM Affinity 3 Birthing Bed - RENTAL 25
Hill-Rom P8000 - RENTAL 10
Hill-Rom P8050 OB-GYN Stretcher - RENTAL 11
Hill-Rom TotalCare Bed - RENTAL 25

L

LIFEPAK 20 - RENTAL 32
Lifepak 20E Defibrillator - RENTAL 32

M

Medfusion 3500 V4 Syringe Infusion Pump - RENTAL 19
Midmark 204 Exam Table - RENTAL 29
Midmark 510 Stretcher - RENTAL 10
Midmark 647 Power Podiatry Chair - RENTAL 27

Midmark Ritter 222 Barrier-Free Exam Table - RENTAL 28
Midmark Ritter 75E - RENTAL 27
Midmark Ritter 75L - RENTAL 28
Midmark Ritter M11 - RENTAL 35
Mindray Accutorr V Vital Signs Monitor - RENTAL 3
Mindray DPM 4 - RENTAL 5
Mindray Passport V - RENTAL 6
Mindray Passport V with CO2 - RENTAL 5
Mindray Spectrum Patient Monitor - RENTAL 4
Mortara ELI 250 ECG Machine - RENTAL 31

R

Ritter M11 Autoclave, M11-001 - RENTAL 34
Ritter Midmark M9, M9-001 - RENTAL 34

S

Skytron 3500 Elite - RENTAL 17
Skytron 6500 - RENTAL 18
Steris Amsco 3080 - RENTAL 18
Steris Amsco Century V148 - RENTAL 8
Steris Amsco Century V160 - RENTAL 8
Steris Amsco Eagle 3043 - RENTAL 7
Steris Amsco Eagle 3053 - RENTAL 7
Steris Amsco Model 3085 - RENTAL 19
Steris Analog Tabletop Warmer - RENTAL 14
Steris Caviwave Single Chamber Ultrasonic Cleaning System - RENTAL 13
Steris Eagle 3000 Small Series - RENTAL 9
Steris Harmony LC 500 OR Lights - RENTAL 11
Steris Powerpack Chomalox Boiler - RENTAL 33
Steris QDJ04 Dual Compartment Blanket Warmer - RENTAL 15
Steris Reliance 444 - RENTAL 16
Steris Synergy Washer - RENTAL 16
Steris V116 Steam Sterilizer - RENTAL 9

T

Tuttnauer 3870M - RENTAL 35

V

Valleylab Force 2 - RENTAL 29
Valleylab Force FX & Force FXC - RENTAL 30

W

Welch Allyn Spot Vital Signs Model 420 NTB - RENTAL 4

Z

Zimmer ATS 2000 - RENTAL 12
Zimmer ATS 3000 - RENTAL 13
Zimmer ATS1200 Single Channel Tourniquet System - RENTAL 12
Zoll AED Plus - RENTAL 37
Zoll M Series 3 lead Biphasic, Pacing AED - RENTAL 33