TABLE OF CONTENTS

Vaporizers 2
- Vetland D20 Veterinary Isoflurane Vaporizer 2
- Mindray Desflurane Vaporizer 3
- Mindray Sevoflurane Vaporizer with quick fill adapter 3
- Datex Ohmeda Tec 6 Plus Vaporizer 4
- Datex Ohmeda Tec 7 Sevoflurane Vaporizer 4
- Datex Ohmeda Tec 7 Isoflurane Vaporizer 5
- Datex Ohmeda Tec 5 Sevoflourane Vaporizer 5
- Drager 19.1 Sevoflurane Vaporizer 6
- Drager 19.1 Isoflurane Vaporizer 6
- Drager Vapor 2000 Sevoflurane Vaporizer 7
- Drager Vapor 2000 Isoflurane Vaporizer 7



VAPORIZERS

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

BRAND NEW

RECONDITIONED

REFURBISHED



VETLAND D20 VETERINARY ISOFLURANE VAPORIZER

The **Vetland D20 Veterinary Isoflurane Vaporizer** is a high-precision, agent-specific device designed for continuous-flow veterinary anesthesia systems. Calibrated exclusively for isoflurane, it offers a concentration range of 0–5%, with automatic temperature and flow compensation to ensure consistent delivery regardless of ambient conditions or varying gas flow rates. Its universal Selectatec-style mounting allows for easy integration with a variety of anesthesia machines, with options for funnel-fill or keyed-fill systems—both featuring a large locking cap to prevent spills and incorrect fills. Built for durability, the rugged, corrosion-resistant construction ensures lifelong calibration, backed by a five-year warranty. Compact and lightweight at approximately 5.5 pounds with dimensions of 8” H × 4” W × 6” D, the D20 provides a reliable, precise, and versatile anesthesia vaporizing solution for veterinary use.
Read More
SKU: AM-520-11
Price: Request Quote



MINDRAY DESFLURANE VAPORIZER

Desflurane Vaporizer for Mindray A4 Anesthesia Machine

Read More
SKU: AM-040-001997-00
Price: Request Quote



MINDRAY SEVOFLURANE VAPORIZER WITH QUICK FILL ADAPTER

Includes Quick Fill Bottle Adapter 0004-00-0100
For Mindray A4 Anesthesia Machine

Read More
SKU: AM-0992-00-0148
Price: Request Quote

DATEX OHMEDA TEC 6 PLUS VAPORIZER



The Tec 6 Plus is an electronic vaporizers that heats Desflurane to maintain constant temperature and vapour pressure for consistent output. Tec 6 Plus incorporates an LED display indicating vapouriser status: no output, low agent, warm-up, operational, and alarm battery low.

Read More
SKU: AM-TEC6P-DES-C
Price: Request Quote



DATEX OHMEDA TEC 7 SEVOFLURANE VAPORIZER

The Tec 7 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models. Features Easy-fil, which is designed to simplify the filling process while minimizing leaks. Sevoflurane only. Other agents require a separate vaporizer.

Read More
SKU: AM-TEC7-SEVO
Price: Request Quote



DATEX OHMEDA TEC 7 ISOFLURANE VAPORIZER

The Tec 7 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models. Features Easy-fil, which is designed to simplify the filling process while minimizing leaks. Isoflurane only. Other agents require a separate vaporizer.

Read More
SKU: AM-TEC4-ISO
Price: Request Quote



DATEX OHMEDA TEC 5 SEVOFLOURANE VAPORIZER

The Tec 5 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models.

Features Key Fill or Funnel Fill, which is designed to simplify the filling process while minimizing leaks. Sevoflurane only. Other agents require a separate vaporizer

Read More
SKU: AM-TEC5-SEVO
Price: Request Quote



DRAGER 19.1 SEVOFLURANE VAPORIZER

The vaporizers are designed for out of circuit use in continuous flow techniques of inhalation anaesthesia. Temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow. Each vaporer is clearly labeled with the name of the anaesthetic agent for which it is designed.

Read More
SKU: AM-DR-19.1-SEV
Price: Request Quote



DRAGER 19.1 ISOFLURANE VAPORIZER

The vaporizers are designed for out of circuit use in continuous flow techniques of inhalation anaesthesia. Temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow. Each vaporer is clearly labeled with the name of the anaesthetic agent for which it is designed.

Read More
SKU: AM-DR-19.1-ISO
Price: Request Quote



DRAGER VAPOR 2000 SEVOFLURANE VAPORIZER

The Vapor 2000 is a vaporizer for volatile anesthetic agents with unsurpassed performance. Originally developed for the extreme requirements of mobile army surgical units, these demanding ambient specifications bene t the normal, daily clinical routine in anesthesia.

Read More
SKU: AM-DR2000-SEVO
Price: Request Quote



DRAGER VAPOR 2000 ISOFLURANE VAPORIZER

The Vapor 2000 is a vaporizer for volatile anesthetic agents with unsurpassed performance. Originally developed for the extreme requirements of mobile army surgical units, these demanding ambient specifications bene t the normal, daily clinical routine in anesthesia.

Read More
SKU: AM-DR2000-ISO
Price: Request Quote

INDEX

D

Datex Ohmeda Tec 5 Sevoflourane Vaporizer 5
Datex Ohmeda Tec 6 Plus Vaporizer 3
Datex Ohmeda Tec 7 Isoflurane Vaporizer 4
Datex Ohmeda Tec 7 Sevoflurane Vaporizer 4
Drager 19.1 Isoflurane Vaporizer 6
Drager 19.1 Sevoflurane Vaporizer 5
Drager Vapor 2000 Isoflurane Vaporizer 7
Drager Vapor 2000 Sevoflurane Vaporizer 6

M

Mindray Desflurane Vaporizer 2
Mindray Sevoflurane Vaporizer with quick fill adapter 3

V

Vetland D20 Veterinary Isoflurane Vaporizer 2